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**5.1.2.e**

	<b>SOP - Laboratory Cleaning &amp; Disinfection</b>	Page 1 of 1
---	---	-------------

**Purpose:**

To provide a daily clean workbench in the laboratory.

**Interested Parties:**

Team leaders and lab assistants

**Definitions:**

- a. IPA: Isopropyl alcohol
- b. BSC: biological safety cabinet
- c. Non-disposable equipment: laptop, multichannel pipette, 3D-printed moulds.

**Responsibility:**

Every laboratory worker is responsible for cleaning his/her workbench after the work is done. Clean and disinfect work surfaces and non-disposable equipment before and after any activity, and when a possible contamination is detected.

**Procedure:**

- a. Wear gloves before cleaning and disinfecting any surface.
- b. Prepare and organize the workplace to facilitate cleaning and disinfection
- c. Take paper towels and spray a generous amount of IPA (80%) to disinfect the work area. **Note:** Alternatively, spray directly on contaminated surface and use multiple paper towels to clean and remove contaminants.
- d. Clean the surface of the BSC/workbench
- e. Use dry paper towels to remove excess liquid

**Non-disposable Equipment**

- f. Spray IPA (80%) on a dry paper towel and wipe the high-touch surfaces on the device. **Note:** Do not spray IPA directly on the piece of equipment!
- g. Dry the device using a dry paper towel
- h. Inspect the equipment for remaining smudges and repeat the cleaning procedure

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
PhoenixDx® 2019-nCoV  
RNA Detection Kit

Procomcure GmbH

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PROBE	HERKUNFT	PhoenixDX IVD KONTROLLE
1		PASS
2		PASS
3		PASS
4		PASS
5		PASS
6	RINGVERSUCH	FAIL
7	RINGVERSUCH	FAIL
8	RINGVERSUCH	FAIL
9	RINGVERSUCH	FAIL
10	RINGVERSUCH	FAIL
11	RINGVERSUCH	FAIL
12	RINGVERSUCH	PASS
13	RINGVERSUCH	PASS
14		PASS
15		PASS
16		PASS
17		PASS
18		PASS
19		PASS
20		PASS
21	RNA delivered	PASS
22	RNA delivered	PASS
23	RNA delivered	PASS
24		PASS
25		PASS
26		PASS
27		PASS
28		PASS
29		PASS
30		PASS

PROBE	HERKUNFT	PhoenixDX IVD KONTROLLE
30		PASS
31		PASS
32	RINGVERSUCH	FAIL
33	RINGVERSUCH	FAIL
34		PASS
35		PASS
36		PASS
37		PASS
38	RNA delivered	PASS
39	RNA delivered	PASS
40	RNA delivered	PASS
41	RNA delivered	PASS
42		PASS
43		PASS
44		PASS
45		PASS
46		PASS
47		PASS
48		PASS
49		PASS
50	RNA delivered	PASS
51	RNA delivered	PASS
52		PASS
53		PASS
54		PASS
55		PASS
56		PASS
57		PASS
58		PASS
59		PASS
60		PASS

PROBE	HERKUNFT	PhoenixDX IVD KONTROLLE
61		PASS
62		PASS
63		PASS
64		PASS
65		PASS
66		PASS
67	RNA delivered	FAIL
68	RNA delivered	FAIL
69	RNA delivered	PASS
70		PASS
71		PASS
72		PASS
73		PASS
74		PASS
75		PASS
76		PASS
77		PASS
78	RINGVERSUCH	FAIL
79	RINGVERSUCH	FAIL
80	RINGVERSUCH	PASS
81		PASS
82		PASS
83		PASS
84		PASS
85		PASS
86		PASS
87		PASS
88		PASS
89		PASS
90		PASS
91		PASS

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[Novogenia.com](http://Novogenia.com)

5.1.2.e

## In-house validation and performance characteristics of a high throughput COVID-19 qRT-PCR detection assay using 384-well format and 10µl reaction volume

### 5.1.2.e

Laboratory of the Novogenia GmbH  
Strass 19, 5301 Eugendorf, AUSTRIA

20.12.2020

#### Abstract

The laboratory diagnostic gold standard for the diagnosis of infections with coronavirus SARS-CoV-2 is the direct virus detection from respiratory secretions using a real-time reverse transcriptase (RT) polymerase chain reaction (PCR) test for the detection of SARS-CoV-2. We developed an automated high throughput SARS-CoV-2 detection workflow using samples collected in transport media, a magnetic bead RNA isolation, and a 384-well format RT-qPCR using the FTD SARS-CoV-2 Assay from Fast Track Diagnostics. Performance characteristics were measured using 4 different panels showing the analytical and clinical performance of the workflow. The following data was collected: Analytical sensitivity, analytical specificity, cross-reactivity, and clinical performance.

#### Methods

##### Sample collection

Sensitivity panel: SARS-CoV-2 sample from a 1:100 dilution derived from stock (hCoV-19/Netherlands/NoordBrabant\_10003/2020) in MEM with Hanks' salts; heat-inactivated at 60°C for 2 hours. The sample contains SARS-CoV-2 with a  $5.62 \times 10^5$  TCID<sub>50</sub>/ml and up to  $4.56 \times 10^7$  digital copies of SARS-CoV-2 positive-strand RNA/ml. For the validation experiment, we used a 10-fold dilution series ranging from  $4.56 \times 10^7$  viral copies per millilitre (cp/ml) to 4.56 cp/ml.

Specificity panel: The panel consists of 298 negative samples (111 non-template controls (NTC) controls, 187 extracted confirmed negative clinical samples).

Cross-reactivity panel: 10 simulated clinical samples (1ml) containing either human respiratory seasonal viruses, heat-inactivated SARS-CoV-2 virus, or no virus (a negative control) in MEM with Hanks' salts and Hep2 cells (10.000 cells/ml).

Clinical performance panel: 91 confirmed positive (19) and negative (72) clinical samples collected with nasopharyngeal swabs, oropharyngeal swabs, or a gargling kit. The results have been confirmed by two other CE IVD labelled RT-qPCR assays.

##### RNA extraction

RNA was extracted using the SphaeraMag DNA/RNA Isolation Kit 96-Universal (Procomcure Biotech) on Allsheng Auto-Pure96 automated systems according to the manufacturer protocol with the exceptions of using 300 µL of starting input and eluting in 60 µL elution buffer.

##### RT-qPCR

RT-qPCR was performed using the FTD SARS-CoV-2 Assay (Fast Track Diagnostics) which targets the highly conserved regions within the ORF1ab and N gene. All reactions were performed on a Quantstudio 5 (Thermo Fisher Scientific, Applied Biosystems) according to the manufacturers'

instructions in 10µl reaction volumes on 384-well PCR plates. Clinical samples were confirmed using the SARS-CoV-2 RealFast Assay (Vienna lab) and the PhoenixDX SARS-CoV-2 assay (Procomcure).

## Results

### Analytical sensitivity

The experiment aimed to prove acceptable sensitivity, reproducibility, and limit of detection of the highly automated SARS-CoV-2 workflow using the FTD SARS-CoV-2 Assay. The analytical sensitivity is the ability of the assay to consistently detect a given target sequence in the tested biological specimens. The lowest concentration detectable in greater than or equal to 95% of tested specimens is then defined as the “Limit of Detection” (LoD). We created a serial dilution of the SARS-CoV-2 stock sample achieving final viral copies of 9120000 down to 0 digital copies per reaction. All PCR reactions were performed in quadruplicates and the results including mean Ct values were summarized in table 1.

Table 1:

Dilution	digital copies / ml sample	digital copies per reaction (4µl)			
	$4.56 \times 10^9$	91200000			
10-1	$4.56 \times 10^8$	9120000	Result	Mean Ct	Positive(s)
10-2	$4.56 \times 10^7$	912000	Pos	21.44	(4/4)
10-3	$4.56 \times 10^6$	91200	Pos	23.31	(4/4)
10-4	$4.56 \times 10^5$	9120	Pos	26.72	(4/4)
10-5	$4.56 \times 10^4$	912	Pos	30.26	(4/4)
10-6	$4.56 \times 10^3$	91	Pos	33.24	(4/4)
10-7	$4.56 \times 10^2$	9	Pos	37.48	(4/4)
10-8	$4.56 \times 10^1$	1	Pos	38.87	(1/4)
10-9	4.56	0	Neg	Undeterm.	(0/4)

The data demonstrate robust performance and acceptable sensitivity of the laboratory workflow. Mean Ct values ranging from 21.44 for the highest ( $4.56 \times 10^7$  copies/ml transport medium or  $9.12 \times 10^5$  copies/reaction) to 38.87 for the lowest (45.6 copies/ml transport medium or 1 copy per reaction) concentration of SARS-CoV-2 virus. The method shows a reproducible low limit of detection (LOD) with  $4.56 \times 10^2$  copies per ml transport media or 9 copies per reaction.

### Analytical specificity

Analytical specificity measures the proportion of samples without SARS-CoV-2 that are correctly classified as negative. Therefore, 111 non-template controls (NTC) controls and 187 confirmed negative clinical samples have been tested to evaluate the occurrence of potential non-specific amplification when using the FTD SARS-CoV-2 kit. Analytical specificity is displayed in Table 2.

Table 2:

	Total reactions	Positive results	Analytical sensitivity %	Confidence interval %
Negative clinical samples	187	0	100	98.05 to 100
Non-template control	111	0	100	96.73 to 100
Total	298	0	100	98.77 to 100

Overall, an analytical specificity of 100% was reached.

### Cross-reactivity

The analysis was tested for potential cross-reactions with common human seasonal viruses causing respiratory infections. The panel contained either human respiratory seasonal viruses, heat-

inactivated SARS-CoV-2 virus or no virus. The qPCRs were performed in quadruplicates. Results including mean Ct values are shown in table 3.

Table 3:

Definition	Virus	Concentration	Result	Mean Ct
EQA6_CoV20_01	CoV-NL63		Neg	Undetermined
EQA6_CoV20_02	CoV-229E		Neg	Undetermined
EQA6_CoV20_03	SARS-CoV-2	$4.56 * 10^1$	Neg	Undetermined
EQA6_CoV20_04	SARS-CoV-2	$4.56 * 10^2$	Pos	37.86
EQA6_CoV20_05	CoV-OC43		Neg	Undetermined
EQA6_CoV20_06	SARS-CoV-2	$4.56 * 10^2$	Pos	37.07
EQA6_CoV20_07	Influenza virus A(H3N2)		Neg	Undetermined
EQA6_CoV20_08	No virus		Neg	Undetermined
EQA6_CoV20_09	SARS-CoV-1		Neg	Undetermined
EQA6_CoV20_10	SARS-CoV-2 (d1)	$4.56 * 10^4$	Pos	31.16

No unspecific signal could be detected indicating high specificity and no cross-reactivity of the assay. Sample EQA6\_CoV20\_03 was not detected confirming the limit of detection shown in table 1.

#### Clinical performance

As analytical sensitivity and specificity needs to be differentiated from diagnostic sensitivity and specificity, clinical performance was evaluated by analysing 91 confirmed positive (19) or negative (72) samples. The results are shown in table 4.

Table 4:

	Total number	Percentage	Confidence interval %
<b>Sensitivity</b>	19/19	100	82.35 to 100
<b>Specificity</b>	72/72	100	95.01 to 100

The results showed an overall diagnostic sensitivity of 100 % and an overall diagnostic specificity of 100 % for the detection of SARS-CoV-2.

#### **Conclusion**

We established and validated a high throughput SARS-CoV-2 qRT-PCR detection workflow using samples collected in transport media, a magnetic bead RNA isolation, and a 384-well format RT-qPCR assay. Our data demonstrates a highly specific and sensitive workflow ensuring a robust analytical and clinical performance.

Internal Audit Planning PGL 2021-2022																
Norm	ISO	Section	Auditor		may-21	jun-21	jul-21	aug-21	sep-21	oct-21	nov-21	dec-21	jan-22	feb-22	mar-22	apr-22
Management	4.1	5.1.2.e	5.1.2.e													
Quality management, document management and control of records	4.2, 4.3, 4.13	5.1.2.e	5.1.2.e	System audit												
				System audit												
Service agreements / referral laboratory	4.4, 4.5, 4.6, 4.7	5.1.2.e	5.1.2.e	System audit												
				System audit												
Complaints, evaluations and audits	4.8, 4.14	5.1.2.e	5.1.2.e	System audit												
				System audit												
Non-conformities, corrective and preventive measures and continous improvement	4.9, 4.10, 4.11, 4.12	5.1.2.e	5.1.2.e	System audit												
				System audit												
Management review	4.15	5.1.2.e	5.1.2.e													
Personnel, accommodation and enviromental requirements	5.1, 5.2	5.1.2.e	5.1.2.e	System audit												
				System audit												
Laboratory equipment, reagents and consumables	5.3	Office / storage	Interne auditor	System audit												
				System audit												
Pre-analysis	5.4	Sample preparation	Interne Auditor	System audit												
				System audit												
Examination processes, quality examination results, post-examination processes	5.5, 5.6, 5.7	Molecular diagnostics	Interne auditor	System audit												
				System audit												
Reporting results and release of results	5.8, 5.9	Back office	5.1.2.e	System audit												
				System audit												
Laboratory information management	5.10	5.1.2.e	5.1.2.e	System audit												
				System audit												

planning
audit done
audit report

	<b>Internal Audit Procedure</b>	Page 1 of 3
---	---------------------------------	-------------

## PURPOSE

Assessing whether the regulations/procedures of the various sections within the laboratory meet the requirements of ISO 15189 and whether the agreements are known and effective in practice. This way continuous improvement will be established.

## DEFINITIONS AND ABBREVIATIONS

- Audit: a systematic and independent investigation to determine whether the quality activities and results thereof correspond to established regulations and whether the latter have been effectively implemented, as well as be suitable for achieving the targets.
- Auditor: a trained employee who carries out the audit
- Auditee: person to whom the audit is performed

## RESPONSIBILITY

Management:

- determining the planning of the audits
- ultimately responsible for follow-up of actions in response to the audit results

Quality officer:

- drawing up an annual audit plan
- assign audits to auditors
- control of audits in accordance with planning
- monitor progress improvement process
- archiving audit plan, reporting and action plan
- education and training internal auditors

Auditee:

- undergoes the audit

Auditor:

- conducting audit in an objective manner
- timely planning, execution and reporting of an audit
- review the required documents and method descriptions

Supervisor:

- create preconditions for conducting audits
- drawing up an action plan to resolve the findings from the report
- discussing the findings and actions with employees
- take corrective and preventive measures within the stipulated time

	<b>Internal Audit Procedure</b>	Page 2 of 3
---	---------------------------------	-------------

## PROCEDURE

1. The quality officer makes an annual audit plan with a scope of 1 year. In this plan, one or two auditors are linked to each process or department to be tested, who are sufficiently independent with regard to the processes that are audited.

The frequency of an audit depends on the nature of the process/department and the results of previous audits (1-3 years). In principle, the cycle of an internal audit must be completed within one year. Each individual element of ISO 15189 is tested annually. Every year, each individually accredited activity of ISO 15189 must be screened. A lower audit frequency is only allowed if there were no significant deviations during 2 consecutive audit periods.

2. The audit planning is submitted to the management for authorization. The auditors receive an invitation from the quality officer stating the procedure/department to which the audit relates and the period in which the audit must be performed. The supervisor must be available for questions during the audit.

3. The auditors contact the supervisor at least two weeks in advance to make agreements about the date and time for the audit. Prior to the audit, the designated auditors collect and study the necessary information such as quality manual, procedures, protocols, work instructions, and possibly. prior audit reports of the relevant project/procedure. The auditors can request this information from the quality officer. After this preparation, a checklist with points for attention is drawn up for the audit.

The auditors conduct the audit through interviews and work floor observations.

At the start of a visit, the auditors indicate the purpose, method and estimated duration of the audit. During the audit, identified shortcomings and opportunities for improvement are discussed with the auditees and summarized as much as possible.

4. During the audit, notes are made to reconstruct the information obtained afterwards. The notes will be kept until the audit report has been finalized. The auditors complete the checklist in full.

5. The auditors draw up a written report within two weeks of the audit. A template is available for this. The report contains the following parts:

- Standard form with general data such as the audited section, audit date, report date, auditors, relevant procedure and auditees.
- Short report of the conversations with general findings and (non-)conformities

The auditors send the completed checklist and the audit report to the quality officer.

6. The quality officer sends a copy of the audit report to the concerned manager. The manager checks the audit report for factual inaccuracies and reports these to the auditors within two weeks so that the report can be adjusted if necessary. The modified report is resubmitted via the quality officer.

	<b>Internal Audit Procedure</b>	Page 3 of 3
---	---------------------------------	-------------

7. If after 2 weeks the auditee does not contact the auditors, the report is considered final. The report is signed for approval by the manager.

8. The original report is archived by the quality officer. A copy of the final report goes to the supervisor and the head of the department. The supervisor determines whether and how the audit report is handled in the department.

9. The quality officer supplies the action plan by the audit findings. The supervisor determines which actions must be taken in response to the reported findings. The supervisor fills in the actions to be taken in the action plan, supplied by the quality officer, as added to the audit report. Responsible persons and deadlines are determined.

10. The action plan must be drawn up within two weeks of the final report. The management must agree to this and determine this by means of the date and initials. The manager sends a copy of the adopted plan of action to the quality officer and keeps the original.

11. The quality officer is responsible for discussing the audit report and the actions with the employees, and the supervisor is responsible for carrying out the actions in the plan of action within the specified period. The progress of the actions from the action plan are discussed in the weekly staff meeting.

12. When the actions have been carried out, they are signed off with date and initials on the action plan by the manager(s). When all actions have been initialed, the action plan is archived with the relevant audit report by the quality officer.

13. A follow-up audit can be planned and performed by the quality officer in a relatively short period of time. Depends on the goal and the activities to be tested.

#### 14. Taking Corrective Actions

Measures are taken to relieve the irregular and/or incorrect activities. The 4 O's system " Cause, Size, Solution and Operationality" is used, to provide insight into the corrected measures taken.

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### **Formation and duration of the Contract**

The Contract will run from the date of signature of this Contract until 25<sup>th</sup> April 2021 and will then automatically expire without any notice being required, unless it is extended by mutual written agreement 1 (one) month in advance the latest by period(s) of 1 (one) month. The Minister may decide on any such extension in his sole discretion and will not be under any obligation to agree to an extension.

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### **Formation and duration of the Contract**

The Contract will run from the date of signature of this Contract until July 15 2021, or if later, 21 days after the tender award has been granted and will then automatically expire without any notice being required latest by period(s) of 1 (one) month. The Minister may decide on any such extension in his sole discretion and will not be under any obligation to agree to an extension.

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[As of Regarding up to long run]

The Parties have agreed to build up to the Target capacity as described in table below and the Minister agrees to a Guaranteed volume in the following way:

5.1.2.e

Predicated percentage guaranteed volume			
# month of the runtime of the addendum	Target Capacity per day*	Percentage guaranteed	Guaranteed volume per day*
Start 26 April 2021	20.000	20%	4.000
May 2021	20.000	15%	3.000
June 2021	20.000	10%	2.000
Till 15 July 2021	20.000	10%	2.000
<b>Ultimate Capacity* : 50.000 tests per day</b>			
Maximum contracted volume (Target) = 20.000 test per day*			

)\* Measured as average per calendar week

#### 4. Wijziging bijlage A "Schedule 1: Services, Laboratories, Volumes and Price" onder punt 7 "Pricing Table":

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**Pricing Table**

Pricing table (excl. VAT)		
NUMBER OF REPORTED TESTS PER DAY*		PRICE PER REPORTED TEST
Minimum	Maximum	
0	No maximum	€ 50,00

Contractor will be paid 50 euro, per reported tests. [for the sake of clarity: For – test: to be deleted]

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**Artikel 5.2 luidt tot 25 april 2021**

5.2 During the utilization of the Laboratory the initial quality requirements as referred to in article 5.1 have to be maintained or improved. As soon as Contractor foresees that the Laboratory cannot meet these minimum requirements anymore Contractor shall immediately inform the Minister. In that case the Minister is entitled to terminate the contract with immediate effect without any prior notice being required, by registered letter and without being due any compensation for lost income or other damage by the Contractor.

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5.2 During the utilization of the Laboratory the initial quality requirements as referred to in article 5.1 have to be maintained or improved. As soon as Contractor foresees that the Laboratory cannot meet these minimum requirements anymore Contractor shall immediately inform the Minister. In that case the Minister is entitled to withdraw or suspend its approval for that particular Laboratory resulting in the fact that Contractor is not any longer allowed to utilize that Laboratory for the Diagnostic test. In case no Laboratories are allowed to provide the Services, the Contract may be terminated with immediate effect by (electronic) letter without any compensation for lost of income or other damage by the Contractor

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5.3 If the Minister has suspended its approval as mentioned in the previous the Minister will not allocate any Sample to the Laboratory. If the Minister terminates the suspension after resolving the quality problem underlying the suspension, the allocation of the Samples will resume two weeks after the suspension has ended.

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5.5 If during such regular assessment or audit it appears that a Laboratory is no longer meeting the quality requirements as referred to in 5.1, the Minister is entitled to withdraw the approval for that particular Laboratory resulting in the fact that Contractor is not any longer allowed to utilize that Laboratory for the Diagnostic tests.

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- Cooperation in relation to specific research commissioned by the Minister

Bijlage B "Schedule 2: Operational agreements" onder punt 1 "Organization" Toevoeging extra subbullit bij bullit 5 "contractor is responsible for":

- o collection of positive Samples and transport of these Samples to a location for research

Bijlage B "Schedule 2: Operational agreements" toevoeging van punt 7 "Collection of Samples and Transport to location for research":

- In case of research commissioned by the Minister, the Contractor will select and pack the Samples and transport them to a location designated by the Minister.

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- Negatively and undetermined tested Samples will be appropriately destroyed by Contractor 3 (three) days after publication of the test result of the CoronIT.

5.1.2.e

- Positively tested Samples will be appropriate and traceable stored by Contractor for 2 (two) months after publication of the test result of the Samples in CoronIT.

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- Negatively and undetermined tested Samples must be appropriately destroyed by Contractor 3 (three) days after publication of the test result of the Samples in CoronIT.
- Positively tested Samples must be appropriate and traceable stored by Contractor for 2 (two) months after publication of the test result of the Samples in CoronIT, after which the Samples must be appropriately destroyed by Contractor

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**Timely information (< 2 hours after arrival of Samples in Lab) of Sampling hubs/LCDK about disturbances leading to delays in publication of test results**

- bonus 5 points in case of less than two – according the Minister well-founded – notifications by the Sampling hubs per week. These notifications concern alleged deviation in the performance of the Agreement and must be communicated in writing ( email) to the Contractor and Minister ( via the Testing Service Ops [ops@lcdk.nl](mailto:ops@lcdk.nl)).

**Timely information (> 2 hours in advance) of Sampling hubs about delays in supplies or Sample collection**

- bonus 5 points in case of less than two – according the Minister well-founded – notifications by the Sampling hubs per week. These notifications concern alleged deviation in the performance of the Agreement and must be communicated in writing ( email) to the Contractor and Minister ( via the Testing Service Ops [ops@lcdk.nl](mailto:ops@lcdk.nl)).

**No interruption of testing activities caused by shortage in agreed safety stocks**

- bonus 5 points if no disturbance of sample taking/testing process

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**Contract number: .....**

**PUBLIC SERVICE CONTRACT ("CONTRACT") FOR THE  
PROVISION OF LOGISTIC AND COVID-19 TESTING SERVICES**

**BETWEEN**

**THE STATE OF THE NETHERLANDS  
THE MINISTRY OF HEALTH, WELFARE AND SPORTS**

**AND**

**PGL Finance BV**

**The undersigned:**

1. The State of the Netherlands, which has its seat in The Hague, represented by the Minister of Health, Welfare and Sport, Directorate-General for Public Health, legally represented in this matter by [REDACTED] Head of Testing Services hereinafter referred to as the "Minister",  
  
and
2. The private limited company PGL Finance BV, which has its registered office at Bieslandseweg 2, 2645BM Delfgauw, the Netherlands, registered with the Chamber of Commerce under number 27352865, legally represented in this matter by [REDACTED], hereinafter referred to as "Contractor",

Each party also referred to as "Party", and both parties together as "Parties",


**WHEREAS:**

- a) the Minister has a facilitating role in infectious disease control on the basis of the Public Health Act (Wpg). Pursuant to Article 3(1) of the Wpg, the Minister has the task to promote the quality and efficiency of public health care. In addition, he is responsible for the maintenance and improvement of the so-called support structure. On the basis of article 7(1) of the Wpg, the Minister is also in charge of combating an infectious disease from Group A, such as Covid-19;
- b) The Minister therefore has engaged a number of parties to provide logistic- and diagnostic services in order to urgently provide temporarily additionally needed laboratory capacity to test samples at (community) sampling hubs across the Netherlands for the COVID-19 virus;
- c) Contractor is an innovative diagnostic knowledge center, and offers a wide range of diagnostics and laboratory services in the field of clinical chemistry, immunology and microbiology, including COVID-19;

- d) Contractor is prepared to do everything reasonable to urgently provide temporary additional laboratory capacity in an efficient and cost-effective way;
- e) As it is strictly necessary for reasons of extreme urgency the Minister, and in accordance with article 2.32, section 1, paragraph c of the Dutch Public Procurement Act 2012 ('Dutch PPA' or 'Aanbestedingswet 2012'), the Minister awards this Contract to Contractor. The Minister has awarded this Contract by a negotiated procedure without publication, as it is strictly necessary where, for reasons of extreme urgency brought about by events related to the COVID 19 pandemic which were unforeseeable by the Minister, the time limits for the open or restricted procedures or competitive procedures with negotiation cannot be complied with (also see EC 2020/C 108 I/01);
- f) the circumstances invoked to justify extreme urgency shall not in any event be attributable to the Minister. The events were unforeseeable by the Minister, because the estimated need for test capacity was reached months earlier than expected;
- g) The extreme urgency makes compliance with general deadlines impossible because the time that is needed for compliance will result in a growing shortage of test capacity. There is a causal link between the unforeseen event and the extreme urgency because the gap between the official estimation of the needed test capacity and the current rise in demand for test capacity is growing every day. This procedure is only used by the Minister in order to cover the gap until more stable solutions can be found;
- h) Parties have discussed the main conditions for the provision of test volumes by Contractor;
- i) Parties wish to lay down their agreements in this Contract. With the signature of this Contract all previously made arrangements are no longer valid and are replaced with the arrangements in this Contract.

## **AGREE AS FOLLOWS:**

### **Definitions**

The following terms are written with initial capitals in this Contract and defined as follows:


Contract	this contract as defined in article 1.7
CoronIT	CoronIT is the name of registration software
Diagnostic test	detection of SARS-CoV-2 RNA by RT-PCR
First production date:	the date that all processes for providing the Services are in place (week n)
Guaranteed volume	the minimum number of Samples as percentage of the Target capacity, guaranteed by the Contractor
Laboratories	laboratories meeting the quality requirements of this Contract, performing Covid-19 Diagnostics tests
Operational agreements	the agreements as laid down in Schedule 2
RIVM	Rijksinstituut voor Volksgezondheid en Milieu, meaning the National Institute for Public Health and the Environment
Services	the services as defined in Schedule 1, "Services"
Samples	sample material of visitors of the Sample hubs
Sampling hubs	centres where Samples are taken
Sampling kits	kits containing all materials needed for the taking of Samples
Target capacity	the volume the Minister and Contractor agreed upon as reserved volume for use by the Minister
Ultimate capacity	the ultimate volume Contractor can offer to the Minister in the long run
Volumes	amount of Services as defined in Schedule 1, 'Volumes of Services'

## **1. Object of the Contract**

- 1.1 Parties agree on the provision by Contractor of temporarily additionally needed diagnostic and logistic services ("Services") in certain volumes ("Volumes") for the duration of this Contract including the supply to Sampling hubs of Sampling Kits and labels, tube racks and containers.
- 1.2 Contractor guarantees to provide to the Minister the Target capacity and the Services under the conditions as set out in this Contract and in specific Schedule 1 (Services, Laboratory, Volumes and Price), attached to this Contract as Annex A.
- 1.3 The Minister is free to order lesser Volumes than the Target capacity and/or only part of the Services as described in Schedule 1.
- 1.4 The Minister guarantees payment to Contractor of 30% of the Target capacity (the Guaranteed Volume), by way of compensation to the Contractor for the reservation of the Target capacity per day as provided for in this article, on the condition that Contractor performs its other obligations in conformity with this Contract.
- 1.5 On request of the Minister Contractor is prepared to explore an increase of the Volumes to the Ultimate capacity as mentioned in Schedule 1 Services, Laboratory, Volumes and Price) attached to this Contract as Annex A. If Parties agree to increase the Target capacity towards the Ultimate capacity, the new Volumes will be the new Target capacity.
- 1.6 Parties will comply with the operational organisation arrangements as mentioned in Schedule 2 Operational Arrangements attached to this Contract as Annex B and the Key Performance indicators as mentioned in Schedule 3 attached to this Contract as Annex C.
- 1.7 The following documents are an integral part of this Contract. In the event of inconsistencies, a higher ranked document takes precedence over a lower ranked document:
  1. this document;
  2. Annex A: Schedule 1 – Services, Laboratory, Volumes and Price;
  3. Annex B: Schedule 2 – Operational agreements;
  4. Annex C: Schedule 3 – Key Performance Indicators;
  5. Annex D: General Government Purchasing Conditions 2018 (ARIV 2018);
  6. Annex E: Financial Brochure e-invoicing government ("Financiele Bijsluiter").

## **2. Formation and duration of the Contract**

- 2.1 The Contract will run from the date of signature of this Contract until 25<sup>th</sup> April 2021 and will then automatically expire without any notice being required, unless it is extended by mutual written agreement 1 (one) month in advance the latest by period(s) of 1 (one) month. The Minister may decide on any such extension in his sole discretion and will not be under any obligation to agree to an extension.
- 2.2. The Guarantee as mentioned in article 1.6 will only be invoked after the mutually agreed First production date.

## **3. Price and invoicing**

- 3.1 The agreed price to be paid to Contractor for all Services to be performed by Contractor ("Price") and the Volumes to be met for the Services are laid down in Schedule 1.
- 3.2 The Price is exclusive of VAT.
- 3.3 The Price includes all Services as mentioned in Schedule 1.
- 3.4 Contractor will invoice every 2 (two) weeks for the attention of the Minister electronically in the manner prescribed in annex D ("Financiële Bijsluiter").
- 3.5 Contractor will send a copy of the invoice together with a print from CoronIT showing the Volumes tested during the invoiced period, together as attached pdf-files via 1 (one) e-mail message to OGZ-Covid@rivm.nl.
- 3.6 Contractor is allowed to invoice only Volumes actually tested (Real volume) covered by the invoice.
- 3.7 In case the Real volume per day is less than the Guaranteed volume per day, as defined in this Contract and Annexes, Contractor is allowed to invoice the Price for the Real volume as well as the difference between the Guaranteed volume and the Real volume. In this case Contractor will inform the Minister of this situation via e-mail message to secretariaat@lcdk.nl before sending in the invoice as mentioned in article 3.5
- 3.8 Invoices will be paid with a maximum payment term of 30 (thirty) days from the date of receipt of the invoice.
- 3.9 The Minister will check the invoices on a regular basis. Incorrect invoices will be discussed with Contractor and in mutual consent be corrected immediately.

#### **4. Contacts**

##### **4.1 The Minister's contact is:**

[REDACTED]

Head of Testing Services ("Dienst Testen") – Ministry of Health, The Netherlands

Via e-mail: [secretariaat@lcdk.nl](mailto:secretariaat@lcdk.nl)

##### **4.2 Contractor's contact is:**

[REDACTED]

CEO

E-mail: [REDACTED]@puregeneticlifestyle.com

#### **5. Quality**

5.1 Contractor guarantees to the Minister that before utilization of a Laboratory for Diagnostic tests, this Laboratory meets the Quality Requirements for Laboratories as mentioned in Schedule 1.

5.2 During the utilization of the Laboratory the initial quality requirements as referred to in article 5.1 have to be maintained or improved. As soon Contractor foresees that the Laboratory cannot meet these minimum requirements anymore Contractor shall immediately inform the Minister. In that case the Minister is entitled to terminate the contract with immediate effect without any prior notice being required, by registered letter and without being due any compensation for lost income or other damage by the Contractor.

5.3 During utilization of the Laboratories the quality of the Laboratories will be regularly assessed and/or audited on behalf of the Minister.

5.4 If during such regular assessment or audit it appears that a Laboratory is no longer meeting the quality requirements as referred to in 5.1, the Minister is entitled to terminate the contract with immediate effect without any prior notice being required, by registered letter and without being due any compensation for lost income or other damage by the Contractor.

5.5 Contractor guarantees that the quality of Services is in accordance with prevailing laws and regulations, and practices, which would be considered as good practice of a skilled and experienced organization of a similar scale and professional standing, engaged in the provision of comparable and related services.

[REDACTED]

[REDACTED]

**6. Safety stock**

- 6.1 Contractor guarantees at the Laboratories a safety stock of reagents and disposables for the assigned Volumes of Diagnostic tests for at least 10 (ten) days based on the average Target volume.
- 6.2 Contractor guarantees at a central place a safety stock of reagents and disposables for the assigned Volumes of Diagnostic tests for at least 30 (thirty) days of the average Target capacity.
- 6.3 Contractor guarantees at the Sample hubs a safety stock of Sampling kits for the assigned Volumes for 7 (seven) days based on the average Target volume.
- 6.4 Contractor guarantees at a central place a safety stock of Sampling kits and containers for the assigned Volumes for 30 (thirty) days of the average Target volume.

**7. Throughput time, timely information**

- 7.1 Contractor guarantees a maximum throughput time of the Service of 24 (twenty-four) hours, and an average throughput time of the Service of 22 (twenty-two) hours.
- 7.2 Contractor shall timely inform Sampling hubs and the Minister about any disturbances leading to delays in publication of test results, but in any case, within 2 (two) hours of arrival of Samples in the Laboratory and shall endeavour to avoid structural complaints from Sampling hubs.
- 7.3 Contractor shall timely inform Sampling hubs about any delays in supplies or Sample collection, but in any case, at least 2 (two) hours in advance.

**8. Use of third Parties**

- 8.1 In performing the Services of the Contract, Contractor may make use of a third party or a party in which Contractor has a minority interest only if this party has accepted the same quality requirements as laid down in this Contract and Annexes in writing and after Contractor has obtained the prior Minister's consent in writing.
- 8.2 Notwithstanding the Minister's consent, Contractor remains fully responsible for discharging the obligations imposed under this Contract.

## **9. Confidentiality and Media**

- 9.1 Unless required by law or necessary in light of the objectives of this Contract, the Parties shall keep the terms and the nature of all discussions between them relating to this Contract and relating to affairs pursuant to the Contract confidential and shall not publish, disclose or provide to third parties any information from this Contract and Annexes. Parties will ensure that their respective employees, agents and sub-contractors fully comply herewith.
- 9.2 Nothing in this Agreement shall prevent the Minister from any disclosure required under the Dutch Act Dutch Act 'Wet Openbaarheid van Bestuur' (*Freedom of information act*) or to comply with the parliamentary rights on information. Such disclosure shall not constitute a default under this Contract.
- 9.3 The fact that this Agreement has been concluded as well as the essence of this Contract may be disclosed by the Minister to the Staten-Generaal (parliament).
- 9.4 The Parties shall consult each other timely in case they are approached by Media or consider to react to publications or news in the Media, or to approach Media, in relation to the subject of this Contract. In their contacts with Media, the Parties will respect each other's reasonable interests.

## **10. Data protection**

- 10.1 In so far as the Contractor, as a processor within the meaning of the General Data Protection Regulation, processes personal data for the Sampling hub(s), the Contractor guarantees the application of appropriate technical and organisational measures, so that processing meets the requirements of the General Data Protection Regulation and the protection of the data subjects is ensured.
- 10.2 Parties agree that personal data will be processed in accordance with the agreements as made between Contractor and the Sampling hub(s) in relation to the connection and use of CoronIT.

## **11. Limitation of liability**

- 11.1 The total aggregate liability of each Party under this Contract (for breach of contract, tort or otherwise) shall not exceed EUR 2.500.000 (two and a half million EUR), in no event shall any party be liable to the other party for any indirect or consequential damages (*gevolgschade*) out of or in connection with this Contract. Parties do not exclude their liability for gross negligence or wilful misconduct.

**12. Insurance**

- 12.1 The Contractor has taken out and will retain insurance cover for the total aggregate liability as mentioned in article 11.1.
- 12.2 At the Ministers request, the Contractor will immediately present either the original or certified copies of the policies and proof of the payment of premium for the insurances referred to in paragraph 1 or a statement from the insurer to the effect that these policies exist and the premium has been paid. The insurance premiums owed by the Contractor are deemed to be included in the agreed prices and fees.

**13. Cancellation and termination**

- 13.1 Without prejudice to the other provisions of this Contract, either Party may cancel the Contract in full or in part with immediate effect by registered letter, without recourse to the courts, if the other Party is in default or compliance is permanently or temporarily impossible.
- 13.2 In case Contractor is not able to meet the Target Capacity and/or the minimum agreed total score as defined in Annex C / Schedule 3 during a period of 3 (three) consecutive weeks or 3 (three) times in a total time period of 13 (thirteen) weeks, the Minister is entitled to cancel the Contract in full with immediate effect by registered letter, without recourse to the courts and without any compensation of Contractors loss of profit incurred.
- 13.3 The Minister may cancel the Contract forthwith out of court by registered mail, without being required to send any demand or notice of default, if the Contractor does not meet the Target Capacity and/or the minimum agreed total score as defined in Annex C / Schedule 3 as referred to in 13.2, if the Contractor applies for or is granted a provisional or definitive suspension of payment, files for bankruptcy, is declared bankrupt, if its business is wound up, if it ceases trading, if a substantial proportion of its assets are seized, if it undergoes a merger or a division or is dissolved which reasonably effects the compliance to the Contract, or if it should reasonably be deemed on any other grounds to be no longer capable of discharging its obligations under the Contract.


#### **14. Force Majeure**

- 14.1 In the event either Party is unable to perform its obligations under the terms of this Contract, because of causes reasonably beyond its control (including but not limited to government restrictions or lack of access to equipment, unforeseen at the moment of signing this Contract) or damage reasonably beyond its control, such Party shall not be liable against the other Party for any damages resulting from such failure to perform or otherwise from such causes.
- 14.2 The term 'Force Majeure' is in any event understood not to include: allocation of testing capacity by Contractor for the benefit of other clients, staff shortages, strikes, staff illness, shortages of materials (with exception of an acknowledged worldwide shortage), transport problems, late delivery or the unsuitability of items required for the performance of the Services, liquidity or solvency problems on the part of the Contractor or failures on the part of parties engaged by the Contractor.
- 14.3 If a Force Majeure event prevents the Contractor meeting its Target capacity as defined in Schedule 2, the Contractor will give notice thereof within 48 (forty-eight) hours of the occurrence such event, failing which it will not be entitled to rely on paragraph 1 of this Article 14. In this case Contractor will not be relieved from providing its Services in respect of the volumes it can reasonably supply despite such Force Majeure.
- 14.4 A notice of Force Majeure shall be well motivated and, in the event of foreign Government restrictions, be accompanied by an authorized statement from the relevant government evidencing such fact.

#### **15. Other Terms and Conditions**

- 15.1 This Contract is subject solely to the General Government Purchasing Conditions 2018 (ARIV 2018), in so far as the Contract does not depart from them. Any general and special terms and conditions drawn up by the Contractor do not apply. These General Terms are attached to this Contract in Annex D.
- 15.2 Parties will discuss in good faith and will mutually decide to any amendments and/or additional conditions to this Contract and/or Annexes in case unforeseen circumstances affect the performance of the agreed contractual obligations. Amendments to this Contract and/or Annexes shall not be binding unless they are agreed upon and made in writing and signed by both parties.

15.3 Contractor will provide all relevant information and after timely pre-announcement grant access to utilized locations on the Minister's reasonable request for verification purposes of (the quality of) the utilized Locations and the agreed Services.

**16. Governing law and jurisdiction**

16.1 This Contract is governed by the laws of the Netherlands.

16.2 Any dispute that arises between the parties in relation to the explanation or application of this Contract, or any agreement resulting from this Contract, shall be settled by the competent court of the Hague, the Netherlands. Before starting court litigation, Parties agree to jointly attempt to come to a solution through mediation unless one of the Parties is of the opinion that this is not effective under the circumstances. In the case of mediation, Parties will jointly appoint an MfN-registered mediator.

**17. Final provisions**

17.1 This Contract sets out the entire agreement between the Parties. Any derogations from, amendments of or additions to this Contract and/or Annexes, including changes in the scope of the Services are binding only if expressly agreed in writing by the Parties.

17.2 Any written or oral agreements previously made by the Parties about the Services that are the object of this Contract are nullified by the signature of this Contract.

Done on the date stated below and signed by Parties.

The Hague, the Netherlands,  
For the State, Minister of Health

31/12/2020

Head of Testing Services

Delfgauw, the Netherlands  
For PGL Finance BV

31/12/2020

CEO

## **Schedule 1 Services, Laboratories, Volumes and Price**

### **1. Services**

Contractor will provide the following services ("Services"):

- Set-up and distribution of an appropriate introduction of relevant persons to be contacted and an instruction letter to the Sampling hubs;
- Self-sufficient procurement of all materials needed for the execution of the Services;
- Supply of Sampling kits and containers and/or tube racks to Sampling hubs;
- Collection and transportation of Samples from Sampling hubs to Laboratories;
- Diagnostic test of Samples;
- Registration of the results of the Diagnostic test in CoronIT;
- Appropriate and traceable storage of positive tested Samples;
- Destruction of tested Samples;
- Dedicated, appropriate and where relevant pro-active Support to the Services.

For operational details of the Services see Schedule 2 'Operational details of Services'.

For key performance indicators of the Services see Schedule 3 'Key Performance Indicators'.

### **2. Sampling hubs**

Contractor will provide Sampling kits and containers to and will pick up Samples from the Sampling hubs in the Netherlands as assigned on behalf of the Minister (via "Dienst Testen").

### **3. Quality requirements laboratories**

Before utilization by the Contractor Laboratories first have to meet the following initial quality requirements:

- The Laboratory is capable in processing potential SAR-CoV-2 positive specimens at BSL-2 (ML-2) level;
- The Laboratory has at least one dedicated Medical Microbiologist medically responsible for the quality of the Diagnostic services and available for Sampling hubs for answering questions about the test results;
- The laboratory has a valid certificate for ISO15189, ISO17025 with RT-PCR in scope provided by a conform European standard accredited certification institute.

After a Laboratory has met all initial quality requirements the Laboratory has to be validated by the RIVM in the following way:

- Successful execution of a diagnostic sensitivity test of samples provided by the RIVM;
- Successful execution of a diagnostic specificity test of samples provided by the RIVM;
- Both tests are only successful if at least the core/reference Samples are diagnosed.

After the validation the Implementation phase starts in the following way:

- A between the Parties agreed implementation plan;
- A site visit;
- An audit;
- A between the Parties agreed plan how to build up the Target volume.

After the successful fulfilment of all quality requirements Contractor accepts for the purpose of monitoring the quality of the Laboratory during utilization regular proficiency tests.

#### 4. Approved Laboratories

At the moment of concluding this Contract the following Laboratories are used and approved:

Location	% of shares	Initial requirements RIVM Approval	Max test capacity per laboratory	Remarks
Novogenia, Eugendorf, Austria	15%	Yes	50.000	
Novogenia, Vlaardingen	15%	N/A	N/A	Labhub

If Contractor after the signature of this Contract wishes to use another Laboratory Contractor will see to it that before utilization the Laboratory in accordance with article 5 of this Contract meets the Quality requirements and that the Minister has given approval before utilization.

#### 5. Allowed third parties

Party	% of shares	Activities	Remarks
n/a			

If Contractor after the signature of this Contract wishes to use another third party, Contractor will inform the Minister before utilization in accordance with this Contract.

#### 6. Volumes of Services (tests per day)

Regarding the Volumes of Services four kinds of volumes are identified:

1. Target capacity      the Volume the Minister and Contractor agreed upon as reserved Volume for use by the Minister.
2. Real volume          the total actual tested Samples for a particular day
3. Guaranteed volume   the minimum guaranteed number of tested Samples as a percentage of the Target capacity.
4. Ultimate capacity    the ultimate Volume Contractor can offer to the Minister in the long run.

The Parties have agreed to build up to the Target capacity as described in table below and the Minister agrees to a Guaranteed volume (30% of Target capacity) in the following way:

Target capacity build up and Garanteed volume		
Milestones in calender weeks	Target Capacity per day	Guaranteed volume per day*
Start contract Week 1	500	N/A
Week 2	3.000	900
Week 3	10.000	3.000
Week 4	18.000	5.400
Week 5	20.000	6.000
Week 6	20.000	6.000
Week 7	20.000	6.000
Week 8	20.000	6.000
Week 9	20.000	6.000
Week 10	20.000	6.000
Week 11	20.000	6.000
Week 12	20.000	6.000
Week 13	20.000	6.000
Week 14	20.000	6.000
Week 15	20.000	6.000
Week 16	20.000	6.000
<b>Ultimate Capacity: 50.000 tests per day</b> <b>Maximum contracted volume (Target) = 20.000</b>		

)\* Measured as average per calendar week

## 7. Pricing Table

<b>Pricing table (excl. VAT) based on bandwidths of Volumes*</b>			
Till 10.000 tests per day: 65 Euro, after 10K each extra 5.000 tests price decreases till 50 Euro			
NUMBER OF REPORTED TESTS PER DAY		PRICE PER REPORTED TEST	
Minimum	Maximum		
0	9.999	€	65,00
10.000	14.999	€	62,00
15.000	19.999	€	57,50
20.000	24.999	€	53,00
>25.000	>25000	€	50,00

)\* Measured as average per calendar week

For the first 10.000 reported tests Contractor will be paid 65 euro, for the next 5.000 reported tests 62 euro, for the next 5.000 reported tests 57,50 euro, for the next 5.000 reported tests 53 euro and the number of reported tests thereafter 50 euro per reported test.

## 8. Safety stock

### *At Laboratories*

- Contractor guarantees at the Laboratories a safety stock of reagents and disposables for the assigned Volumes of Diagnostic tests for at least 10 days based on the average Target volume.
- Contractor guarantees at a central place a safety stock of reagents and disposables for the assigned Volumes of Diagnostic tests for at least 30 days of the average Target capacity.

### *At Sample hubs*

- Contractor guarantees at the Sample hubs a safety stock of Sampling kits for the assigned Volumes for 7 days based on the average Target volume.
- Contractor guarantees at a central place a safety stock of Sampling kits and containers for the assigned Volumes for 30 days of the average Target volume.


## 9. Example

This paragraph contains **an example** of how a situation in practice may work out on a specific day. **For instance** Parties agree to **106.000 test** per day as Ultimate capacity and a maximum Target capacity of 80.000 tests per day.

### EXAMPLE

Milestones	Target Capacity per day	Garanteed volume per day	Real Volume per day
Start contract Week n			
Week n+1	8.000	2.400	6.392
Week n+2	20.000	6.000	14.348
Week n+3	40.000	12.000	18.356
Week n+4	40.000	12.000	16.964
Week n+5	60.000	18.000	29.408
Week n+6	60.000	18.000	43.572
Week n+7	80.000	24.000	48.980
Week n+8	80.000	24.000	55.568
Week n+9	80.000	24.000	44.100
Week n+10	80.000	24.000	53.072
Week n+11	80.000	24.000	68.956
Week n+12	80.000	24.000	54.752
Ultimate Capacity = 106.000 per day; Maximum contracted volume (Target) = 80.000			

### EXAMPLE

## **Schedule 2 Operational agreements**

### **1. Organization**

- The Minister is overall responsible for the execution and the management of the Contract and the payment of the Contractor.
- RIVM is responsible for the validation and approval of the quality prior to the utilization of the Laboratories of Contractor.
- The Minister (via "Dienst Testen") is responsible for the assignment of Sampling hubs and Volumes to Contractor and the monitoring of the quality of Contractors Laboratories throughout utilization.
- Sampling hubs are responsible for stock keeping of the by Contractor provided Sampling Kits, sample taking and registration of the samples and citizens in CoronIT.
- Contractor is responsible for:
  - the assignment of dedicated and knowledgeable teams of medical microbiologists, logistical experts and project management for troubleshooting, (pro-active) analysis and continuous improvement of the quality of the Services and answering questions about these Services;
  - the supply of Sampling Kits, labels, tube racks and containers;
  - collection and transportation of the samples;
  - the testing of the samples and the registration of the testing results in CoronIT;
  - the traceable storage and/or the destruction of Samples;
  - monthly reports about complaints, responses and solutions to the Minister (via "Dienst Testen").

### **2. Supply of Sampling kits, containers and/or tube racks**

- Contractor provides Sampling kits with the following materials:
  - A pre-labelled collection tube with cap filled with buffer (Lyses, UTM or GLy) or an unlabelled collection tube with cap filled with buffer (Lyses, UTM or GLy) and separate labels, to be chosen by the Sampling hub;
  - The label is resistant to methyl alcohol (70%), storage at -80°C, cannot be erased by thermal heat and is provided with a unique barcode, containing information of the Contractor or Laboratory and a tube number, which can be scanned with a maximum failure rate of 0,01% or 100 ppm;

- The buffer is compatible with the Diagnostic test and not dangerous during use by the Sampling Hub staff;
- A CE-marked and patient-friendly nasopharyngeal- and throat swab;
- A seal bag with absorption filter for a safe transportation of the Sample.
- Contractor supplies beside the Sampling kits also appropriate disinfected containers and/or tube racks in an amount so that the sample taking process never must wait and, in a quality, so that collection and transportation of the Samples can be securely executed.
- Contractor supplies the Sampling kits, containers and/or tube racks in consultation with the Sampling hubs minimum once a week, during opening times and free of charge.
- Contractor is responsible for an appropriate registration of the supplies.
- Contractor provides for its risk and costs a safety stock of Sampling kits at Sampling hubs for a minimum of 7 days of the assigned Volumes.  
Contractor provides for its risk and costs a safety stock of Sampling kits at Contractor or a central location in the Netherlands for a minimum of 30 days of the assigned number of tests.

Communication	Service window	Telephone Number	Response time
Logistic issues	8:00 – 22:00; 7 days/week	085-1058499	15 minutes

### 3. Storage of Sampling kits

- Sampling hub will store the Sampling kits on its own costs at a temperature between 15 and 25°C.

### 4. Sampling

- The Sampling Hub registers the used collection tube and the information of the sampled citizen in CoronIT
- The Sampling Hub collect samples in accordance with the instructions of Contractor and RIVM ([www.lci.rivm.nl/richtlijnen/COVID-19](http://www.lci.rivm.nl/richtlijnen/COVID-19)).
- The Sampling Hub is responsible for the careful packaging of the samples in the sealing bags or tube racks and the containers.

## 5. Collection and transportation of the Samples to the Laboratories

- Contractor collects in consultation with the Sampling hubs during opening times the Samples minimum twice a day, the last moment of collection just before closure time of the Sample hub.
- Contractor keeps a registration of the collected Samples.
- Contractor consolidates and/or transports the Samples to the Laboratory in accordance to UN3373.
- Contractor chooses a collection frequency and way of consolidation and/or transportation so that the total agreed maximum throughput time will be as short as reasonably possible, but at any time within the agreed maximum throughput time.
- Contractor keeps a registration so that collected and transported Samples can be traced.

Communication	Service window	Telephone Number	Response time
Logistic issues	8:00 – 22:00; 7 days/week	085-1058499	15 minutes

## 6. Diagnostic test of Samples

- Contractor assures that the connection with CoronIT is up and running before starting the testing of the Samples.
- Contractor keeps a registration of the received Samples. Contractor will check the received tubes and will scan the tubes for registration in CoronIT.
- Missing tubes and/or missing or unreadable labels will be immediately (within 2 hours) reported to the Sampling hubs where the batch Samples came from.
- If Contractor foresees that the agreed throughput time cannot be met, Contractor will immediately inform the Sampling hub and timely (< 2 hours after reception of the Samples) inform the Sampling hub and the Minister ("Dienst Testen") about the causes, the corrective actions and when the test results than can be expected;
- Contractor tests all received Samples.
- Contractor reports the results of these tests in CoronIT in an accurate way.
- Contractor will do its utmost to keep the number of 'undetermined' test results as low as reasonable possible, but at any time < 1%.
- Contractor will do its utmost to realize the shortest possible throughput time between the moment of Sampling at the Sampling hubs and the publication of the test results of the Samples in CoronIT, but anyway within the agreed maximum throughput time.

Communication	Service window	Telephone Number	Response time
Logistic issues	8:00 – 22:00; 7 days/week	085-1058499	15 minutes
Microbiological issues	8:00 – 22:00; 7 days/week	085-1058499	15 minutes

Missing tubes and/or missing or unreadable labels reporting	a.s.a.p.
---	----------

Any foreseen delay of publication of the test results reporting	a.s.a.p.
---	----------

Contractor reports the test results of Samples in CoronIT	a.s.a.p.
---	----------

The between Parties agreed average throughput time per week is:	
For Laboratories situated in the Netherlands;	< 16 hours
For Laboratories situated outside the Netherlands.	< 22 hours

The between Parties agreed maximum throughput time per week is:	
For Laboratories situated in the Netherlands;	24 hours
For Laboratories situated outside the Netherlands.	24 hours

## 7. Storage and destruction of Samples

- Negatively and undetermined tested Samples will be appropriately destroyed by Contractor 3 (three) days after publication of the test result of the Samples in CoronIT.
- Positively tested Samples will be appropriate and traceable stored by Contractor for 2 (two) months after publication of the test result of the Samples in CoronIT.

Total days of storage of negative or undetermined Samples	3 days
Total days of storage positive Samples	2 months

**Schedule 3. Key Performance Indicators**

Key Performance Indicators	Score	Score																																																																						
<b>Reliability of average throughput time per week *</b> <table><tr><td>For Laboratories situated in the Netherlands;</td><td>&lt; 16 hours</td></tr><tr><td>For Laboratories situated outside the Netherlands.</td><td>&lt; 22 hours</td></tr></table> - 50 points for agreed throughput time - minus x points (progressive) per hour > agreed throughput time - bonus x points (progressive) per hour < agreed throughput time <table><tr><th colspan="3">AVERAGE THROUGHPUT TIME IN THE NETHERLANDS</th></tr><tr><td>Average TPT &lt;= 12,5 uur</td><td>+20</td><td>70 points</td></tr><tr><td>12,5u &lt; Average TPT &lt;= 13,5u</td><td>+10</td><td>60 points</td></tr><tr><td>13,5u &lt; Average TPT &lt;= 14,5u</td><td>+5</td><td>55 points</td></tr><tr><td>14,5u &lt; Average TPT &lt;= 15,5u</td><td>+2</td><td>52 points</td></tr><tr><td><b>15,5u &lt; Average TPT &lt;= 16,5u</b></td><td><b>0</b></td><td><b>50 points</b></td></tr><tr><td>16,5u &lt; Average TPT &lt;= 17,5u</td><td>-2</td><td>48 points</td></tr><tr><td>17,5u &lt; Average TPT &lt;= 18,5u</td><td>-5</td><td>45 points</td></tr><tr><td>18,5u &lt; Average TPT &lt;= 19,5u</td><td>-10</td><td>40 points</td></tr><tr><td>20,5u &lt; Average TPT &lt;= 21,5u</td><td>-15</td><td>35 points</td></tr><tr><td>21,5u &lt; Average TPT &lt;= 22,5u</td><td>-20</td><td>30 points</td></tr><tr><td>22,5u &lt; Average TPT &lt;= 23,5u</td><td>-30</td><td>20 points</td></tr><tr><td>Average TPT &gt; 23,5u</td><td>-40</td><td>10 points</td></tr></table> <table><tr><th colspan="3">AVERAGE THROUGHPUT TIME OUTSIDE THE NETHERLANDS</th></tr><tr><td>Average TPT &lt;= 18,5 uur</td><td>+20</td><td>70 points</td></tr><tr><td>18,5u &lt; Average TPT &lt;= 19,5u</td><td>+10</td><td>60 points</td></tr><tr><td>19,5u &lt; Average TPT &lt;= 20,5u</td><td>+5</td><td>55 points</td></tr><tr><td>20,5u &lt; Average TPT &lt;= 21,5u</td><td>+2</td><td>52 points</td></tr><tr><td><b>21,5u &lt; Average TPT &lt;= 22,5u</b></td><td><b>0</b></td><td><b>50 points</b></td></tr><tr><td>22,5u &lt; Average TPT &lt;= 23,5u</td><td>-10</td><td>40 points</td></tr><tr><td>23,5u &lt; Average TPT &lt;= 24,5u</td><td>-25</td><td>25 points</td></tr><tr><td>Average TPT &gt; 24,5u</td><td>-40</td><td>10 points</td></tr></table>	For Laboratories situated in the Netherlands;	< 16 hours	For Laboratories situated outside the Netherlands.	< 22 hours	AVERAGE THROUGHPUT TIME IN THE NETHERLANDS			Average TPT <= 12,5 uur	+20	70 points	12,5u < Average TPT <= 13,5u	+10	60 points	13,5u < Average TPT <= 14,5u	+5	55 points	14,5u < Average TPT <= 15,5u	+2	52 points	<b>15,5u &lt; Average TPT &lt;= 16,5u</b>	<b>0</b>	<b>50 points</b>	16,5u < Average TPT <= 17,5u	-2	48 points	17,5u < Average TPT <= 18,5u	-5	45 points	18,5u < Average TPT <= 19,5u	-10	40 points	20,5u < Average TPT <= 21,5u	-15	35 points	21,5u < Average TPT <= 22,5u	-20	30 points	22,5u < Average TPT <= 23,5u	-30	20 points	Average TPT > 23,5u	-40	10 points	AVERAGE THROUGHPUT TIME OUTSIDE THE NETHERLANDS			Average TPT <= 18,5 uur	+20	70 points	18,5u < Average TPT <= 19,5u	+10	60 points	19,5u < Average TPT <= 20,5u	+5	55 points	20,5u < Average TPT <= 21,5u	+2	52 points	<b>21,5u &lt; Average TPT &lt;= 22,5u</b>	<b>0</b>	<b>50 points</b>	22,5u < Average TPT <= 23,5u	-10	40 points	23,5u < Average TPT <= 24,5u	-25	25 points	Average TPT > 24,5u	-40	10 points	<b>10 - 70</b>	
For Laboratories situated in the Netherlands;	< 16 hours																																																																							
For Laboratories situated outside the Netherlands.	< 22 hours																																																																							
AVERAGE THROUGHPUT TIME IN THE NETHERLANDS																																																																								
Average TPT <= 12,5 uur	+20	70 points																																																																						
12,5u < Average TPT <= 13,5u	+10	60 points																																																																						
13,5u < Average TPT <= 14,5u	+5	55 points																																																																						
14,5u < Average TPT <= 15,5u	+2	52 points																																																																						
<b>15,5u &lt; Average TPT &lt;= 16,5u</b>	<b>0</b>	<b>50 points</b>																																																																						
16,5u < Average TPT <= 17,5u	-2	48 points																																																																						
17,5u < Average TPT <= 18,5u	-5	45 points																																																																						
18,5u < Average TPT <= 19,5u	-10	40 points																																																																						
20,5u < Average TPT <= 21,5u	-15	35 points																																																																						
21,5u < Average TPT <= 22,5u	-20	30 points																																																																						
22,5u < Average TPT <= 23,5u	-30	20 points																																																																						
Average TPT > 23,5u	-40	10 points																																																																						
AVERAGE THROUGHPUT TIME OUTSIDE THE NETHERLANDS																																																																								
Average TPT <= 18,5 uur	+20	70 points																																																																						
18,5u < Average TPT <= 19,5u	+10	60 points																																																																						
19,5u < Average TPT <= 20,5u	+5	55 points																																																																						
20,5u < Average TPT <= 21,5u	+2	52 points																																																																						
<b>21,5u &lt; Average TPT &lt;= 22,5u</b>	<b>0</b>	<b>50 points</b>																																																																						
22,5u < Average TPT <= 23,5u	-10	40 points																																																																						
23,5u < Average TPT <= 24,5u	-25	25 points																																																																						
Average TPT > 24,5u	-40	10 points																																																																						
<b>Reliability of maximum throughput time per week*</b> <table><tr><td>For Laboratories situated in the Netherlands;</td><td>24 hours</td></tr><tr><td>For Laboratories situated outside the Netherlands.</td><td>24 hours</td></tr></table> - bonus 15 points only if all weekly tested Samples are reported within the maximum throughput time (Yes = 15 points/No = 0 points)	For Laboratories situated in the Netherlands;	24 hours	For Laboratories situated outside the Netherlands.	24 hours	<b>0 or 15</b>																																																																			
For Laboratories situated in the Netherlands;	24 hours																																																																							
For Laboratories situated outside the Netherlands.	24 hours																																																																							
<b>Timely information (&lt; 2 hours after arrival of Samples in Lab) of Sampling hubs/LCDK about disturbances leading to delays in publication of test results</b> - bonus 5 points if no structural complaints** from Sampling hubs	<b>0 or 5</b>																																																																							

<b>Timely information (&gt; 2 hours in advance) of Sampling hubs about delays in supplies or Sample collection</b> - bonus 5 points if no structural complaints** from Sampling hubs	<b>0 or 5</b>	
<b>No non-performance caused by shortage in agreed safety stocks</b> - bonus 5 points if no disturbance of sample taking/testing process	<b>0 or 5</b>	
<b>Total points</b>	<b>100</b>	

<b>The minimum total score per week is:</b>	<b>70</b>
---	-----------

\* based on information derived from CoronIT)

\*\* Structural complaints mean more non-conformances from two or more Sampling hubs reported by e-mail to Contractor [complaints@contractor.com](mailto:complaints@contractor.com) with a copy to [contractmanagement@lcdk.nl](mailto:contractmanagement@lcdk.nl) and accepted by Contractor.

Dienst Testen will weekly set up and share with Contractor the score on KPI's. The score on KPI's, Strategic and operational issues will be discussed in meetings every week.

**Annex D**

**General Government Purchasing Conditions 2018 (ARIV 2018)**

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5.1.2.e

## Validation Report

### KingFisherFlex Purification system MagMAX Viral/PathogenII Nucleic Acid Isolation Kit - Thermo Fisher TaqPath Covid-19 CE-IVD RT-PCR Kit

#### Workflow 3

In-house validation and performance characteristics of a high throughput COVID-19 qRT-PCR detection assay

Author: 5.1.2.e and 5.1.2.e

Date: 22-7-2021

Pure Medical  
Olivier van Noortlaan 120  
3133 AT Vlaardingen, The Netherlands

#### Aim of the Validation

The aim of the validation is to determine the analytical sensitivity and specificity of the analysis system for the SARS-CoV-2 virus.

The KingFisherFlex Purification system MagMAX Viral/PathogenII Nucleic Acid Isolation Kit is intended as an accessory for the isolation of viral RNA and DNA (depending on the viral particle to be detected) from nasopharyngeal and oropharyngeal swab samples. The excellent purity of the isolated viral nucleic acid isolated with the KingFisherFlex Purification system MagMAX Viral/PathogenII Nucleic Acid Isolation Kit allows immediate use in downstream applications such as PCR assays.

The verification is done with two different analytical runs by two different technicians.

#### Sample Collection and Transport

Human oropharyngeal and nose samples swabs arrive in viral transport medium (EchoSafe). All samples should be processed as soon as possible.

The oropharyngeal and nose samples were taken by GGD personnel.

Samples were collected at the GGD and transported to Pure Medical laboratory. All samples were individually divided into two parts: half of the sample specimen was analysed by reference laboratory Novogenia and the other half was analysed at our laboratory.

### RNA extraction

The RNA extraction takes place with the MagMAX™ Viral/Pathogen II Nucleic Acid Isolation Kit (Cat. No. A48383) with the King Fisher FLEX purification system. This kit is designed to recover RNA and DNA from viral particles contained in viral transport media (VTM). The kit utilizes MagMAX™ magnetic-bead technology, ensuring reproducible recovery of high-quality nucleic acid.

### RT-qPCR

RT-qPCR was performed using the Thermo Fisher TaqPath Covid-19 CE-IVD RT-PCR Kit which targets the highly conserved regions within the viral genes N-Gene, ORF-1 ab and S-Gene. All reactions were performed on a Quantstudio 5 (Thermo Fisher Scientific, Applied Biosystems) according to the manufacturers' instructions in 10µl reaction volumes on 96-well PCR plates. Clinical samples were confirmed using the PhoenixDx® 2019-nCoV RNA Detection Kit / Analytica REALQUALITY RQ-2019-nCoV.

### **WORKFLOW 3 KingFisherFlex Purification system MagMAX Viral/PathogenII Nucleic Acid Isolation Kit - Thermo Fisher TaqPath Covid-19 CE-IVD RT-PCR Kit**

We refer to our SOP Workflow 3 for the sample extraction and PCR determination.

We compare our results with the results of the reference laboratory Novogenia, see the table.

### WORKFLOW INFORMATION

Workflow	Sample Preparation Kit	PCR Instrument	PCR Kit	Sample Preparation Duration	PCR Duration
Comparator (Reference) Novogenia	Qiagen Biorobot Universal/Quick viral kit ( Zymo research)	Viaa 7 - Thermofisher	Phoenix-DX-SARS-COV 2	Setup: 15 min.	Setup: 15 min.
				Run: 45 min.	Run: 138 min.
					Analyses: 5 min.
Thermo Fisher TaqPath Covid-19 CE-IVD RT-PCR Kit / KingFisherFlex Purification system MagMAX Viral/PathogenII Nucleic Acid Isolation Kit	MagMAX Viral/PathogenII Nucleic Acid Isolation Kit	Quant Studio 5	Thermo Fisher TaqPath Covid-19 CE-IVD RT-PCR Kit	Setup: 15 min.	Setup: 5 min.
				Run: 23 min.	Run: 64 min.
					Analyses: 5 min.

SP Platform / SP Kit	Sample/Elution volume
Qiagen Biorobot Universal/Quick viral kit ( Zymo research)	200 µL / 60 µL
MagMAX Viral/PathogenII Nucleic Acid Isolation Kit	200µL / 50 µL

PCR Instrument	PCR Kit	Method
Viaa 7 - Thermofisher	Phoenix-DX-SARS-COV 2	Qual.
Quant Studio 5	Thermo Fisher TaqPath Covid-19 CE-IVD RT-PCR Kit	Qual.

## CLINICAL SAMPLE INFORMATION

139 human respiratory swab samples for analysis

## RESULTS

MagMAX Viral/Pathogen II Nucleic Acid Isolation Kit

		Phoenix-DX-SARS-COV 2		
$\Sigma$ Samples	139	POS	NEG	$\Sigma$ Samples
Thermo Fisher TaqPath Covid-19 CE-IVD RT-PCR Kit	POS	83	0	83
	NEG	0	56	56
	$\Sigma$ Samples	83	56	139

Thermo Fisher TaqPath Covid-19 CE-IVD RT-PCR Kit	AVE (Ct IC)	SD (Ct IC)	MIN (Ct IC)	MAX (Ct IC)
Swabs	31	1	29	31

Sensitivity 100%  
Specificity 100%

## Conclusion

We established and validated a high throughput SARS-CoV-2 qRT-PCR detection workflow using samples collected in transport media, a magnetic bead RNA isolation, and a 96-well format RT-qPCR assay. Our data demonstrates a highly specific and sensitive workflow ensuring a robust analytical and clinical performance. We conclude that the results are according to standards.

**Van:** 5.1.2.e  
**Aan:** 5.1.2.e @gmail.com; 5.1.2.e @puregeneticlifestyle.com; 5.1.2.e @puregeneticlifestyle.com  
**Cc:** 5.1.2.e  
**Onderwerp:** 5.1.2.e inspectiebezoek Pure Medical; Vorderen van informatie  
**Datum:** donderdag 18 november 2021 11:05:00  
**Bijlagen:** [image002.jpg](#)  
[20211104\\_Pure\\_Medical bezoeksverslag 4 november.pdf](#)

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The screenshot shows the iProva search interface. The search bar at the top contains the word 'infuus'. Below the search bar, the results are listed as '190 resultaten'. On the left side, there are several filter categories: TITEL, TYPE, VERZAMELING, LOCATIE, and ONDERWERP. The 'VERZAMELING' category is highlighted with a red arrow pointing to it, and a red text box next to it says 'Met deze velden kun je de zoekresultaten verder filteren'. The search results are listed on the right side, starting with 'Toedienen medicijnen via subcutane infuuscanule spuitpomp'. A red arrow points to the search bar with the text 'Lijst met documenten op basis van het ingegeven trefwoord'.

The screenshot shows the iProva search interface after filtering. The search bar still contains 'infuus'. The results are now '4 resultaten'. The 'VERZAMELING' filter is selected, and the results are filtered to show only documents related to 'Cliënt processen'. The results are listed on the right side, starting with 'Handhygiëne - nieuw formaat- def'. A red arrow points to the 'VERZAMELING' filter, and another red arrow points to the 'Cliënt processen' option under it.

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
The screenshot shows the Zenya document viewer interface. At the top, there is a toolbar with several icons and labels: a back arrow, 'Favoriet maken' (Make favorite), 'Opmerking toevoegen' (Add comment) which is highlighted with a red box, 'Afdrukken' (Print), 'Volledig scherm' (Full screen), 'Versie vergelijken' (Compare versions), and 'Meer' (More). Below the toolbar, the document title 'Meldcode Huiselijk Geweld (Versie 2)' is displayed. On the left side, there is a vertical sidebar with the label 'Inhoudsopgave' (Table of contents). The main content area features the 'Accolade Zorg' logo with the tagline 'met elkaar'. Below the logo is a table with the following information:

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- MagMAX™ Viral/Pathogen II Nucleic Acid Isolation & King Fisher FLEX purification system
- Thermo Fisher TaqPath Covid-19 CE-IVD RT-PCR Kit
- Quantstudio 5

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# TaqPath™ COVID-19 CE-IVD RT-PCR Kit

## INSTRUCTIONS FOR USE

Multiplex real-time RT-PCR test intended for the qualitative detection of nucleic acid from SARS-CoV-2

**Catalog Number** A48067

**Publication Number** MAN0019215

**Revision** G.0

For In Vitro Diagnostic Use.

**ThermoFisher**  
S C I E N T I F I C

Life Technologies Corporation |  
6055 Sunol Blvd |  
Pleasanton, California 94566 USA

Life Technologies Europe B.V.  
Kwartsweg 2, 2665 NN Bleiswijk  
The Netherlands

For descriptions of symbols on product labels or product documents, go to [thermofisher.com/symbols-definition](https://thermofisher.com/symbols-definition).

The customer is responsible for compliance with regulatory requirements that pertain to their procedures and uses of the instrument.

The information in this guide is subject to change without notice.

**DISCLAIMER:** TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

**Revision history:** Pub. No. MAN0019215 (English)

Revision	Date	Description
G.0	29 March 2021	<ul style="list-style-type: none"><li>Removed the licensing statement covering Limited Use Label Licenses.</li><li>Updated the Intended use, the product description, the assay limitations, and reactivity (inclusivity).</li><li>Updated the contents and storage to include the software.</li><li>Provided new options for RT-PCR plates.</li><li>Added clinical performance evaluation.</li></ul>
F.0	1 March 2021	<ul style="list-style-type: none"><li>Added Applied Biosystems™ COVID-19 Interpretive Software v1.5 and Applied Biosystems™ COVID-19 Interpretive Software v2.5.</li><li>Removed Applied Biosystems™ COVID-19 Interpretive Software v1.3 and Applied Biosystems™ COVID-19 Interpretive Software v2.3.</li><li>Added general laboratory recommendations.</li><li>Updated template files for the following instruments:<ul style="list-style-type: none"><li>7500 Fast Dx Real-Time PCR Instrument</li><li>7500 Fast Real-Time PCR Instrument with SDS Software v1.5.1</li></ul></li><li>Updated instructions to obtain the COVID-19 Interpretive Software.</li><li>Added an optional step to review the amplification curve for each sample in the COVID-19 Interpretive Software.</li><li>Updated wording for the interpretation of results (Table 6 on page 57) and specified that the Ct for MS2 can be &gt;32 if any of the viral targets is positive (Table 15 on page 66).</li></ul>
E.0	15 July 2020	<ul style="list-style-type: none"><li>Added Applied Biosystems™ COVID-19 Interpretive Software v1.3 and Applied Biosystems™ COVID-19 Interpretive Software v2.3.</li><li>Removed Applied Biosystems™ COVID-19 Interpretive Software v1.2, Applied Biosystems™ COVID-19 Interpretive Software v2.1, and Applied Biosystems™ COVID-19 Interpretive Software v2.2.</li><li>Added new guidelines for sample collection, storage, packaging, and shipping to “Assay limitations” on page 13.</li><li>Added C<sub>t</sub> cutoff information (Appendix A, “C<sub>t</sub> cutoff values for assay targets”).</li><li>Updated instructions to seal, vortex, and centrifuge the RT-PCR plates.</li><li>Provided additional details for naming negative control wells in 384-well plates with unique names.</li><li>Updated reactivity (inclusivity) (page 60).</li><li>Added more information about Limit of Detection experiments to “Limit of detection (LoD)” on page 58.</li></ul>
D.0	2 June 2020	<ul style="list-style-type: none"><li>Removed instructions to mix by pipetting up and down 10 times when preparing RT-PCR plates. Added instructions to vortex the plates for 10–20 seconds to ensure proper mixing.</li><li>Added instructions to unseal and reseal one extraction plate at a time when preparing 384-well RT-PCR plates.</li><li>Updated instructions to create a unique name for each well in the physical plate, not only the wells with a patient sample.</li><li>Added Applied Biosystems™ QuantStudio™ 5 Real-Time PCR Instrument, 384-well block.</li></ul>

Revision	Date	Description
D.0 (continued)	2 June 2020	<ul style="list-style-type: none"> <li>Added QuantStudio™ 5 Dx IVD Software for the Applied Biosystems™ QuantStudio™ 5 Dx Real-Time PCR Instrument.</li> <li>Added COVID-19 Interpretive Software v2.2.</li> <li>Updated retesting requirements for samples with inconclusive results.</li> </ul>
C.0	1 May 2020	<ul style="list-style-type: none"> <li>Updated specimen types to upper respiratory specimens (such as nasopharyngeal, oropharyngeal, nasal, and mid-turbinate swabs, and nasopharyngeal aspirate) and bronchoalveolar lavage (BAL) specimens ("Intended Use" on page 7, "Assay limitations" on page 13, and "Workflow" on page 15).</li> <li>Added procedures to perform RT-PCR with the following instruments: Applied Biosystems™ 7500 Real-Time PCR Instrument, Applied Biosystems™ QuantStudio™ 5 Real-Time PCR Instruments (96-well 0.1-mL block and 96-well 0.2-mL block), Applied Biosystems™ QuantStudio™ 5 Dx Real-Time PCR Instrument, and QuantStudio™ 7 Flex Real-Time PCR Instrument (384-well block).</li> <li>Added new software versions: Applied Biosystems™ COVID-19 Interpretive Software v1.2 and Applied Biosystems™ COVID-19 Interpretive Software v2.1.</li> <li>Changed the shelf life of TaqPath™ COVID-19 Control Dilution Buffer from 9 months to 12 months.</li> <li>Updated guidelines extracted RNA (Chapter 3, "Guidelines for the RNA extracted from the sample").</li> <li>Updated requirements for controls based on the addition of 384-well real-time RT-PCR plates ("Guidelines for RT-PCR" on page 20).</li> <li>Added reaction plates without a barcode as an option for real-time RT-PCR (MicroAmp™ Fast Optical 96-Well Reaction Plate, 0.1 mL, MicroAmp™ Optical 96-Well Reaction Plate, 0.2 mL, and MicroAmp™ Optical 384-Well Reaction Plate).</li> <li>Created separate procedures to prepare the RT-PCR reactions based on the volume of sample input (<math>\leq 200 \mu\text{L}</math> or <math>&gt; 200 \mu\text{L}</math>), the RT-PCR plate (96-well or 384-well), and added specific instructions to vortex and centrifuge the reaction plate.</li> <li>Specified that retesting must be done with the original sample ("Interpretation of the results" on page 57).</li> <li>Updated control requirements in interpretation of results, based on addition of 384-well real-time RT-PCR plates ("Interpretation of the results" on page 57).</li> <li>Added "Interfering substances" on page 60.</li> <li>Removed references to COVID-19 Interpretive Software v1.1.</li> </ul>
B.0	27 March 2020	<ul style="list-style-type: none"> <li>Corrected description of MS2 Phage Control in 'Product Description' from DNA Control to RNA Control.</li> <li>Expanded description of the service and technical support information available at <a href="https://www.thermofisher.com/contactus">https://www.thermofisher.com/contactus</a> in 'Customer and technical support'.</li> <li>Changed the shelf life of the TaqPath™ COVID-19 Control Dilution Buffer.</li> </ul>
A.0	24 March 2020	New document.

**TRADEMARKS:** All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. TaqMan is a registered trademark of Roche Molecular Systems, Inc., used under permission and license. Nasacort is a trademark of AVENTISUB LLC. Dymista is a trademark of Meda Pharmaceuticals Inc. NeilMed and Nasogel are trademarks of NeilMed Products, Inc. Chloraseptic is a trademark of Medtech Products Inc. Bactroban is a trademark of GLAXOSMITHKLINE LLC. Similasan is a trademark of Similasan AG Corporation Switzerland. Spectrum Solutions is a trademark of Spectrum Solutions, LLC. Roche is a trademark of Hoffmann-LaRoche Inc. Cobas is a trademark of Roche Diagnostics Operations, Inc. Illumina is a trademark of Illumina, Inc.

©2021 Thermo Fisher Scientific Inc. All rights reserved.

# Contents

■	<b>CHAPTER 1</b>	<b>TaqPath™ COVID-19 CE-IVD RT-PCR Kit</b>	
		product information .....	7
		Intended Use .....	7
		Product description .....	7
		Contents and storage .....	8
		Required materials not supplied .....	9
		Instrument and software compatibility .....	11
		General laboratory recommendations .....	12
		Assay limitations .....	13
		Samples and controls .....	14
		Workflow .....	15
■	<b>CHAPTER 2</b>	<b>Before you begin .....</b>	<b>16</b>
		Warnings and precautions .....	16
		Sample collection, transport, and storage .....	17
■	<b>CHAPTER 3</b>	<b>Guidelines for the RNA extracted from the sample .....</b>	<b>18</b>
		Sample input volumes .....	18
		Guidelines for RNA extraction .....	18
		Add MS2 Phage Control .....	19
■	<b>CHAPTER 4</b>	<b>Prepare RT-PCR reactions .....</b>	<b>20</b>
		Guidelines for RT-PCR .....	20
		Prepare the RT-PCR reactions (≤200-μL sample input volume, 96-well reaction plate) ...	21
		Prepare the RT-PCR reactions (≤200-μL sample input volume, 384-well reaction plate) ..	22
		Prepare the RT-PCR reactions (>200-μL sample input volume, 96-well reaction plate) ...	25
		Prepare the RT-PCR reactions (>200-μL sample input volume, 384-well reaction plate) ..	27

■	<b>CHAPTER 5</b>	Perform RT-PCR using the Applied Biosystems™ 7500 Fast Dx Real-Time PCR Instrument .....	29
		Dye calibration for the 7500 Real-Time PCR Instrument series .....	29
		Transfer the template (SDT) file for the 7500 Fast Dx Real-Time PCR Instrument .....	29
		Set up and run the 7500 Fast Dx Real-Time PCR Instrument .....	30
■	<b>CHAPTER 6</b>	Perform RT-PCR using the Applied Biosystems™ 7500 Fast Real-Time PCR Instrument .....	32
		Dye calibration for the 7500 Real-Time PCR Instrument series .....	32
		Transfer the template (SDT or EDT) file for the 7500 Fast Real-Time PCR Instrument ....	33
		Set up and run the 7500 Fast Real-Time PCR Instrument (SDS Software v1.5.1) .....	33
		Set up and run the 7500 Fast Real-Time PCR Instrument (7500 Software v2.3) .....	35
■	<b>CHAPTER 7</b>	Perform RT-PCR using the Applied Biosystems™ 7500 Real-Time PCR Instrument .....	37
		Dye calibration for the 7500 Real-Time PCR Instrument series .....	37
		Transfer the template (EDT) file for the 7500 Real-Time PCR Instrument .....	38
		Set up and run the 7500 Real-Time PCR Instrument .....	38
■	<b>CHAPTER 8</b>	Perform RT-PCR using the Applied Biosystems™ QuantStudio™ 5 Real-Time PCR Instrument .....	40
		Dye calibration for the QuantStudio™ 5 Real-Time PCR Instrument .....	40
		Transfer the template (EDT) file for the QuantStudio™ 5 Real-Time PCR Instrument ....	41
		Set up and run the QuantStudio™ 5 Real-Time PCR Instrument (96-well plates) .....	41
		Set up and run the QuantStudio™ 5 Real-Time PCR Instrument (384-well plates) .....	43
■	<b>CHAPTER 9</b>	Perform RT-PCR using the Applied Biosystems™ QuantStudio™ 5 Dx Real-Time PCR Instrument .....	45
		Dye calibration for the QuantStudio™ 5 Dx Real-Time PCR Instrument .....	45
		QuantStudio™ 5 Dx TD Software .....	46
		Transfer the template (EDT) file for the QuantStudio™ 5 Dx Real-Time PCR Instrument (QuantStudio™ 5 Dx TD Software) .....	46
		Set up and run the QuantStudio™ 5 Dx Real-Time PCR Instrument (QuantStudio™ 5 Dx TD Software) .....	46
		QuantStudio™ 5 Dx IVD Software .....	49
		Transfer the template (EDT) file for the QuantStudio™ 5 Dx Real-Time PCR Instrument (QuantStudio™ 5 Dx IVD Software) .....	49
		Install the template file in the QuantStudio™ 5 Dx IVD Software .....	49
		Set up and run the QuantStudio™ 5 Dx Real-Time PCR Instrument (QuantStudio™ 5 Dx IVD Software) .....	50

■	<b>CHAPTER 10</b>	Perform RT-PCR using the Applied Biosystems <sup>™</sup> QuantStudio <sup>™</sup> 7 Flex Real-Time PCR Instrument (384-well block) .....	51
		Dye calibration for the QuantStudio <sup>™</sup> 7 Flex Real-Time PCR Instrument .....	51
		Transfer the template (EDT) file for the QuantStudio <sup>™</sup> 7 Flex Real-Time PCR Instrument (384-well block) .....	52
		Set up and run the QuantStudio <sup>™</sup> 7 Flex Real-Time PCR Instrument (384-well block) ...	52
■	<b>CHAPTER 11</b>	Analysis and results .....	55
		Obtain the Applied Biosystems <sup>™</sup> COVID-19 Interpretive Software .....	55
		Analyze the data .....	55
		Interpretation of the results .....	57
■	<b>CHAPTER 12</b>	Performance characteristics .....	58
		Limit of detection (LoD) .....	58
		Reactivity (Inclusivity) .....	60
		Interfering substances .....	60
		Cross-reactivity .....	62
		Clinical evaluation .....	63
		Clinical performance evaluation .....	64
■	<b>APPENDIX A</b>	C <sub>t</sub> cutoff values for assay targets .....	66
■	<b>APPENDIX B</b>	Documentation and support .....	67
		Related documentation .....	67
		Customer and technical support .....	68
		Limited product warranty .....	68

# TaqPath™ COVID-19 CE-IVD RT-PCR Kit product information

## Intended Use

The TaqPath™ COVID-19 CE-IVD RT-PCR Kit contains the reagents and controls for a real-time reverse transcription polymerase chain reaction (RT-PCR) test intended for the qualitative detection of nucleic acid from SARS-CoV-2 in upper respiratory specimens (such as nasopharyngeal, oropharyngeal, nasal, and mid-turbinate swabs, and nasopharyngeal aspirate), bronchoalveolar lavage (BAL) specimens, and saliva specimens from individuals suspected of COVID-19. For *in vitro* diagnostic use.

Results are for the identification of SARS-CoV-2 RNA. The Applied Biosystems™ COVID-19 Interpretive Software CE-IVD Edition must be used for interpretation of the results generated by the assay. The SARS-CoV-2 RNA is generally detectable in upper respiratory, bronchoalveolar lavage (BAL), and saliva specimens during the acute phase of infection. Positive results are indicative of the presence of SARS-CoV-2 RNA; clinical correlation with patient history and other diagnostic information is necessary to determine patient infection status.

Positive results do not rule out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease. Laboratories may be required to report all positive results to the appropriate Competent Health Authorities.

Negative results do not preclude SARS-CoV-2 infection and should not be used as the sole basis for patient management decisions. Negative results must be combined with clinical observations, patient history, and epidemiological information.

Testing with the TaqPath™ COVID-19 CE-IVD RT-PCR Kit is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of real-time PCR and *in vitro* diagnostic procedures.

## Product description

TaqPath™ COVID-19 CE-IVD RT-PCR Kit includes the following components:

- TaqPath™ COVID-19 RT-PCR Kit, 1000 reactions
  - TaqPath™ COVID-19 Assay Multiplex—Multiplexed assays that contain three primer/probe sets specific to different SARS-CoV-2 genomic regions and primers/probes for bacteriophage MS2.
  - MS2 Phage Control—RNA control to verify the efficacy of the sample preparation and the absence of inhibitors in the PCR reaction. To perform the control, add MS2 Phage Control to the samples before extraction of the RNA.
- TaqPath™ COVID-19 Control—RNA control that contains targets specific to the SARS-CoV-2 genomic regions targeted by the assays.

- TaqPath™ COVID-19 Control Dilution Buffer—Dilution buffer to formulate the TaqPath™ COVID-19 Control.
- TaqPath™ 1-Step Multiplex Master Mix (No ROX)—Ready-to-use PCR mix, including polymerase, dNTP, and salt buffer.
- Package insert—Provides the instructions and the link to download the Instructions for Use and the Applied Biosystems™ COVID-19 Interpretive Software.

## Contents and storage

Table 1 TaqPath™ COVID-19 CE-IVD RT-PCR Kit, 1,000 reactions (Cat. No. [A48067](#))

Component	Contents	Quantity	Volume per tube or bottle	Unit part number	Storage	Shelf life <sup>[1]</sup>
TaqPath™ COVID-19 RT-PCR Kit	TaqPath™ COVID-19 Assay Multiplex (Gene ORF1ab, N Protein, S Protein, MS2)	1 tube	1,500 µL (1,000 reactions)	100093311	–30°C to –10°C	12 months
	MS2 Phage Control	10 tubes	1 mL	100093312	–30°C to –10°C	12 months
TaqPath™ COVID-19 Control		10 tubes	10 µL (1 x 10 <sup>4</sup> copies / µL)	100093314	≤ –70°C	12 months
TaqPath™ COVID-19 Control Dilution Buffer		10 tubes	250 µL	100093291	–30°C to –10°C	12 months
TaqPath™ 1-Step Multiplex Master Mix (No ROX)		1 bottle	10 mL	A48111	–30°C to –10°C	12 months
Applied Biosystems™ COVID-19 Interpretive Software CE-IVD Edition v1.5		<ul style="list-style-type: none"><li>• “Instrument and software compatibility” on page 11</li><li>• “Obtain the Applied Biosystems™ COVID-19 Interpretive Software” on page 55</li></ul>		100098243	n/a <sup>[2]</sup>	n/a <sup>[2]</sup>
Applied Biosystems™ COVID-19 Interpretive Software CE-IVD Edition v2.5				100098246	n/a <sup>[2]</sup>	n/a <sup>[2]</sup>

<sup>[1]</sup> The shelf life of the kit is determined by the component with the shortest shelf life.

<sup>[2]</sup> n/a: not applicable

## Required materials not supplied

Unless otherwise indicated, all materials are available through [thermofisher.com](https://www.thermofisher.com). "MLS" indicates that the material is available from [fisherscientific.com](https://www.fisherscientific.com) or another major laboratory supplier.

Catalog numbers that appear as links open the web pages for those products.

Item	Source
<b>Real-time PCR instrument</b>	
Applied Biosystems™ 7500 Fast Dx Real-Time PCR Instrument (used with SDS Software v1.4.1)	4406984 (with laptop computer) 4406985 (with tower computer)
Applied Biosystems™ 7500 Fast Real-Time PCR Instrument (used with SDS Software v1.5.1 or 7500 Software v2.3)	<a href="#">4351106</a> (with laptop computer) <a href="#">4351107</a> (with desktop computer)
Applied Biosystems™ 7500 Real-Time PCR Instrument (used with 7500 Software v2.3)	<a href="#">4351104</a> (with laptop computer) <a href="#">4351105</a> (with desktop computer)
Applied Biosystems™ QuantStudio™ 5 Real-Time PCR Instrument, 0.2-mL block (used with QuantStudio™ Design and Analysis Desktop Software v1.5.1)	A28569 (with laptop computer) A28574 (with desktop computer) A28139 (instrument only)
Applied Biosystems™ QuantStudio™ 5 Real-Time PCR Instrument, 0.1-mL block (used with QuantStudio™ Design and Analysis Desktop Software v1.5.1)	A28568 (with laptop computer) A28573 (with desktop computer) A28138 (instrument only)
Applied Biosystems™ QuantStudio™ 5 Real-Time PCR Instrument, 384-well block (used with QuantStudio™ Design and Analysis Desktop Software v1.5.1)	A28570 (with laptop computer) A28575 (with desktop computer) A28140 (instrument only)
Applied Biosystems™ QuantStudio™ 5 Dx Real-Time PCR Instrument (used with QuantStudio™ 5 Dx TD Software v1.0.2 or QuantStudio™ 5 Dx IVD Software v1.0.2)	A32005 (with laptop computer) A32006 (with tower computer)
Applied Biosystems™ QuantStudio™ 7 Flex Real-Time PCR Instrument, 384-well block (used with QuantStudio™ Real-Time PCR Software v1.3)	4485695 (with laptop computer) 4485701 (with desktop computer)
<b>Equipment</b>	
Laboratory freezers <ul style="list-style-type: none"> <li>• -30°C to -10°C</li> <li>• ≤ -70°C</li> </ul>	MLS
Centrifuge, with a rotor for microplates	MLS

(continued)

Item	Source
Microcentrifuge	MLS
Laboratory mixer, vortex or equivalent	MLS
Single and multichannel adjustable pipettors (1.00 µL to 1,000.0 µL)	MLS
Cold block (96-well or 384-well) or ice	MLS
<b>Kits and reagents</b>	
Nuclease-free Water (not DEPC-Treated)	MLS
<b>Calibration plates (7500 Real-Time PCR Instrument series)</b>	
ABY™ Dye Spectral Calibration Plate for Multiplex qPCR, Fast 96-well (0.1-mL)	A24734
JUN™ Dye Spectral Calibration Plate for Multiplex qPCR, Fast 96-well (0.1-mL)	A24735
ABY™ Dye Spectral Calibration Plate for Multiplex qPCR, 96-well (0.2-mL)	A24738
JUN™ Dye Spectral Calibration Plate for Multiplex qPCR, 96-well (0.2-mL)	A24737
<b>Calibration plates (QuantStudio™ 7 Flex Real-Time PCR Instrument)</b>	
ABY™ Dye Spectral Calibration Plate for Multiplex qPCR, 384-well	A24736
JUN™ Dye Spectral Calibration Plate for Multiplex qPCR, 384-well	A24733
<b>Tubes, plates, and other consumables</b>	
MicroAmp™ Fast Optical 96-Well Reaction Plate with Barcode, 0.1 mL	4346906, 4366932
MicroAmp™ Fast Optical 96-Well Reaction Plate, 0.1 mL	<a href="#">4346907</a>
MicroAmp™ Optical 96-Well Reaction Plate with Barcode, 0.2 mL	<a href="#">4306737</a> , <a href="#">4326659</a>
MicroAmp™ EnduraPlate™ Optical 96-Well Fast Clear Reaction Plates with Barcode	4483485, 4483494
MicroAmp™ EnduraPlate™ Optical 96-Well Fast Clear Reaction Plates	<a href="#">A36930</a>
MicroAmp™ Optical 96-Well Reaction Plate, 0.2 mL	<a href="#">N8010560</a> , <a href="#">4316813</a>

(continued)

Item	Source
MicroAmp™ EnduraPlate™ Optical 96-Well Clear Reaction Plates with Barcode	<a href="#">4483354</a> , <a href="#">4483352</a>
MicroAmp™ EnduraPlate™ Optical 96-Well Clear Reaction Plates	<a href="#">A36924</a>
MicroAmp™ Optical 384-Well Reaction Plate with Barcode	4309849, 4326270, 4343814
MicroAmp™ Optical 384-Well Reaction Plate	4343370
MicroAmp™ EnduraPlate™ Optical 384-Well Clear Reaction Plates with Barcode	<a href="#">4483285</a> , <a href="#">4483273</a>
MicroAmp™ EnduraPlate™ Optical 384-Well Clear Reaction Plates	<a href="#">A36931</a>
MicroAmp™ Optical Adhesive Film	<a href="#">4311971</a> , <a href="#">4360954</a>
MicroAmp™ Adhesive Film Applicator	<a href="#">4333183</a>
Nonstick, RNase-free microcentrifuge tubes (1.5 mL and 2.0 mL)	<a href="https://thermofisher.com/plastics">thermofisher.com/plastics</a>
Sterile aerosol barrier (filtered) pipette tips	<a href="https://thermofisher.com/pipettetips">thermofisher.com/pipettetips</a>

## Instrument and software compatibility

The following table lists the version of the Applied Biosystems™ COVID-19 Interpretive Software that is compatible with your instrument and its associated analysis software.

For information on how to obtain the Applied Biosystems™ COVID-19 Interpretive Software, see “Obtain the Applied Biosystems™ COVID-19 Interpretive Software” on page 55.

To obtain the analysis software or firmware for use with your real-time PCR instrument, go to [thermofisher.com/qpcrsoftware](https://thermofisher.com/qpcrsoftware), then select your instrument in the **Real-Time PCR** section.

Instrument	Analysis software used with the instrument	Compatible COVID-19 Interpretive Software version
7500 Fast Dx Real-Time PCR Instrument	SDS Software v1.4.1	v1.5
7500 Fast Real-Time PCR Instrument	SDS Software v1.5.1 or 7500 Software v2.3	v1.5
7500 Real-Time PCR Instrument	7500 Software v2.3	v1.5

(continued)

Instrument	Analysis software used with the instrument	Compatible COVID-19 Interpretive Software version
QuantStudio™ 5 Real-Time PCR Instrument with instrument firmware v1.3.3 96-well, 0.2-mL block 96-well, 0.1-mL block	QuantStudio™ Design and Analysis Desktop Software v1.5.1	v2.5
QuantStudio™ 5 Real-Time PCR Instrument with instrument firmware v1.3.3 384-well block	QuantStudio™ Design and Analysis Desktop Software v1.5.1	v2.5
QuantStudio™ 5 Dx Real-Time PCR Instrument with instrument firmware v1.0.3	QuantStudio™ 5 Dx TD Software v1.0.2	v2.5
QuantStudio™ 5 Dx Real-Time PCR Instrument with instrument firmware v1.0.3	QuantStudio™ 5 Dx IVD Software v1.0.2	v2.5
QuantStudio™ 7 Flex Real-Time PCR Instrument with instrument firmware v1.0.4 384-well block	QuantStudio™ Real-Time PCR Software v1.3	v2.5

## General laboratory recommendations

- Implement standard operating procedures in your laboratory to prevent contamination, such as the following:
  - Frequent glove changes
  - Frequent decontamination of surfaces, equipment, and pipettes with 10% bleach or decontamination solution, followed by 70% ethanol
  - Use of ultraviolet light during biosafety cabinet decontamination (when available)
- To prevent degradation, keep eluted sample RNA, master mixes, assays, and controls on ice or in cold blocks while in use.
- Limit freeze-thaw cycles.
- Aliquot reagents to prevent stock contamination and reduce the number of freeze-thaw cycles.
- After each run, review the amplification curves in the interpretive software for signs of inadequate vortexing or centrifugation. Contact your Applications Support team for additional information or training on data QC in your instrument software.

## Assay limitations

- This assay is intended to be used for *in vitro* diagnostic purposes only. Follow good laboratory practices and all precautions and guidelines in these user guides to avoid cross-contamination between samples.
- The TaqPath™ COVID-19 CE-IVD RT-PCR Kit performance on upper respiratory specimen was established using nasopharyngeal and oropharyngeal swab, nasopharyngeal aspirate, and bronchoalveolar lavage samples only. Nasal swabs and mid-turbinate swabs are considered acceptable specimen types according to the *Interim Guidelines for Collecting, Handling, and Testing Clinical Specimens from Persons for Coronavirus Disease 2019 (COVID-19)* published by the Centers for Disease Control and Prevention, but performance of the TaqPath™ COVID-19 CE-IVD RT-PCR Kit with these specimen types has not been established. Testing of nasal and mid-turbinate nasal swabs (self-collected under supervision of or collected by a healthcare provider) is limited to patients with symptoms of COVID-19.
- The performance on saliva specimen was established with the TaqPath™ COVID-19 CE-IVD RT-PCR Kit and the MagMAX™ Viral/Pathogen II Nucleic Acid Isolation Kit from saliva collected with the Spectrum Solutions™ SDNA Molecular Diagnostic Saliva Collection Devices.
- Saliva must be collected by passive drooling without coughing or throat clearing. Certain drugs, physiological conditions, or pathologies can reduce or increase salivary flow. The presence of the virus in saliva can be the result of both viral excretion from the salivary glands and the upper and lower respiratory tract. Improper saliva production methodology may hinder the ability of the assay to detect the target sequences. Though self-collected saliva may be acceptable, specimen collection without supervision may result in poor specimen quality, and may prevent the certification of the patient identity.
- Samples must be collected, transported, and stored using appropriate procedures and conditions. Improper collection, transport, or storage of specimens may hinder the ability of the assay to detect the target sequences.
- This kit uses purified RNA as a sample for the analysis. The quality of the RNA recovered from biological samples is essential for the quality of the results generated with this kit.
- False-negative results may arise from:
  - Improper sample collection
  - Degradation of the viral RNA during shipping/storage
  - Specimen collection after nucleic acid can no longer be found in the specimen matrix
  - Using poor extraction method
  - The presence of RT-PCR inhibitors
  - Mutation in the SARS-CoV-2 virus
  - Failure to follow instructions for use
  - Absence of specimen

- False-positive results may arise from:
  - Cross contamination during specimen handling or preparation
  - Cross contamination between patient samples
  - Specimen mix-up
  - RNA contamination during product handling
  - Fluorescence artefact in early amplification cycles caused by improper mixing or bubbles in the reaction mix
- The impacts of vaccines, antiviral therapeutics, antibiotics, chemotherapeutic or immunosuppressant drugs have not been evaluated. The TaqPath™ COVID-19 CE-IVD RT-PCR Kit cannot rule out diseases caused by other bacterial or viral pathogens.
- Negative results do not preclude infection with SARS-CoV-2 virus, and should not be the sole basis of a patient management decision.
- Laboratories may be required to report all positive results to the appropriate Competent Health Authorities.
- Use of saliva may be authorized under specific conditions defined by the Competent Health Authorities (e.g. with collection under supervision only / self-collection, for iterative testing, as secondary intention, under provision of test confirmation in case of positive saliva) and may be subject to restrictions at borders entry point in some Countries (travel). Always follow the applicable regulations of the concerned jurisdiction.

## Samples and controls

Patient samples must be collected according to appropriate laboratory guidelines. Positive and negative test controls must be included to accurately interpret patient test results.

Include the following controls:

Control	Used to monitor	Assays
Positive Control (TaqPath™ COVID-19 Control Kit)	RT-PCR reaction setup and reagent integrity	All three SARS-CoV-2 assays
MS2 Phage Control	RNA extraction and absence of inhibitors in the qPCR reaction	MS2 assay
Negative Control	Cross-contamination during RNA extraction and reaction setup	All three SARS-CoV-2 assays
		MS2 assay

## Workflow

The workflow begins with nucleic acid purified from upper respiratory specimens (such as nasopharyngeal, oropharyngeal, nasal, and mid-turbinate swabs, and nasopharyngeal aspirate) and bronchoalveolar lavage (BAL) specimens.

The purified nucleic acid is reverse transcribed into cDNA and amplified using the TaqPath™ COVID-19 CE-IVD RT-PCR Kit and one of the following real-time PCR instruments:

- Applied Biosystems™ 7500 Fast Dx Real-Time PCR instrument
- Applied Biosystems™ 7500 Fast Real-Time PCR Instrument
- Applied Biosystems™ 7500 Real-Time PCR Instrument
- Applied Biosystems™ QuantStudio™ 5 Real-Time PCR Instrument, 96-well, 0.2-mL block
- Applied Biosystems™ QuantStudio™ 5 Real-Time PCR Instrument, 96-well, 0.1-mL block
- Applied Biosystems™ QuantStudio™ 5 Real-Time PCR Instrument, 384-well block
- Applied Biosystems™ QuantStudio™ 5 Dx Real-Time PCR Instrument
- Applied Biosystems™ QuantStudio™ 7 Flex Real-Time PCR Instrument, 384-well block

In the process, probes anneal to three (3) target sequences that are specific to SARS-CoV-2. Each target is located between unique forward and reverse primers for the following genes:

- ORF1ab
- N Protein
- S Protein

During the extension phase of the PCR cycle, the 5' nuclease activity of Taq polymerase degrades the probe. This degradation causes the reporter dye to separate from the quencher dye, which generates a fluorescent signal. With each cycle, additional reporter dye molecules are cleaved from their respective probes, which increases the fluorescence intensity. Fluorescence intensity is monitored at each PCR cycle by the real-time PCR instrument.

The data are analyzed and interpreted by the Applied Biosystems™ COVID-19 Interpretive Software.

The workflow dedicated to test from saliva is detailed in the *TaqPath™ COVID-19 CE-IVD RT-PCR Kit with saliva Instructions for Use* (Pub. No. MAN0024936).

## Warnings and precautions

The TaqPath™ COVID-19 CE-IVD RT-PCR Kit workflow should be performed by qualified and trained staff to avoid the risk of erroneous results. Use separate areas for the preparation of patient samples and controls to prevent false positive results. Samples and reagents must be handled under a laminar airflow hood or biological safety cabinet.

- Specimens should always be treated as if infectious and/or biohazardous in accordance with safe laboratory procedures.
- Follow necessary precautions when handling specimens. Use personal protective equipment (PPE) consistent with current guidelines for the handling of potentially infectious samples.
- Always use pipette tips with aerosol barriers. Tips that are used must be sterile and free from DNases and RNases.
- Do not eat, drink, smoke, or apply cosmetic products in the work areas.
- Modifications to assay reagents, assay protocol, or instrumentation are not permitted, and are in violation of the *In Vitro Diagnostic Directive 98/79/EC*.
- Do not use the kit after the expiry date.
- Dispose of waste in compliance with the local regulations.
- Safety Data Sheets are available upon request.
- Laboratories may be required to report all positive results to the appropriate Competent Health Authorities.

The following countries require the CE-marked In Vitro Diagnostics: Austria, Belgium, Bulgaria, Croatia, Cyprus, Czechia, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Poland, Portugal, Romania, Slovakia, Slovenia, Spain, Sweden, UK, Norway, Iceland, Liechtenstein, Switzerland, Turkey.

- Positive results are indicative of the presence of SARS-CoV-2 RNA.
- Handle all samples and controls as if they are capable of transmitting infectious agents.
- Reagents must be stored and handled as specified in Table 1.
- The quality of the sample preparation (purified RNA) may influence the quality of the qPCR test. Laboratories shall only use the purification method they have selected. Laboratories that do not have any selected method may use the MagMAX™ Viral/Pathogen Nucleic Acid Isolation Kit or the MagMAX™ Viral/Pathogen II Nucleic Acid Isolation Kit.
- The qPCR kit includes an RNA phage control to verify the efficacy of the sample preparation and the absence of inhibitors in the PCR reaction. To perform the control, add MS2 Phage Control to the samples before extraction of the RNA, according to the recommendations in Chapter 3, “Guidelines for the RNA extracted from the sample”.

## Sample collection, transport, and storage

---

**Note:** Handle all samples and controls as if they are capable of transmitting infectious agents.

---

- Refer to the *Interim Guidelines for Collecting, Handling, and Testing Clinical Specimens from Persons for Coronavirus Disease 2019 (COVID-19)*, published by the Centers for Disease Control and Prevention, for specimen collection and storage guidelines.
- Specimens must be packaged, shipped, and transported according to the current edition of the *International Air Transport Association (IATA) Dangerous Goods Regulations* ([iata.org/en/programs/cargo/dgr](https://www.iata.org/en/programs/cargo/dgr)).

# Guidelines for the RNA extracted from the sample

## Sample input volumes

To meet the need for a broad range of sample input volumes, 2 RT-PCR protocols are provided, depending on the sample input volume.

The following sample input volumes were tested:

- Low input: 200 µL
- High input: 400 µL

Sample input volumes > 400 µL have not been tested with the TaqPath™ COVID-19 CE-IVD RT-PCR Kit.

The following RNA extraction kits were tested:

- MagMAX™ Viral/Pathogen Nucleic Acid Isolation Kit
- MagMAX™ Viral/Pathogen II Nucleic Acid Isolation Kit

## Guidelines for RNA extraction

Extracted RNA can be prepared with any standard RNA extraction procedure or RNA extraction kit. The minimum recommended elution volume is 50 µL.

The following guidelines are defined as optimal for the use of the TaqPath™ COVID-19 CE-IVD RT-PCR Kit.

- The MS2 Phage Control provided in the kit must be used to verify the efficacy of the sample preparation and the absence of inhibitors in the RT-PCR reaction. It is added to the sample before extraction (see “Add MS2 Phage Control” on page 19).
- Use Nuclease-free Water (not DEPC-Treated) containing MS2 Phage Control as the Negative Control. The purified Negative Control is used as the Negative Control for RT-PCR.
- We recommend that the volume of MS2 Phage Control is 2.5% of the sample input volume.

## Add MS2 Phage Control

The MS2 Phage Control must be used to verify the efficacy of the sample preparation and the absence of inhibitors in the RT-PCR reaction.

Add the appropriate volume of MS2 Phage Control to each sample well and to the Negative Control well just before lysis during RNA extraction.

We recommend that the volume of MS2 Phage Control is 2.5% of the sample input volume. Examples are provided in the table below:

Sample input volume	MS2 Phage Control volume
200 µL	5 µL MS2 Phage Control
400 µL	10 µL MS2 Phage Control

# 4

## Prepare RT-PCR reactions

■ Guidelines for RT-PCR .....	20
■ Prepare the RT-PCR reactions ( $\leq 200$ - $\mu$ L sample input volume, 96-well reaction plate) .....	21
■ Prepare the RT-PCR reactions ( $\leq 200$ - $\mu$ L sample input volume, 384-well reaction plate) .....	22
■ Prepare the RT-PCR reactions ( $> 200$ - $\mu$ L sample input volume, 96-well reaction plate) .....	25
■ Prepare the RT-PCR reactions ( $> 200$ - $\mu$ L sample input volume, 384-well reaction plate) .....	27

---

**Note:** The procedure used to prepare the RT-PCR reactions depends on the original sample input volume that was used during the preparation of RNA ( $\leq 200$   $\mu$ L or  $> 200$   $\mu$ L). The procedure used to prepare the RT-PCR reactions also depends on whether a 96-well plate or a 384-well plate is used for RT-PCR.

---

### Guidelines for RT-PCR

---

#### IMPORTANT!

- For each RT-PCR reaction plate, include the following controls:
    - One Positive Control
    - One Negative Control from each extraction run.  
For example, if RNA samples from 4 extraction runs are combined on one 384-well RT-PCR reaction plate, then 4 Negative Control wells must be run on that 384-well reaction plate.
  - Prepare the RT-PCR reaction plate on ice and keep it on ice until it is loaded into the real-time PCR instrument.
  - Run the plate immediately after preparation. Failure to do so could result in degraded RNA samples.
  - To prevent contamination, prepare reagents in a PCR workstation or equivalent amplicon-free area. Do not use the same pipette for controls and RNA samples, and always use aerosol barrier pipette tips.
  - Maintain an RNase-free environment.
  - Protect assays from light.
  - Keep RNA samples and components on ice during use.
-

## Prepare the RT-PCR reactions (≤200-μL sample input volume, 96-well reaction plate)

Use this procedure under the following conditions:

- Original sample input volume of up to 200 μL was used for extraction
  - Instrument is compatible with 96-well RT-PCR reaction plates
1. If frozen, thaw the reagents on ice.
  2. Gently vortex the reagents, then centrifuge briefly to collect liquid at the bottom of the tube.
  3. Dilute TaqPath™ COVID-19 Control ( $1 \times 10^4$  copies/μL) to a working stock of 25 copies/μL:
    - a. Pipet 98 μL of TaqPath™ COVID-19 Control Dilution Buffer into a microcentrifuge tube, then add 2 μL of TaqPath™ COVID-19 Control. Mix well, then centrifuge briefly.
    - b. Pipet 87.5 μL of TaqPath™ COVID-19 Control Dilution Buffer into a second microcentrifuge tube, then add 12.5 μL of the dilution created in substep 3a. Mix well, then centrifuge briefly.

---

**Note:** The TaqPath™ COVID-19 Control does not contain the MS2 template.

---

4. Prepare the Reaction Mix:
  - a. For each run, combine the following components sufficient for the number of RNA samples to be tested plus one Positive Control and one Negative Control.

All volumes include 10% overage for pipette error.

---

**IMPORTANT!** The volumes in this table assume that you extracted sample RNA using an original sample input volume of up to 200 μL.

---

Component	Volume per RNA sample or control	Volume for $n$ RNA samples plus 2 controls	Volume for 94 RNA Samples plus 2 controls
TaqPath™ 1-Step Multiplex Master Mix (No ROX™) (4X)	6.25 μL	$6.875 \times (n + 2)$ μL	660 μL
COVID-19 Real Time PCR Assay Multiplex	1.25 μL	$1.375 \times (n + 2)$ μL	132 μL
Nuclease-free Water	7.50 μL	$8.25 \times (n + 2)$ μL	792 μL
<b>Total Reaction Mix volume</b>	<b>15.0 μL</b>	—	<b>1584 μL</b>

5. Set up the reaction plate:
  - a. Pipet 15.0 μL of the Reaction Mix prepared in step 4 into each well of a MicroAmp™ Fast Optical 96-Well Reaction Plate with Barcode, 0.1 mL or a MicroAmp™ Optical 96-Well Reaction Plate with Barcode, 0.2 mL.

For other reaction plates that can be used, see “Required materials not supplied” on page 9.

- b. Gently vortex the sealed plate containing the purified sample RNA and Negative Control from the RNA extraction procedure, then centrifuge briefly to collect liquid at the bottom of the plate.
- c. Unseal the plate containing the purified sample RNA and Negative Control from the RNA extraction procedure. Add either sample RNA, Negative Control, or Positive Control to each well of the reaction plate according to Table 2 on page 22.
- d. Seal the plate thoroughly with MicroAmp™ Optical Adhesive Film.

**IMPORTANT!** When applying the MicroAmp™ Optical Adhesive Film, ensure that pressure is applied across the entire plate and that there is a tight seal across every individual well. Failure to do so runs the risk of an improperly sealed well, leading to potential well-to-well contamination during vortexing and evaporation during PCR.

- e. Vortex the plate at the highest setting speed for 10–30 seconds with medium pressure. Move the plate around to ensure equal contact on the vortex mixer platform.

**IMPORTANT!** Vortex for 10–30 seconds to ensure proper mixing. Failure to do so might result in false classification of samples.

- f. Centrifuge the reaction plate for 1–2 minutes at  $\geq 650 \times g$  ( $\geq 650$  RCF) to remove bubbles and to collect the liquid at the bottom of the reaction plate.

Table 2 Reaction plate volumes

Component	Volume per reaction		
	RNA Sample reaction	Positive Control reaction	Negative Control reaction
Reaction Mix	15.0 $\mu\text{L}$	15.0 $\mu\text{L}$	15.0 $\mu\text{L}$
Purified sample RNA (from RNA extraction)	10.0 $\mu\text{L}$	—	—
Positive Control (diluted TaqPath™ COVID-19 Control, from step 3)	—	2.0 $\mu\text{L}$	—
Nuclease-free Water	—	8.0 $\mu\text{L}$	—
Purified Negative Control (from RNA extraction)	—	—	10.0 $\mu\text{L}$
<b>Total volume</b>	<b>25.0 <math>\mu\text{L}</math></b>	<b>25.0 <math>\mu\text{L}</math></b>	<b>25.0 <math>\mu\text{L}</math></b>

## Prepare the RT-PCR reactions ( $\leq 200\text{-}\mu\text{L}$ sample input volume, 384-well reaction plate)

Use this procedure under the following conditions:

- Original sample input volume of up to 200  $\mu\text{L}$  was used for extraction
- Instrument is compatible with 384-well RT-PCR reaction plates

1. If frozen, thaw the reagents on ice.
2. Gently vortex the reagents, then centrifuge briefly to collect liquid at the bottom of the tube.
3. Dilute TaqPath™ COVID-19 Control ( $1 \times 10^4$  copies/ $\mu\text{L}$ ) to a working stock of 25 copies/ $\mu\text{L}$ :
  - a. Pipet 98  $\mu\text{L}$  of TaqPath™ COVID-19 Control Dilution Buffer into a microcentrifuge tube, then add 2  $\mu\text{L}$  of TaqPath™ COVID-19 Control. Mix well, then centrifuge briefly.
  - b. Pipet 87.5  $\mu\text{L}$  of TaqPath™ COVID-19 Control Dilution Buffer into a second microcentrifuge tube, then add 12.5  $\mu\text{L}$  of the dilution created in substep 3a. Mix well, then centrifuge briefly.

---

**Note:** The TaqPath™ COVID-19 Control does not contain the MS2 template.

---

4. Prepare the Reaction Mix.
  - a. For each run, combine the following components sufficient for the number of RNA samples, plus one Positive Control per 384-well real-time RT-PCR plate, and one Negative Control from each extraction run.  
  
For example, if RNA samples from 4 extraction runs are being combined on one 384-well real-time RT-PCR plate, then 4 Negative Control wells need to be run on that 384-well real-time RT-PCR plate.

All volumes include 10% overage for pipette error.

---

**IMPORTANT!** The volumes in this table assume that you extracted sample RNA using an original sample input volume of up to 200  $\mu\text{L}$ .

---

Component	Volume per RNA Sample or control	Volume for $n$ RNA Samples plus $y$ Negative Controls plus 1 Positive Control	Volume for 379 RNA Samples plus 4 Negative Controls plus 1 Positive Control
TaqPath™ 1-Step Multiplex Master Mix (No ROX™) (4X)	5.00 $\mu\text{L}$	$5.50 \times (n + y + 1) \mu\text{L}$	2112.0 $\mu\text{L}$
COVID-19 Real Time PCR Assay Multiplex	1.00 $\mu\text{L}$	$1.10 \times (n + y + 1) \mu\text{L}$	422.4 $\mu\text{L}$
Nuclease-free Water	4.00 $\mu\text{L}$	$4.40 \times (n + y + 1) \mu\text{L}$	1690.0 $\mu\text{L}$
<b>Total Reaction Mix volume</b>	<b>10.0 <math>\mu\text{L}</math></b>	—	<b>4224.4 <math>\mu\text{L}</math></b>

5. Set up the reaction plate:
  - a. Pipet 10.0  $\mu\text{L}$  of the Reaction Mix prepared in step 4 into each well of a MicroAmp™ Optical 384-Well Reaction Plate with Barcode.  
For other reaction plates that can be used, see “Required materials not supplied” on page 9.
  - b. Gently vortex the sealed plate containing the purified sample RNA and Negative Control from the RNA extraction procedure, then centrifuge briefly to collect liquid at the bottom of the plate.

- c. Unseal the plate containing the purified sample RNA and Negative Control from the RNA extraction procedure. Add either sample RNA, Negative Control, or Positive Control to each well of the reaction plate according to Table 3 on page 24.

---

**IMPORTANT!** To prevent sample contamination, unseal one extraction plate at a time, then reseal it after adding the samples to the RT-PCR reaction plate.

---

- d. Seal the plate thoroughly with MicroAmp™ Optical Adhesive Film.

---

**IMPORTANT!** When applying the MicroAmp™ Optical Adhesive Film, ensure that pressure is applied across the entire plate and that there is a tight seal across every individual well. Failure to do so runs the risk of an improperly sealed well, leading to potential well-to-well contamination during vortexing and evaporation during PCR.

---

- e. Vortex the plate at the highest setting speed for 10–30 seconds with medium pressure. Move the plate around to ensure equal contact on the vortex mixer platform.

---

**IMPORTANT!** Vortex for 10–30 seconds to ensure proper mixing. Failure to do so might result in false classification of samples.

---

- f. Centrifuge the reaction plate for 1–2 minutes at  $\geq 650 \times g$  ( $\geq 650$  RCF) to remove bubbles and to collect the liquid at the bottom of the reaction plate.

**Table 3** Reaction plate volumes

Component	Volume per reaction		
	RNA Sample reaction	Positive Control reaction	Negative Control reaction
Reaction Mix	10.0 $\mu\text{L}$	10.0 $\mu\text{L}$	10.0 $\mu\text{L}$
Purified sample RNA (from RNA extraction)	10.0 $\mu\text{L}$	—	—
Positive Control (diluted TaqPath™ COVID-19 Control from step 3)	—	2.0 $\mu\text{L}$	—
Nuclease-free Water	—	8.0 $\mu\text{L}$	—
Purified Negative Control (from RNA extraction)	—	—	10.0 $\mu\text{L}$
<b>Total volume</b>	<b>20.0 <math>\mu\text{L}</math></b>	<b>20.0 <math>\mu\text{L}</math></b>	<b>20.0 <math>\mu\text{L}</math></b>

## Prepare the RT-PCR reactions (>200-μL sample input volume, 96-well reaction plate)

Use this procedure under the following conditions:

- Original sample input volume of >200 μL, but ≤400 μL, was used for extraction

**Note:** A sample input volume >400 μL has not been tested with the TaqPath™ COVID-19 CE-IVD RT-PCR Kit.

- Instrument is compatible with 96-well RT-PCR reaction plates

1. If frozen, thaw the reagents on ice.

2. Gently vortex the reagents, then centrifuge briefly to collect liquid at the bottom of the tube.

3. Dilute TaqPath™ COVID-19 Control ( $1 \times 10^4$  copies/μL) to a working stock of 25 copies/μL:

- Pipet 98 μL of TaqPath™ COVID-19 Control Dilution Buffer into a microcentrifuge tube, then add 2 μL of TaqPath™ COVID-19 Control. Mix well, then centrifuge briefly.
- Pipet 87.5 μL of TaqPath™ COVID-19 Control Dilution Buffer into a second microcentrifuge tube, then add 12.5 μL of the dilution created in substep 3a. Mix well, then centrifuge briefly.

**Note:** The TaqPath™ COVID-19 Control does not contain the MS2 template.

4. Prepare the Reaction Mix:

- For each run, combine the following components sufficient for the number of RNA samples to be tested plus one Positive Control and one Negative Control.

All volumes include 10% overage for pipette error.

**IMPORTANT!** The volumes in this table assume that you extracted sample RNA using an original sample input volume between 201 μL and 400 μL.

Component	Volume per RNA Sample or control	Volume for $n$ RNA Samples plus 2 controls	Volume for 94 RNA Samples plus 2 controls
TaqPath™ 1-Step Multiplex Master Mix (No ROX™) (4X)	6.25 μL	$6.875 \times (n + 2)$ μL	660 μL
COVID-19 Real Time PCR Assay Multiplex	1.25 μL	$1.375 \times (n + 2)$ μL	132 μL
Nuclease-free Water	12.50 μL	$13.75 \times (n + 2)$ μL	1320 μL
<b>Total Reaction Mix volume</b>	<b>20.0 μL</b>	—	<b>2112 μL</b>

5. Set up the reaction plate:

- Pipet 20.0 μL of the Reaction Mix prepared in step 4 into each well of a MicroAmp™ Fast Optical 96-Well Reaction Plate with Barcode, 0.1 mL or a MicroAmp™ Optical 96-Well Reaction Plate with Barcode, 0.2 mL.

For other reaction plates that can be used, see “Required materials not supplied” on page 9.

- b. Gently vortex the sealed plate containing the purified sample RNA and Negative Control from the RNA extraction procedure, then centrifuge briefly to collect liquid at the bottom of the plate.
- c. Unseal the plate containing the purified sample RNA and Negative Control from the RNA extraction procedure. Add either sample RNA, Negative Control, or Positive Control to each well of the reaction plate according to Table 4 on page 26.
- d. Seal the plate thoroughly with MicroAmp™ Optical Adhesive Film.

---

**IMPORTANT!** When applying the MicroAmp™ Optical Adhesive Film, ensure that pressure is applied across the entire plate and that there is a tight seal across every individual well. Failure to do so runs the risk of an improperly sealed well, leading to potential well-to-well contamination during vortexing and evaporation during PCR.

---

- e. Vortex the plate at the highest setting speed for 10–30 seconds with medium pressure. Move the plate around to ensure equal contact on the vortex mixer platform.

---

**IMPORTANT!** Vortex for 10–30 seconds to ensure proper mixing. Failure to do so might result in false classification of samples.

---

- f. Centrifuge the reaction plate for 1–2 minutes at  $\geq 650 \times g$  ( $\geq 650$  RCF) to remove bubbles and to collect the liquid at the bottom of the reaction plate.

Table 4 Reaction plate volumes

Component	Volume per reaction		
	RNA Sample reaction	Positive Control reaction	Negative Control reaction
Reaction Mix	20.0 µL	20.0 µL	20.0 µL
Purified sample RNA (from RNA extraction)	5.0 µL	—	—
Positive Control (diluted TaqPath™ COVID-19 Control from step 3)	—	2.0 µL	—
Nuclease-free Water	—	3.0 µL	—
Purified Negative Control (from RNA extraction)	—	—	5.0 µL
<b>Total volume</b>	<b>25.0 µL</b>	<b>25.0 µL</b>	<b>25.0 µL</b>

## Prepare the RT-PCR reactions (>200-μL sample input volume, 384-well reaction plate)

Use this procedure under the following conditions:

- Original sample input volume of >200 μL, but ≤400 μL, was used for extraction

**Note:** A sample input volume >400 μL has not been tested with the TaqPath™ COVID-19 CE-IVD RT-PCR Kit.

- Instrument is compatible with 384-well RT-PCR reaction plates

1. If frozen, thaw the reagents on ice.

2. Gently vortex the reagents, then centrifuge briefly to collect liquid at the bottom of the tube.

3. Dilute TaqPath™ COVID-19 Control ( $1 \times 10^4$  copies/μL) to a working stock of 25 copies/μL:

- Pipet 98 μL of TaqPath™ COVID-19 Control Dilution Buffer into a microcentrifuge tube, then add 2 μL of TaqPath™ COVID-19 Control. Mix well, then centrifuge briefly.
- Pipet 87.5 μL of TaqPath™ COVID-19 Control Dilution Buffer into a second microcentrifuge tube, then add 12.5 μL of the dilution created in substep 3a. Mix well, then centrifuge briefly.

**Note:** The TaqPath™ COVID-19 Control does not contain the MS2 template.

4. Prepare the Reaction Mix.

- For each run, combine the following components sufficient for the number of RNA samples, plus one Positive Control per 384-well real-time RT-PCR plate, and one Negative Control from each extraction run.

For example, if RNA samples from 4 extraction runs are being combined on one 384-well real-time RT-PCR plate, then 4 Negative Control wells need to be run on that 384-well real-time RT-PCR plate.

All volumes include 10% overage for pipette error.

**IMPORTANT!** The volumes in this table assume that you extracted sample RNA using an original sample input volume between 201 μL and 400 μL.

Component	Volume per RNA Sample or control	Volume for $n$ RNA Samples plus $y$ Negative Controls plus 1 Positive Control	Volume for 379 RNA Samples plus 4 Negative Controls plus 1 Positive Control
TaqPath™ 1-Step Multiplex Master Mix (No ROX™) (4X)	5.00 μL	$5.50 \times (n + y + 1)$ μL	2112.0 μL
COVID-19 Real Time PCR Assay Multiplex	1.00 μL	$1.10 \times (n + y + 1)$ μL	422.4 μL
Nuclease-free Water	9.00 μL	$9.90 \times (n + y + 1)$ μL	3802.0 μL
<b>Total Reaction Mix volume</b>	<b>15.0 μL</b>	—	<b>6336.4 μL</b>

## 5. Set up the reaction plate:

- a. Pipet 15.0 μL of the Reaction Mix prepared in step 4 into each well of a MicroAmp™ Optical 384-Well Reaction Plate with Barcode.

For other reaction plates that can be used, see “Required materials not supplied” on page 9.

- b. Gently vortex the sealed plate containing the purified sample RNA and Negative Control from the RNA extraction procedure, then centrifuge briefly to collect liquid at the bottom of the plate.
- c. Unseal the plate containing the purified sample RNA and Negative Control from the RNA extraction procedure. Add either sample RNA, Negative Control, or Positive Control to each well of the reaction plate according to Table 5 on page 28.

**IMPORTANT!** To prevent sample contamination, unseal one extraction plate at a time, then reseal it after adding the samples to the RT-PCR reaction plate.

- d. Seal the plate thoroughly with MicroAmp™ Optical Adhesive Film.

**IMPORTANT!** When applying the MicroAmp™ Optical Adhesive Film, ensure that pressure is applied across the entire plate and that there is a tight seal across every individual well. Failure to do so runs the risk of an improperly sealed well, leading to potential well-to-well contamination during vortexing and evaporation during PCR.

- e. Vortex the plate at the highest setting speed for 10–30 seconds with medium pressure. Move the plate around to ensure equal contact on the vortex mixer platform.

**IMPORTANT!** Vortex for 10–30 seconds to ensure proper mixing. Failure to do so might result in false classification of samples.

- f. Centrifuge the reaction plate for 1–2 minutes at  $\geq 650 \times g$  ( $\geq 650$  RCF) to remove bubbles and to collect the liquid at the bottom of the reaction plate.

Table 5 Reaction plate volumes

Component	Volume per reaction		
	RNA Sample reaction	Positive Control reaction	Negative Control reaction
Reaction Mix	15.0 μL	15.0 μL	15.0 μL
Purified sample RNA (from RNA extraction)	5.0 μL	—	—
Positive Control (diluted TaqPath™ COVID-19 Control from step 3)	—	2.0 μL	—
Nuclease-free Water	—	3.0 μL	—
Purified Negative Control (from RNA extraction)	—	—	5.0 μL
<b>Total volume</b>	<b>20.0 μL</b>	<b>20.0 μL</b>	<b>20.0 μL</b>

# Perform RT-PCR using the Applied Biosystems™ 7500 Fast Dx Real-Time PCR Instrument

■ Dye calibration for the 7500 Real-Time PCR Instrument series .....	29
■ Transfer the template (SDT) file for the 7500 Fast Dx Real-Time PCR Instrument .....	29
■ Set up and run the 7500 Fast Dx Real-Time PCR Instrument .....	30

## Dye calibration for the 7500 Real-Time PCR Instrument series

A maintained instrument will be calibrated for many dyes. In addition to those dyes, the instrument operator must calibrate the instrument for ABY™ dye and JUN™ dye that are used with this kit. For all other assays, refer to the standard calibration process.

## Transfer the template (SDT) file for the 7500 Fast Dx Real-Time PCR Instrument

The template (SDT) file contains the settings for the instrument run. It is installed on the computer with Applied Biosystems™ COVID-19 Interpretive Software, and must be transferred via a USB drive or other method to the computer on which SDS Software v1.4.1 is installed.

---

**IMPORTANT!** Be careful to select the appropriate template file for your instrument. Failure to do so can cause errors in the analysis.

---

1. On the computer with Applied Biosystems™ COVID-19 Interpretive Software, navigate to the following directory (where <...> is the installation directory):  
`<...>\Applied Biosystems\COVID-19 Interpretive Software\Client\docs\User Documents`
2. Select the SDT file:  
`TaqPath COVID-19 Kit Template 7500fastDx sds1_4_1 v1-4.sdt`
3. Transfer the SDT file to the computer with SDS Software v1.4.1, using a USB drive or other method.

## Set up and run the 7500 Fast Dx Real-Time PCR Instrument

For more information about the instrument, see the documents listed in “Related documentation” on page 67.

1. Using SDS Software v1.4.1, open the SDT file that you transferred in “Transfer the template (SDT) file for the 7500 Fast Dx Real-Time PCR Instrument” on page 29.

---

**IMPORTANT!** Be careful to select the appropriate template file for your instrument. Failure to do so can cause errors in the analysis.

---

2. Confirm the run settings in the template and adjust as necessary.

- **Assay: Standard Curve (Absolute Quantitation)**
- **Run mode: Standard 7500**
- **Passive reference: None**
- **Sample volume: 25 µL**

---

**IMPORTANT!** The passive reference must be set to **None**.

---

3. Confirm that the reporter dye and the detector pairs are correct in the **Detector Manager** in the **Tools** menu.

Reporter dye	Detector
FAM	ORF1ab
VIC	N gene
ABY	S gene
JUN	MS2

4. Confirm that the targets above are assigned to each well in the plate layout.
5. Confirm the labeling of the control wells.
  - The template has one positive control and one negative control assigned to wells for reference.
  - Move the control well assignments by copying the existing control wells and pasting them according to their location on the physical plate.
6. For wells with a positive control, confirm that **Task** is set to **Standard**.
7. For wells with a negative control, confirm that **Task** is set to **NTC**.
8. Edit the plate layout to assign a unique sample name to each well in the physical plate. For wells with a patient sample, confirm that **Task** is set to **Unknown** for all detectors.

---

**Note:** Wells that do not have a sample name will not be analyzed by the software.

---

9. Confirm the thermal protocol.

Step	Temperature	Time	Number of cycles
UNG incubation	25°C	2 minutes	1
Reverse transcription	53°C	10 minutes	1
Activation	95°C	2 minutes	1
Denaturation	95°C	3 seconds	40
Anneal / extension	60°C	30 seconds	

10. Click **Save As**, enter a file name, then click **Save**.
11. Reopen the file to connect the computer to the instrument, load the plate, then start the run on the real-time PCR instrument.
12. After the instrument run is complete, open the SDS file in SDS Software v1.4.1. Analyze, then save the file.

# Perform RT-PCR using the Applied Biosystems™ 7500 Fast Real-Time PCR Instrument

■ Dye calibration for the 7500 Real-Time PCR Instrument series .....	32
■ Transfer the template (SDT or EDT) file for the 7500 Fast Real-Time PCR Instrument .....	33
■ Set up and run the 7500 Fast Real-Time PCR Instrument (SDS Software v1.5.1) .....	33
■ Set up and run the 7500 Fast Real-Time PCR Instrument (7500 Software v2.3) .....	35

## Dye calibration for the 7500 Real-Time PCR Instrument series

A maintained instrument will be calibrated for many dyes. In addition to those dyes, the instrument operator must calibrate the instrument for ABY™ dye and JUN™ dye that are used with this kit. For all other assays, refer to the standard calibration process.

## Transfer the template (SDT or EDT) file for the 7500 Fast Real-Time PCR Instrument

The template (SDT or EDT) file contains the settings for the instrument run. It is installed on the computer with Applied Biosystems™ COVID-19 Interpretive Software, and must be transferred via a USB drive or other method to the computer on which instrument data collection software is installed.

---

**IMPORTANT!** Be careful to select the appropriate template file for your instrument and data collection software. Failure to do so can cause errors in the analysis.

---

1. On the computer with Applied Biosystems™ COVID-19 Interpretive Software, navigate to the following directory (where <...> is the installation directory):

<...>\Applied Biosystems\COVID-19 Interpretive Software\Client\docs\User Documents

2. Select the correct SDT or EDT file for the version of the instrument software that you are using:

Data collection software version	Template file
SDS Software v1.5.1	TaqPath COVID-19 Kit Template 7500fast sds1_5_1 v1-4.sdt
7500 Software v2.3	TaqPath COVID-19 Kit Template 7500fast sds2_3 v1-2.edt

3. Transfer the appropriate SDT or EDT file to the computer with your data collection software, using a USB drive or other method.

## Set up and run the 7500 Fast Real-Time PCR Instrument (SDS Software v1.5.1)

This procedure is specific for the 7500 Fast Real-Time PCR Instrument using SDS Software v1.5.1.

For more information about the instrument, see the documents listed in “Related documentation” on page 67.

1. Using SDS Software v1.5.1, access the appropriate template file.
  - a. Create a new experiment.
  - b. In the **Template** field, browse to, then open the SDT file that you transferred in “Transfer the template (SDT or EDT) file for the 7500 Fast Real-Time PCR Instrument” on page 33.

---

**IMPORTANT!** Be careful to select the appropriate template file for your instrument and data collection software. Failure to do so can cause errors in the analysis.

---

2. Confirm the run settings in the template and adjust as necessary.

- **Assay: Standard Curve (Absolute Quantitation)**
- **Run mode: Standard 7500**
- **Passive reference: None**
- **Sample volume: 25 µL**

---

**IMPORTANT!** The passive reference must be set to **None**.

---

3. Confirm that the reporter dye and the detector pairs are correct in the **Detector Manager** in the **Tools** menu.

Reporter dye	Detector
FAM	ORF1ab
VIC	N gene
ABY	S gene
JUN	MS2

4. Confirm that the targets above are assigned to each well in the plate layout.
5. Confirm the labeling of the control wells.
- The template has one positive control (PC) and one negative control (NC) assigned to wells for reference.
  - Move the control well assignments by copying the existing control wells and pasting them according to their location on the physical plate.
6. For wells with a positive control, confirm that **Task** is set to **Standard**.
7. For wells with a negative control, confirm that **Task** is set to **NTC**.
8. Edit the plate layout to assign a unique sample name to each well in the physical plate.  
For wells with a patient sample, ensure that **Task** is set to **Unknown** for all detectors.

---

**Note:** Wells that do not have a sample name will not be analyzed by the software.

---

9. Confirm the thermal protocol.

Step	Temperature	Time	Number of cycles
UNG incubation	25°C	2 minutes	1
Reverse transcription	53°C	10 minutes	1
Activation	95°C	2 minutes	1
Denaturation	95°C	3 seconds	40
Anneal / extension	60°C	30 seconds	

10. Click **Save As**, enter a file name, then click **Save**.

11. Reopen the file, load the plate, then start the run on the instrument.
12. After the instrument run is complete, open the SDS file in SDS Software v1.5.1. Analyze, then save the file.

## Set up and run the 7500 Fast Real-Time PCR Instrument (7500 Software v2.3)

This procedure is specific for the 7500 Fast Real-Time PCR Instrument using SDS Software v2.3.

For more information about the instrument, see the documents listed in “Related documentation” on page 67.

1. In the 7500 Software v2.3 home page, click **Template**.
2. Browse to, then open the EDT file that you transferred in “Transfer the template (SDT or EDT) file for the 7500 Fast Real-Time PCR Instrument” on page 33.

---

**IMPORTANT!** Be careful to select the appropriate template file for your instrument and data collection software. Failure to do so can cause errors in the analysis.

---

3. In the **Experiment Properties** window, enter or confirm the following information:
  - **Experiment name:** Enter a unique name
  - **Instrument type:** 7500 Fast (96 wells)
  - **Type of experiment:** Quantitation - Standard Curve
  - **Reagents:** TaqMan™
  - **Ramp Speed:** Standard
4. In the **Plate Setup** window, in the **Define Targets and Samples** tab and the **Define Targets** pane, confirm that the targets, reporter dyes, and quenchers are listed correctly.

Target	Reporter dye	Quencher
ORF1ab	FAM	None
N gene	VIC	None
S gene	ABY	None
MS2	JUN	None

5. In the **Plate Setup** window, in the **Define Targets and Samples** tab and the **Define Samples** pane, confirm the labeling of the control wells.
  - The template has one positive control (PC) and one negative control (NC) assigned to wells for reference.
  - Move the control well assignments by copying the existing control wells and pasting them according to their location on the physical plate.
  - To include additional controls, select **Add New Sample**.

6. Select **Add New Sample** to assign a unique sample name to each well in the physical plate.
7. In the **Plate Setup** window, in the **Assign Targets and Samples** tab, confirm that four targets are assigned to each well in the plate layout. To designate a target or sample to a well, select the well, then check the **Assign** checkbox.

---

**Note:** Wells that do not have a sample name will not be analyzed by the software.

---

8. For wells with a positive control, confirm that **Task** is set to **S** for Standard.
9. For wells with a negative control, confirm that **Task** is set to **N** for Negative.
10. For wells with a patient sample, confirm that **Task** is set to **U** for Unknown.
11. Confirm that the **Passive Reference** is set to **None**.
12. In the **Run Method** window, confirm that the **Reaction Volume Per Well** is 25 µL, then confirm the thermal protocol.

Step	Temperature	Time	Number of cycles
UNG incubation	25°C	2 minutes	1
Reverse transcription	53°C	10 minutes	1
Activation	95°C	2 minutes	1
Denaturation	95°C	3 seconds	40
Anneal / extension	60°C	30 seconds	

13. Select **Start Run**, enter a file name, then click **Save**.
14. After the instrument run is complete, click **Analyze**, then save the file.

# Perform RT-PCR using the Applied Biosystems™ 7500 Real-Time PCR Instrument

■ Dye calibration for the 7500 Real-Time PCR Instrument series .....	37
■ Transfer the template (EDT) file for the 7500 Real-Time PCR Instrument .....	38
■ Set up and run the 7500 Real-Time PCR Instrument .....	38

## Dye calibration for the 7500 Real-Time PCR Instrument series

A maintained instrument will be calibrated for many dyes. In addition to those dyes, the instrument operator must calibrate the instrument for ABY™ dye and JUN™ dye that are used with this kit. For all other assays, refer to the standard calibration process.

---

**IMPORTANT!** Use only the calibration plates listed in “Required materials not supplied” on page 9.

---

## Transfer the template (EDT) file for the 7500 Real-Time PCR Instrument

The template (EDT) file contains the settings for the instrument run. It is installed on the computer with Applied Biosystems™ COVID-19 Interpretive Software, and must be transferred via a USB drive or other method to the computer on which 7500 Software v2.3 is installed.

---

**IMPORTANT!** Be careful to select the appropriate template file for your instrument. Failure to do so can cause errors in the analysis.

---

1. On the computer with Applied Biosystems™ COVID-19 Interpretive Software, navigate to the following directory (where <...> is the installation directory):

```
<...>\Applied Biosystems\COVID-19 Interpretive Software\Client\docs\User Documents
```

2. Select the EDT file:

```
TaqPath COVID-19 Kit Template 7500std sds2_3 v1-2.edt
```

3. Transfer the EDT file to the computer with your data collection software, using a USB drive or other method.

## Set up and run the 7500 Real-Time PCR Instrument

This procedure is specific for the 7500 Real-Time PCR Instrument using 7500 Software v2.3.

For more information about the instrument, see the documents listed in “Related documentation” on page 67.

1. In the 7500 Software v2.3 home page, click **Template**.
2. Browse to, then open the template file that you transferred in “Transfer the template (EDT) file for the 7500 Real-Time PCR Instrument” on page 38.

---

**IMPORTANT!** Be careful to select the appropriate template file for your instrument. Failure to do so can cause errors in the analysis.

---

3. In the **Experiment Properties** window, enter or confirm the following information:
  - **Experiment name:** Enter a unique name
  - **Instrument type:** 7500 (96 wells)
  - **Type of experiment:** Quantitation - Standard Curve
  - **Reagents:** TaqMan™
  - **Ramp Speed:** Standard

4. In the **Plate Setup** window, in the **Define Targets and Samples** tab and the **Define Targets** pane, confirm that the targets, reporter dyes, and quenchers are listed correctly.

Target	Reporter dye	Quencher
ORF1ab	FAM	None
N gene	VIC	None
S gene	ABY	None
MS2	JUN	None

5. In the **Plate Setup** window, in the **Define Targets and Samples** tab and the **Define Samples** pane, confirm the labeling of the control wells.
- The template has one positive control (PC) and one negative control (NC) assigned to wells for reference.
  - Move the control well assignments by copying the existing control wells and pasting them according to their location on the physical plate.

To include additional controls, select **Add New Sample**.

6. Select **Add New Sample** to assign a unique sample name to each well in the physical plate.
7. In the **Plate Setup** window, in the **Assign Targets and Samples** tab, confirm that four targets are assigned to each well in the plate layout. To designate a target or sample to a well, select the well, then check the **Assign** checkbox.

---

**Note:** Wells that do not have a sample name will not be analyzed by the software.

---

8. For wells with a positive control, confirm that **Task** is set to **S** for Standard.
9. For wells with a negative control, confirm that **Task** is set to **N** for Negative.
10. For wells with a patient sample, confirm that **Task** is set to **U** for Unknown.
11. Confirm that the **Passive Reference** is set to **None**.
12. In the **Run Method** window, confirm that the **Reaction Volume Per Well** is 25 µL, then confirm the thermal protocol.

Step	Temperature	Time	Number of cycles
UNG incubation	25°C	2 minutes	1
Reverse transcription	53°C	10 minutes	1
Activation	95°C	2 minutes	1
Denaturation	95°C	3 seconds	40
Anneal / extension	60°C	30 seconds	

13. Select **Start Run**, enter a file name, then click **Save**.
14. After the instrument run is complete, click **Analyze**, then save the file.

# Perform RT-PCR using the Applied Biosystems™ QuantStudio™ 5 Real-Time PCR Instrument

■ Dye calibration for the QuantStudio™ 5 Real-Time PCR Instrument . . . . .	40
■ Transfer the template (EDT) file for the QuantStudio™ 5 Real-Time PCR Instrument . . . . .	41
■ Set up and run the QuantStudio™ 5 Real-Time PCR Instrument (96-well plates) . . . . .	41
■ Set up and run the QuantStudio™ 5 Real-Time PCR Instrument (384-well plates) . . . . .	43

## Dye calibration for the QuantStudio™ 5 Real-Time PCR Instrument

A maintained instrument will be calibrated for all dyes that are used with this kit. Ensure that the calibrations for FAM™ dye, VIC™ dye, ABY™ dye, and JUN™ dye are current. If calibration is required, refer to the standard calibration process in the instrument user guide.

## Transfer the template (EDT) file for the QuantStudio™ 5 Real-Time PCR Instrument

The template (EDT) file contains the settings for the instrument run. It is installed on the computer with Applied Biosystems™ COVID-19 Interpretive Software, and must be transferred via a USB drive or other method to the computer on which QuantStudio™ Design and Analysis Desktop Software v1.5.1 is installed.

**IMPORTANT!** Be careful to select the appropriate template file for your instrument and block type. Failure to do so can cause errors in the analysis.

1. On the computer with Applied Biosystems™ COVID-19 Interpretive Software, navigate to the following directory (where <...> is the installation directory):

<...>\Applied Biosystems\COVID-19 Interpretive Software\Client\docs\User Documents

2. Select the correct EDT file for your instrument.

Instrument	Template file
QuantStudio™ 5 Real-Time PCR Instrument (96-well, 0.1-mL block)	TaqPath COVID-19 Kit Template QS5 0_1ml_da1_5_1 v2-2.edt
QuantStudio™ 5 Real-Time PCR Instrument (96-well, 0.2-mL block)	TaqPath COVID-19 Kit Template QS5 0_2ml_da1_5_1 v2-2.edt
QuantStudio™ 5 Real-Time PCR Instrument (384-well block)	TaqPath COVID-19 Kit Template QS5 384_da1_5_1 v2-2.edt

3. Transfer the EDT file to the computer with QuantStudio™ Design and Analysis Desktop Software v1.5.1, using a USB drive or other method.

## Set up and run the QuantStudio™ 5 Real-Time PCR Instrument (96-well plates)

For more information about the instrument, see the documents listed in “Related documentation” on page 67.

1. In the QuantStudio™ Design and Analysis Desktop Software v1.5.1, in the **New Experiment** box, select **Create New Experiment ▶ Template**.
2. Browse to, then open the EDT file that you transferred in “Transfer the template (EDT) file for the QuantStudio™ 5 Real-Time PCR Instrument” on page 41.

**IMPORTANT!** Be careful to select the appropriate template file for your instrument and block type. Failure to do so can cause errors in the analysis.

3. In the **Properties** tab, enter or confirm the following.

- **Name:** Enter a unique name
- **Instrument type:** QuantStudio™ 5 System
- **Block type:** 96-Well 0.2-mL Block or 96-Well 0.1-mL Block
- **Experiment type:** Standard Curve
- **Chemistry:** TaqMan™ Reagents
- **Run Mode:** Standard

4. In the **Method** tab, confirm that the **Volume** is 25 µL, then confirm the thermal protocol.

Step	Temperature	Time	Number of cycles
UNG incubation	25°C	2 minutes	1
Reverse transcription	53°C	10 minutes	1
Activation	95°C	2 minutes	1
Denaturation	95°C	3 seconds	40
Anneal / extension	60°C	30 seconds	

5. In the **Plate** tab, click **Quick Setup**.

6. In the **Plate Attributes** pane, confirm that the **Passive Reference** is set to **None**.

7. In the **Plate** tab, click **Advanced Setup**.

8. In the **Targets** table, confirm that the reporter dye and the target pairs are correct.

Reporter dye	Detector	Quencher
FAM	ORF1ab	None
VIC	N gene	None
ABY	S gene	None
JUN	MS2	None

9. Confirm that the targets above are assigned to each well in the plate layout.

10. In the plate layout pane, confirm the labeling of the control wells.

- The template has one positive control (PC) and one negative control (NC) assigned to wells for reference.
- Move the control well assignments by copying the existing control wells and pasting them according to their location on the physical plate.

11. For all targets in the positive control well, confirm that **Task** is set to **S (Standard)**.

12. For all targets in the negative control well, confirm that **Task** is set to **N (Negative Control)**.

13. In the **Samples** table, click **Add** to define the sample names. Create a unique sample name for each well in the physical plate.

14. To assign a sample to a well, select the well in the plate layout, then select the sample from the **Samples** table.

For all targets in the patient sample wells, confirm that **Task** is set to **U (Unknown)**.

---

**Note:** Wells that do not have a sample name will not be analyzed by the software.

---

15. In the **Run** tab, click **Start Run**, then select your instrument from the drop-down list.
16. Enter a file name in the dialog box that prompts you to save the run file, then save the file.

## Set up and run the QuantStudio™ 5 Real-Time PCR Instrument (384-well plates)

For more information about the instrument, see the documents listed in “Related documentation” on page 67.

1. In the QuantStudio™ Design and Analysis Desktop Software v1.5.1, in the **New Experiment** box, select **Create New Experiment ▶ Template**.
2. Browse to, then open the EDT file that you transferred in “Transfer the template (EDT) file for the QuantStudio™ 5 Real-Time PCR Instrument” on page 41.

---

**IMPORTANT!** Be careful to select the appropriate template file for your instrument and block type. Failure to do so can cause errors in the analysis.

---

3. In the **Properties** tab, enter or confirm the following.

- **Name:** Enter a name
- **Instrument type:** QuantStudio™ 5 System
- **Block type:** 384-well Block
- **Experiment type:** Standard Curve
- **Chemistry:** TaqMan™ Reagents
- **Run Mode:** Standard

4. In the **Method** tab, confirm that the **Volume** is 20 µL, then confirm the thermal protocol.

Step	Temperature	Time	Number of cycles
UNG incubation	25°C	2 minutes	1
Reverse transcription	53°C	10 minutes	1
Activation	95°C	2 minutes	1
Denaturation	95°C	3 seconds	40
Anneal / extension	60°C	30 seconds	

5. In the **Plate** tab, click **Quick Setup**.
6. In the **Plate Attributes** pane, confirm that the **Passive Reference** is set to **None**.

7. In the **Plate** tab, click **Advanced Setup**.

8. In the **Targets** table, confirm that the reporter dye and the target pairs are correct.

Reporter dye	Detector	Quencher
FAM	ORF1ab	None
VIC	N gene	None
ABY	S gene	None
JUN	MS2	None

9. Confirm that the targets above are assigned to each well in the plate layout.

10. In the plate layout pane, confirm the labeling of the control wells.

- The template has one positive control (PC) and one negative control (NC) assigned to wells for reference.
- Move the control well assignments by copying the existing control wells and pasting them according to their location on the physical plate.

---

**Note:** If samples from more than one RNA extraction run are being included on the same RT-PCR 384-well plate, then there must be one negative control for each RNA extraction run included on the RT-PCR 384-well plate (see “Guidelines for RT-PCR” on page 20). Label each negative control with a unique name, for example, **NC1**, **NC2**, **NC3**, and **NC4**.

---

11. For all targets in the positive control well, confirm that **Task** is set to **S (Standard)**.

12. For all targets in the negative control well, confirm that **Task** is set to **N (Negative Control)**.

13. In the **Samples** table, click **Add** to define the names. Create a unique sample name for each well in the physical plate.

The template has one positive control and one negative control assigned to wells for reference. If additional control wells are required, each control well must have a unique name.

14. To assign a sample to a well, select the well in the plate layout, then select the sample from the **Samples** table.

For all targets in the patient sample wells, confirm that **Task** is set to **U (Unknown)**.

---

**Note:** Wells that do not have a sample name will not be analyzed by the software.

---

15. In the **Run** tab, click **Start Run**, then select your instrument from the drop-down list.

16. Enter a file name in the dialog box that prompts you to save the run file, then save the file.

# Perform RT-PCR using the Applied Biosystems™ QuantStudio™ 5 Dx Real-Time PCR Instrument

■ Dye calibration for the QuantStudio™ 5 Dx Real-Time PCR Instrument .....	45
■ QuantStudio™ 5 Dx TD Software .....	46
■ QuantStudio™ 5 Dx IVD Software .....	49

## Dye calibration for the QuantStudio™ 5 Dx Real-Time PCR Instrument

A maintained instrument will be calibrated for all dyes that are used with this kit. Ensure that the calibrations for FAM™ dye, VIC™ dye, ABY™ dye, and JUN™ dye are current. For all other assays, refer to the standard calibration process.

## QuantStudio™ 5 Dx TD Software

Follow the procedures in this section if you are using the QuantStudio™ 5 Dx TD Software.

If you are using QuantStudio™ 5 Dx IVD Software, see “QuantStudio™ 5 Dx IVD Software” on page 49.

### Transfer the template (EDT) file for the QuantStudio™ 5 Dx Real-Time PCR Instrument (QuantStudio™ 5 Dx TD Software)

The template (EDT) file contains the settings for the instrument run. It is installed on the computer with Applied Biosystems™ COVID-19 Interpretive Software, and must be transferred via a USB drive or other method to the computer on which QuantStudio™ 5 Dx TD Software is installed.

1. On the computer with Applied Biosystems™ COVID-19 Interpretive Software, navigate to the following directory (where <...> is the installation directory):

```
<...>\Applied Biosystems\COVID-19 Interpretive Software\Client\docs\User Documents
```

2. Select the EDT file:

```
TaqPath COVID-19 Kit Template QS5 Dx_da1_0_2 v2-2.edt
```

3. Transfer the EDT file to the computer with the QuantStudio™ 5 Dx TD Software, using a USB drive or other method.

---

**IMPORTANT!** Be careful to select the appropriate template file for the instrument and software version that you are using. Failure to do so can cause errors in the analysis.

---

### Set up and run the QuantStudio™ 5 Dx Real-Time PCR Instrument (QuantStudio™ 5 Dx TD Software)

For more information about the QuantStudio™ 5 Dx Real-Time PCR Instrument, see the documents listed in “Related documentation” on page 67.

1. In the QuantStudio™ 5 Dx TD Software, in the **New Experiment** box, select **Create New Experiment ▶ Template**.
2. Browse to, then open the EDT file that you transferred in “Transfer the template (EDT) file for the QuantStudio™ 5 Dx Real-Time PCR Instrument (QuantStudio™ 5 Dx TD Software)” on page 46.

---

**IMPORTANT!** Be careful to select the appropriate template file for the instrument and block type. Failure to do so can cause errors in the analysis.

---

3. In the **Properties** tab, enter or confirm the following.

- **Name:** Enter a unique name
- **Instrument type:** QuantStudio™ 5 Dx System
- **Block type:** 96-Well 0.2-mL Block
- **Experiment type:** Standard Curve
- **Chemistry:** TaqMan™ Reagents
- **Run Mode:** Standard

4. In the **Method** tab, confirm that **Volume** is 25 µL, then confirm the thermal protocol.

Step	Temperature	Time	Number of cycles
UNG incubation	25°C	2 minutes	1
Reverse transcription	53°C	10 minutes	1
Activation	95°C	2 minutes	1
Denaturation	95°C	3 seconds	40
Anneal / extension	60°C	30 seconds	

5. In the **Plate** tab, click **Quick Setup**.
6. In the **Plate Attributes** pane, confirm that the **Passive Reference** is set to **None**.
7. In the **Plate** tab, click **Advanced Setup**.
8. In the **Targets** table, confirm that the reporter dye and the target pairs are correct.

Reporter dye	Detector	Quencher
FAM	ORF1ab	None
VIC	N gene	None
ABY	S gene	None
JUN	MS2	None

9. Confirm that the targets above are assigned to each well in the plate layout.
10. In the plate layout pane, confirm the labeling of the control wells.
- The template has one positive control (PC) and one negative control (NC) assigned to wells for reference.
  - Move the control well assignments by copying the existing control wells and pasting them according to their location on the physical plate.
11. For all targets in the positive control well, confirm that **Task** is set to **S (Standard)**.
12. For all targets in the negative control well, confirm that **Task** is set to **N (Negative Control)**.
13. In the **Samples** table, click **Add** to define the sample names. Create a unique sample name for each well in the physical plate.
14. To assign a sample to a well, select the well in the plate layout, then select the sample from the **Samples** table.
- For all targets in the patient sample wells, confirm that **Task** is set to **U (Unknown)**.

---

**Note:** Wells that do not have a sample name will not be analyzed by the software.

---

15. In the **Run** tab, click **Start Run**, then select your instrument from the drop-down list.
16. Enter a file name in the dialog box that prompts you to save the run file, then save the file.

## QuantStudio™ 5 Dx IVD Software

Follow the procedures in this section if you are using the QuantStudio™ 5 Dx IVD Software.

If you are using QuantStudio™ 5 Dx TD Software, see “QuantStudio™ 5 Dx TD Software” on page 46.

### Transfer the template (EDT) file for the QuantStudio™ 5 Dx Real-Time PCR Instrument (QuantStudio™ 5 Dx IVD Software)

The template (EDT) file contains the settings for the instrument run. It is installed on the computer with Applied Biosystems™ COVID-19 Interpretive Software, and must be transferred via a USB drive or other method to the computer on which QuantStudio™ 5 Dx IVD Software is installed.

1. On the computer with Applied Biosystems™ COVID-19 Interpretive Software, navigate to the following directory (where <...> is the installation directory):

```
<...>\Applied Biosystems\COVID-19 Interpretive Software\Client\docs\User Documents
```

2. Select the EDT file:

```
TaqPath COVID-19 Kit Template QS5 Dx_da1_0_2_CE_IVD v2-2.edt
```

3. Transfer the EDT file to the computer with the QuantStudio™ 5 Dx IVD Software, using a USB drive or other method.

---

**IMPORTANT!** Be careful to select the appropriate template file for the instrument and software version that you are using. Failure to do so can cause errors in the analysis.

---

The template must be installed in the QuantStudio™ 5 Dx IVD Software before the run is started. See “Install the template file in the QuantStudio™ 5 Dx IVD Software” on page 49.

### Install the template file in the QuantStudio™ 5 Dx IVD Software

To install a template, the logged-in user must have the SAE permission of **Manage Installed Templates**. For information about SAE permissions, see the *QuantStudio™ 5 Dx TD Software User Guide* (Pub. No. 100049554).

1. Open the QuantStudio™ 5 Dx IVD Software.
2. When prompted, log in with a user account with the appropriate permissions.
3. In the menu bar, select **Tools ▶ Template Menu**.
4. Click **Install**, then click **Yes** to confirm the **Template Installation Agreement**.
5. Navigate template file that was transferred, then click **Open**.  
The template is now installed and accessible in the **Template Menu**.
6. Click **Close**.

## Set up and run the QuantStudio™ 5 Dx Real-Time PCR Instrument (QuantStudio™ 5 Dx IVD Software)

For more information about the QuantStudio™ 5 Dx Real-Time PCR Instrument, see the documents listed in “Related documentation” on page 67.

---

**IMPORTANT!** When the QuantStudio™ 5 Dx IVD Software is used, the positive control must be placed in well A12 and the negative control must be placed in well A1 on the physical plate.

---

1. In the QuantStudio™ 5 Dx IVD Software, in the **New Experiment** box, select **Experiment Setup**.
2. Select the template file that you installed in “Install the template file in the QuantStudio™ 5 Dx IVD Software” on page 49.
3. Click **Create New Experiment**.
4. In **Samples** table, click **Add** to define the sample names. Create a unique sample name for each well in the physical plate.
5. To assign a sample to a well, select the well in the plate layout, then select sample from the **Samples** table.

---

**Note:** Wells that do not have a sample name will not be analyzed by the software.

---

6. In the **Run** tab, click **Start Run**, then select your instrument from the drop-down list. The **Enter Reason for Change** dialog box will appear.
7. In the **Enter Reason for Change** dialog box, click **OK**.
8. Enter a file name in the dialog box that prompts you to save the run file, then save the file.

# Perform RT-PCR using the Applied Biosystems™ QuantStudio™ 7 Flex Real-Time PCR Instrument (384-well block)

■ Dye calibration for the QuantStudio™ 7 Flex Real-Time PCR Instrument . . . . .	51
■ Transfer the template (EDT) file for the QuantStudio™ 7 Flex Real-Time PCR Instrument (384-well block) . . . . .	52
■ Set up and run the QuantStudio™ 7 Flex Real-Time PCR Instrument (384-well block) . . . . .	52

## Dye calibration for the QuantStudio™ 7 Flex Real-Time PCR Instrument

A maintained instrument will be calibrated for many dyes. In addition to those dyes, the instrument operator must calibrate the instrument for ABY™ dye and JUN™ dye that are used with this kit. For all other assays, refer to the standard calibration process.

## Transfer the template (EDT) file for the QuantStudio™ 7 Flex Real-Time PCR Instrument (384-well block)

The template (EDT) file contains the settings for the instrument run. It is installed on the computer with Applied Biosystems™ COVID-19 Interpretive Software, and must be transferred via a USB drive or other method to the computer on which the QuantStudio™ Real-Time PCR Software v1.3 is installed.

---

**IMPORTANT!** Be careful to select the appropriate template file for your instrument and block type. Failure to do so can cause errors in the analysis.

---

1. On the computer with Applied Biosystems™ COVID-19 Interpretive Software, navigate to the following directory (where <...> is the installation directory):

```
<...>\Applied Biosystems\COVID-19 Interpretive Software\Client\docs\User Documents
```

2. Select the EDT file:

```
TaqPath COVID-19 Kit Template QS7 384 1_3 v2-2.edt
```

3. Transfer the EDT file to the computer with QuantStudio™ Real-Time PCR Software v1.3, using a USB drive or other method.

## Set up and run the QuantStudio™ 7 Flex Real-Time PCR Instrument (384-well block)

For more information about the instrument, see the documents listed in “Related documentation” on page 67.

1. In the QuantStudio™ Real-Time PCR Software v1.3 home screen, click **Template**.
2. Browse to, then open the EDT file that you transferred in “Transfer the template (EDT) file for the QuantStudio™ 7 Flex Real-Time PCR Instrument (384-well block)” on page 52.

---

**IMPORTANT!** Be careful to select the appropriate template file for your instrument and block type. Failure to do so can cause errors in the analysis.

---

3. In the **Experiment Properties** tab, enter or confirm the following.
  - **Experiment Name:** Enter a unique name
  - **Instrument type:** QuantStudio™ 7 Flex System
  - **Block:** 384-well
  - **Type of Experiment:** Standard Curve
  - **Reagents:** TaqMan™
  - **Properties:** Standard

4. In the **Define** tab, in the **Targets** pane, confirm that the targets, reporter dyes, and quenchers are listed correctly.

Target	Reporter dye	Quencher
ORF1ab	FAM	None
N gene	VIC	None
S gene	ABY	None
MS2	JUN	None

5. In the **Define** tab, in the **Samples** pane, define a unique sample name for each well in the physical plate.

---

**Note:** Wells that do not have a sample name will not be analyzed by the software.

---

6. In the **Define** tab, confirm that the **Passive Reference** is set to **None**.
7. In the **Assign** tab, confirm that four targets are assigned to each well in the plate layout.  
To assign a target to a well, select the well, then check the **Assign** checkbox.
8. In the **Assign** tab, in the **Samples** pane, confirm the labeling of the control wells.
  - The template has one positive control (PC) and one negative control (NC) assigned to wells for reference.
  - Move the control well assignments by copying the existing control wells and pasting them according to their location on the physical plate.

---

**Note:** If the 384-well plate includes additional negative controls (from combined extraction runs), define each additional negative control and assign it to a well to match the physical plate. Give each negative control a unique name (for example, **NC1**, **NC2**, **NC3**, and **NC4**).

---

9. In the **Assign** tab, confirm the **Task** assignments.
  - For wells with a Positive Control (PC), confirm that the **Task** is set to **S** for Standard for all of the targets.
  - For wells with a Negative Control (NC), confirm that the **Task** is set to **N** for Negative for all of the targets.
  - For the wells with a patient sample, confirm that the **Task** is set to **U** for Unknown for all of the targets.
10. In the **Assign** tab, assign a sample name to each well to match the physical plate.  
To assign a sample to a well, select the well, then check the **Assign** checkbox.

11. In the **Run Method** tab, confirm that the **Reaction Volume Per Well** is 20 µL, then confirm the thermal protocol.

Step	Temperature	Time	Number of cycles
UNG incubation	25°C	2 minutes	1
Reverse transcription	53°C	10 minutes	1
Activation	95°C	2 minutes	1
Denaturation	95°C	3 seconds	40
Anneal / extension	60°C	30 seconds	

12. In the **Run** tab, click **Start Run**, then select your instrument from the drop-down list.
13. Enter a file name in the dialog box that prompts you to save the run file, then save the file.

# Analysis and results

## Obtain the Applied Biosystems™ COVID-19 Interpretive Software

To perform data analysis and results interpretation, you must use the Applied Biosystems™ COVID-19 Interpretive Software.

To obtain the software, contact your local instrument service team. Go to <https://www.thermofisher.com/contactus>.

## Analyze the data

For detailed instructions about using the software, click the **Help** menu in the COVID-19 Interpretive Software.

1. Using a USB drive or other method, transfer the SDS or EDS files from the computer with the data collection software to the computer with the COVID-19 Interpretive Software.
2. Sign in to the COVID-19 Interpretive Software.
3. In the **Home** screen, click the **Import Samples** button.
4. Select the SDS files or the EDS files to import, then click **Open**.  
After import, the software analyzes the run data, performs Quality Check (QC) analysis, and calculates the interpretive results for each sample and control.
5. In the **Home** screen, under **Batches**, click the **<Batch ID>** link for the batch of interest to display the **Batch Details** screen.
  - In the **Batch Details** screen, view the status and result for each sample in the **Sample Information** tab.
  - **Optional:** To identify signs of inadequate vortexing or centrifugation, select the **Amplification Plot** tab to view the amplification curve for each sample.
  - **Optional:** To enter a comment for the sample, select the **Amplification Plot** tab, select the **Comment** tab, enter a comment, then click **Add**.  
The comment will appear in all exports (CSV files) and reports (PDF files). It will include the user name and date.

6. To generate a Batch Export (CSV file), return to the **Home** screen, select the checkbox for the batch, then click **Export Batch** at the top of the **Home** screen. Click **Open folder location** in the dialog box to locate and open the CSV file.
7. To generate a Batch Report (PDF file), select the checkbox for the batch, then click **Report Batch** at the top of the **Home** screen. Click **Open folder location** in the dialog box to locate and open the PDF file.

## Interpretation of the results

Interpretation of the results is performed by the Applied Biosystems™ COVID-19 Interpretive Software. For information about the  $C_t$  values that are used by the software to interpret results, see Appendix A, “Ct cutoff values for assay targets”.

### Quality control and validity of results

A minimum of one Negative Control and one Positive Control must be present for each run for the COVID-19 assay.

Additional Negative Control wells must be run for each extraction that is represented on a real-time RT-PCR plate. All control wells must pass for the real-time RT-PCR plate to be considered valid.

Validation of results is performed automatically by the Applied Biosystems™ COVID-19 Interpretive Software based on performance of the Positive and Negative Controls.

**Table 6 Patient samples**

ORF1ab	N gene	S gene	MS2	Status	Result	Action
NEG	NEG	NEG	NEG	INVALID	NA	Repeat test by re-extracting the original sample and repeating the RT-PCR. If the repeat result remains invalid, consider collecting a new specimen.
NEG	NEG	NEG	POS	VALID	SARS-CoV-2 Not Detected	Report results to the healthcare provider and appropriate public health authorities. Consider testing for other viruses.
Only one SARS-CoV-2 target = POS			POS or NEG	VALID	SARS-CoV-2 Inconclusive	<ol style="list-style-type: none"> <li>1. Repeat test by re-extracting the original sample and repeating the RT-PCR.</li> <li>2. After retesting one time, report results to the healthcare provider and appropriate public health authorities.</li> </ol> <p><b>IMPORTANT!</b> Samples with a result of SARS-CoV-2 Inconclusive shall be retested one time.</p> <p>If the repeat result remains inconclusive, the healthcare provider should conduct additional confirmation testing with a new specimen, if clinically indicated.</p>
Two or more SARS-CoV-2 targets = POS			POS or NEG	VALID	Positive SARS-CoV-2	Report results to the healthcare provider and appropriate public health authorities.

## Performance characteristics

The analytical and clinical performance of the TaqPath™ COVID-19 CE-IVD RT-PCR Kit was evaluated by determining limit of detection (LoD), characterizing the impact of interfering substances and cross-reactivity, as described in the following sections.

### Limit of detection (LoD)

The LoD study established the lowest SARS-CoV-2 viral concentration (Genomic Copy Equivalents or GCE) that can be detected by the TaqPath™ COVID-19 CE-IVD RT-PCR Kit in a particular specimen type at least 95% of the time. Banked Nasopharyngeal swab (NP) and Bronchoalveolar lavage (BAL) samples were obtained from patients in the United States in the years 2015–2019. The NP and BAL samples were pooled, respectively, and spiked with purified SARS-CoV-2 viral genomic RNA at several concentrations, then processed using the kit workflow. A three-phase approach was used to determine the LoD for each specimen type. In phases I and II, the preliminary LoD was established and confirmed in phase III by testing 20 replicates.

Table 7 LoD determination in BAL

Effective Concentration	Replicate	Mean C <sub>t</sub>				Interpretation	% Positive
		ORF1ab	N	S	MS2		
10 GCE/ reaction	1	29.9	29.1	28.5	23.1	Positive	100%
	2	30.1	29.3	29.7	24.0	Positive	
	3	30.0	29.7	29.3	24.0	Positive	
	4	30.3	29.7	29.1	23.8	Positive	
	5	30.2	29.6	29.6	23.7	Positive	
	6	30.3	29.3	29.7	23.5	Positive	
	7	29.9	29.6	32.8	23.4	Positive	
	8	30.2	29.8	29.2	23.8	Positive	
	9	30.1	29.4	28.6	23.8	Positive	
	10	30.1	29.4	29.1	24.0	Positive	
	11	29.8	29.5	29.4	24.3	Positive	
	12	30.1	29.7	29.1	24.6	Positive	
	13	30.7	30.1	28.4	25.1	Positive	

Table 7 LoD determination in BAL (continued)

Effective Concentration	Replicate	Mean C <sub>t</sub>				Interpretation	% Positive
		ORF1ab	N	S	MS2		
10 GCE/ reaction	14	30.4	29.8	29.1	24.8	Positive	100%
	15	30.2	29.8	29.7	24.9	Positive	
	16	30.3	29.8	29.4	24.6	Positive	
	17	30.4	30.0	31.5	24.7	Positive	
	18	30.4	30.1	29.3	24.9	Positive	
	19	30.9	29.7	29.2	25.4	Positive	
	20	30.3	29.9	29.4	25.7	Positive	

Table 8 LoD determination in Nasopharyngeal specimens

Effective Concentration	Replicate	Mean C <sub>t</sub>				Interpretation	% Positive
		ORF1ab	N	S	MS2		
10 GCE/ reaction	1	30.0	28.9	35.7	25.7	Positive	100%
	2	30.6	28.9	33.6	25.8	Positive	
	3	30.2	28.8	32.0	25.8	Positive	
	4	30.4	28.7	34.2	25.7	Positive	
	5	30.5	29.0	31.4	25.8	Positive	
	6	31.0	29.3	36.6	26.0	Positive	
	7	30.3	29.2	31.1	25.8	Positive	
	8	31.1	29.2	31.8	26.5	Positive	
	9	30.5	28.9	33.0	26.2	Positive	
	10	30.3	28.8	34.7	26.8	Positive	
	11	30.5	29.8	38.7	27.4	Positive	
	12	31.6	29.7	35.0	27.6	Positive	
	13	30.7	29.3	36.4	27.4	Positive	
	14	31.6	28.8	31.3	27.2	Positive	
	15	31.0	29.3	36.0	27.0	Positive	
	16	30.5	29.1	35.7	27.0	Positive	
	17	30.7	29.4	34.8	27.4	Positive	

Table 8 LoD determination in Nasopharyngeal specimens (continued)

Effective Concentration	Replicate	Mean C <sub>t</sub>				Interpretation	% Positive
		ORF1ab	N	S	MS2		
10 GCE/ reaction	18	30.7	29.3	34.6	27.5	Positive	100%
	19	31.0	29.3	35.9	28.7	Positive	
	20	30.4	29.2	32.7	28.4	Positive	

Table 9 LoD results

LoD studies were performed by extracting 400 µL of each specimen followed by elution in 50 µL, then 5 µL of eluate was added to the RT-PCR reaction. The LoD of 10 GCE/reaction is calculated from a starting concentration of 250 GCE/mL of specimen.

Specimen type	Limit of Detection (GCE/reaction)
Bronchoalveolar lavage	10 GCE/reaction
Nasopharyngeal swab	10 GCE/reaction

## Reactivity (Inclusivity)

*In silico* analysis was updated on January 04, 2021 using 309,591 sequenced genomes in GISAID and NCBI GenBank databases which include instances of emerging SARS-CoV-2 variants that are under global surveillance. Based on the BLAST analysis, the TaqPath™ COVID-19 CE-IVD RT-PCR Kit maps with 100% homology to 96.2% of known SARS-CoV-2 strains/isolates in GISAID and 99.8% of known strains/isolates in NCBI GenBank. Mapping was deemed successful for a given isolate if at least two of the three target assays (ORF1ab, S gene, and N gene) showed 100% identity.

## Interfering substances

Pooled SARS-CoV-2-negative nasopharyngeal swab and bronchoalveolar lavage specimens were spiked with purified SARS-CoV-2 viral RNA at 3X the Limit of Detection (30 GCE/reaction) and potential interfering substances at the concentrations above. Each substance was tested with triplicate extractions. The results are presented in Table 10 on page 61.

Pooled SARS-CoV-2-negative nasopharyngeal swab and bronchoalveolar lavage specimens were spiked with potential interfering substances at the concentrations above. Each substance was tested with triplicate extractions. No false positive results were observed for any of the substances at the concentrations tested.

Table 10 Interfering substances

Interfering substance	Final concentration in sample	Agreement with expected results			
		Positive BAL samples	Positive NP samples	Negative BAL samples	Negative NP samples
None	N/A	100% <sup>[1]</sup>	100%	100%	100%
Mucin: bovine submaxillary gland, type I-S	0.1 mg/mL	100% <sup>[2]</sup>	100%	100%	100%
Blood (human)	1% v/v	100% <sup>[3]</sup>	100%	100%	100%
Nasal sprays or drops—Nasacort™	10% v/v	100% <sup>[4]</sup>	100% <sup>[4]</sup>	100%	100%
Nasal corticosteroids—Dymista™	5 µg/mL	100% <sup>[2]</sup>	100%	100%	100%
NeilMed™ Nasogel™	1% w/v	100% <sup>[2]</sup>	100%	100%	100%
Influenza A H1N1 Brisbane/59/07	1 × 10 <sup>5</sup> TCID <sub>50</sub> /mL	100% <sup>[2]</sup>	100%	100%	100%
Throat lozenges, oral anesthetic and analgesic—Chloraseptic™	1% w/v	100% <sup>[3]</sup>	100%	100%	100%
Oseltamivir phosphate	33 µg/mL	100% <sup>[2]</sup>	100%	100%	100%
Antibiotic, nasal ointment—Bactroban™	5 µg/mL	100% <sup>[2]</sup>	100%	100%	100%
Antibacterial, systemic—Tobramycin	0.6 mg/mL	100% <sup>[2]</sup>	100%	100%	100%
Homeopathic allergy relief medicine—Similasan™ Nasal	10% v/v	100%	100%	100%	100%

<sup>[1]</sup> Two of six replicates produced a C<sub>t</sub> >37 or Undetermined for S Gene, but all replicates were called Positive based on the interpretation algorithm.

<sup>[2]</sup> Two of three replicates produced a C<sub>t</sub> >37 or Undetermined for S Gene, but all replicates were called Positive based on the interpretation algorithm.

<sup>[3]</sup> All three replicates produced a C<sub>t</sub> >37 or Undetermined for S Gene but were called Positive based on the interpretation algorithm.

<sup>[4]</sup> One of three replicates produced a C<sub>t</sub> >37 or Undetermined for S Gene, but all replicates were called Positive based on the interpretation algorithm.

## Cross-reactivity

*In silico* analysis of the following forty-three (43) organisms was performed.

**Table 11** Organisms used for *in silico* cross-reactivity analysis

Human coronavirus 229E	Rhinovirus/Enterovirus
Human coronavirus OC43	Parechovirus
Human coronavirus HKU1	<i>Candida albicans</i>
Human coronavirus NL63	<i>Corynebacterium diphtheriae</i>
SARS-coronavirus	<i>Legionella</i> (non-pneumophila)
MERS-coronavirus	<i>Bacillus anthracis</i> (Anthrax)
Adenovirus	<i>Moraxella catarrhalis</i>
Human Metapneumovirus (hMPV)	<i>Neisseria elongata</i> and <i>Neisseria meningitidis</i>
Parainfluenza 1	<i>Pseudomonas aeruginosa</i>
Parainfluenza 2	<i>Staphylococcus epidermidis</i>
Parainfluenza 3	<i>Streptococcus salivarius</i>
Parainfluenza 4	<i>Leptospira</i> sp.
Influenza A	<i>Chlamydomydia pneumoniae</i>
Influenza B	<i>Chlamydomydia psittaci</i>
Influenza C	<i>Coxiella burnetii</i> (Q-Fever)
Enterovirus	<i>Staphylococcus aureus</i>
Respiratory Syncytial Virus A	<i>Haemophilus influenzae</i>
Respiratory Syncytial Virus B	<i>Legionella pneumophila</i>
<i>Bordetella pertussis</i>	<i>Mycobacterium tuberculosis</i>
<i>Mycoplasma pneumoniae</i>	<i>Streptococcus pneumoniae</i>
<i>Pneumocystis jirovecii</i> (PJP)	<i>Streptococcus pyogenes</i>

Among the tested organisms, *Neisseria elongata* showed homology for the forward and reverse primers and probe for the N gene. The forward primer showed  $\geq 80\%$  homology while the reverse primer and probe showed 36% homology. The N gene reverse primer and probe show low homology, therefore the risk of the non-specific amplification is low.

Blast analysis showed  $\geq 80\%$  homology for one assay component (forward primer, reverse primer, or probe) for select isolates. Despite  $\geq 80\%$  homology of one assay component for select isolates, there is no anticipated amplification because hybridization of all three assay components are necessary to generate a signal. We also found multiple instances where different assay components had  $\geq 80\%$

homology to different isolates of the same species. For example, *Bacillus anthracis* strain AFS029987 had  $\geq 80\%$  homology to the ORF1ab forward primer while strain MCCC 1A01412 had  $\geq 80\%$  homology to the ORF1ab reverse primer. Since these are two different organisms, amplification is not likely to occur. The *in silico* analysis indicates that significant amplification of non-target sequences that result in cross-reactivity or potentially interfere with detection of SARS-CoV-2 is not likely to occur.

## Clinical evaluation

A clinical evaluation study was performed to evaluate the performance of the TaqPath™ COVID-19 CE-IVD RT-PCR Kit using nasopharyngeal swab (NP) and bronchoalveolar lavage (BAL) specimens.

A total of sixty (60) contrived positive specimens were tested:

- 30 contrived positive nasopharyngeal swab (NP) specimens
- 30 contrived positive bronchoalveolar lavage (BAL) specimens

Samples were contrived by spiking known concentrations of extracted SARS-CoV-2 viral genomic RNA, relative to the product LoD, into matrices which were determined to be negative by the TaqPath™ COVID-19 CE-IVD RT-PCR Kit prior to spiking in the RNA.

In addition to the contrived positive specimens, sixty (60) negative specimens were tested:

- 30 negative nasopharyngeal swab (NP) specimens
- 30 negative samples bronchoalveolar lavage (BAL) specimens
- All negative samples yielded negative results

Results for positive samples are shown in the tables below:

**Table 12 BAL Clinical Evaluation Study**

Final RNA Concentration in Sample	Number of Positives	Mean C <sub>t</sub>		
		S gene	ORF1ab	N gene
2X LoD	20/20 <sup>[1]</sup>	28.9	29.5	28.7
3X LoD	5/5 <sup>[1]</sup>	28.8	29.2	28.5
5X LoD	5/5	27.4	28.2	27.4

<sup>[1]</sup> Two samples initially gave inconclusive results and were retested. The results were positive after the retest. Mean C<sub>t</sub> values are calculated from the retest results.

Table 13 NP Clinical Evaluation Study

Final RNA Concentration in Sample	Number of Positives	Mean C <sub>t</sub>		
		S gene	ORF1ab	N gene
2X LoD	20/20 <sup>[1]</sup>	30.9	30.6	29.3
3X LoD	5/5	30.0	30.1	28.8
5X LoD	5/5	28.7	29.0	27.9

<sup>[1]</sup> One sample initially gave an inconclusive result and was retested. The result was positive after the retest. Mean C<sub>t</sub> values are calculated from the retest results.

## Clinical performance evaluation

An independent clinical study was performed by MVZ Labor Dr. Limbach & Kollegen GbR (Heidelberg, Germany), an ISO 15189 and ISO 17025 accredited clinical laboratory for the performance of molecular methods (genetic testing) and genome sequencing (DAkkS D-ML-13176-01-00), as well as authorized laboratory for COVID-19 diagnostic activities in Germany.

The study protocol was reviewed and approved by the Ethical Committee of the Baden-Württemberg State Medical Association, and executed in line with the Good Clinical Practices.

The study was conducted on retrospectively collected upper respiratory samples from routine leftover clinical specimen from patients of all ages, collected in Germany in February 2021. All samples had been characterized for SARS-CoV-2 presence/absence by the CE-IVD RT-PCR test from Roche (Roche cobas® SARS-CoV-2 Test).

A total of 450 clinical samples were tested in the study. 175 samples were previously characterized positive, and spanning the dynamic range of the assay C<sub>t</sub> values in equal proportions (range bins: C<sub>t</sub> <20, 20 ≤ C<sub>t</sub> <25, and C<sub>t</sub> ≥30), and 275 clinical samples were previously characterized negative. They were concurrently tested with the comparator and the TaqPath™ COVID-19 CE-IVD RT-PCR Kit, in a blinded and randomized study.

The Roche cobas® SARS-CoV-2 Test was used as the comparator.

NGS sequencing (Illumina® NGS sequencing instrument with EasySeq™ SARS-CoV-2 Whole Genome NGS Sequencing Kit) was used as the resolver for any discordant results between both CE-IVD RT-PCR tests.

### Concordance, Positive Percent Agreement (PPA), and Negative Percent Agreement (NPA):

Positive Percent Agreement (PPA) and Negative Percent Agreement (NPA) and the two-sided 95% Confidence Interval (CI) were calculated as follows:

$$\text{PPA} = 100\% \times (\text{Number of positives by both evaluation method and its comparator} / \text{Number of positives by comparator method})$$

$$\text{NPA} = 100\% \times (\text{Number of negatives by both evaluation method and its comparator} / \text{Number of negatives by comparator method})$$

The 95% *Confidence Interval* was determined using Clopper Pearson Method.

**Table 14** Concordance between the TaqPath™ COVID-19 CE-IVD RT-PCR Kit and the comparator

COVID-19	Comparator: Roche cobas® SARS-CoV-2 Test		
TaqPath™ COVID-19 Kit (CE-IVD)	Positive	Negative	Total
Positive	171	13	184
Negative	1	261	262
Total	172	274	446
Agreement Calculations			95% CI
Positive Percent Agreement (PPA)	171/172	99.4%	96.8% – 100.0%
Negative Percent Agreement (NPA)	261/274	95.3%	92.0% – 97.4%

Four samples with only one target gene as positive by cobas® SARS-CoV-2 Assay were deemed as “Inconclusive” per the study protocol and excluded from the agreement calculation.

The TaqPath™ COVID-19 test and the comparator had fourteen (14) discordant samples, amongst which eight (8) could not be resolved by sequencing due to low concentration of SARS-CoV-2-RNA in the specimens. The remaining six (6) discordant samples were resolved by sequencing, showing presence of SARS-CoV-2, in concordance with the TaqPath™ COVID-19 results.

# C<sub>t</sub> cutoff values for assay targets

The Applied Biosystems™ COVID-19 Interpretive Software uses the following C<sub>t</sub> cutoff values for assay targets during interpretation of the results.

**Table 15 COVID-19 Assay C<sub>t</sub> cutoff values**

Sample or Control	Target	C <sub>t</sub> cutoff
Positive Control	MS2	Valid C <sub>t</sub> values are >37
	Viral targets	Valid C <sub>t</sub> values are ≤37
Negative Control	MS2	Valid C <sub>t</sub> values are ≤32
	Viral targets	Valid C <sub>t</sub> values are >37
Clinical samples	MS2	Valid C <sub>t</sub> values are ≤32 <sup>[1]</sup>
	Viral targets	Positive C <sub>t</sub> values are ≤37

<sup>[1]</sup> If any of the viral targets is positive, the C<sub>t</sub> for MS2 can be >32.

# Documentation and support

## Related documentation

Document	Publication Number
<i>Applied Biosystems™ 7500 Fast Dx Real-Time PCR Instrument Reference Guide</i>	4406991
<i>Applied Biosystems™ 7300/7500/7500 Fast Real-Time PCR System Installation and Maintenance Guide</i>	4347828
<i>QuantStudio™ 3 and 5 Real-Time PCR Systems Installation, Use, and Maintenance Guide</i>	MAN0010407
<i>QuantStudio™ 5 Dx Real-Time PCR Instrument Maintenance and Administration User Guide</i>	100042186
<i>QuantStudio™ 5 Dx TD Software User Guide</i>	100049554
<i>QuantStudio™ 5 Dx IVD Software User Guide</i>	100049556
<i>QuantStudio™ 6 and 7 Flex Real-Time PCR Systems Maintenance and Administration Guide</i>	4489821
<i>QuantStudio™ 6 and 7 Flex Real-Time PCR Systems Quick Reference</i>	4489826
<i>QuantStudio™ Real-Time PCR Software Getting Started Guide</i>	4489822
<i>COVID-19 Interpretive Software CE-IVD Edition Installation Quick Reference</i>	MAN0019230
<i>TaqPath™ COVID-19 CE-IVD RT-PCR Kit with saliva Instructions for Use</i>	MAN0024936

## Customer and technical support

For additional documentation and information about this kit, visit:

<https://www.thermofisher.com/covid19ceivd>

For download instructions for the COVID-19 Interpretive Software, see “Obtain the Applied Biosystems™ COVID-19 Interpretive Software” on page 55.

Refer to the Read Me file provided with the COVID-19 Interpretive Software before contacting support for the software.

Visit: <https://www.thermofisher.com/contactus> for service and support information for this kit, including the following:

- Worldwide contact telephone numbers
- Product support information
- Order and web support
- Product documentation such as:
  - Certificates of Analysis
  - Safety Data Sheets (SDSs; also known as MSDSs)

---

**Note:** For SDSs for reagents and chemicals from other manufacturers, contact the manufacturer.

---

## Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at [www.thermofisher.com/us/en/home/global/terms-and-conditions.html](http://www.thermofisher.com/us/en/home/global/terms-and-conditions.html). If you have any questions, please contact Life Technologies at [www.thermofisher.com/support](http://www.thermofisher.com/support).

# QuantStudio™ 3 and 5 Real-Time PCR Systems

## INSTALLATION, USE, AND MAINTENANCE

### Firmware v1.2.x

Publication Number MAN0010407

Revision C.0

For Research Use Only. Not for use in diagnostic procedures.

**ThermoFisher**  
SCIENTIFIC

The information in this guide is subject to change without notice.

## **DISCLAIMER**

TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

### **Limited Use Label License No. 474: Real-Time PCR System for Research Use Only**

Notice to Purchaser: The purchase of this instrument conveys to the purchaser the limited, non-transferable right to use the purchased instrument only, under intellectual property rights that are owned and/or controlled by Life Technologies and relate specifically to the instrument. Purchase of the instrument includes the right to use the instrument for internal research and to perform services (including the right to report the results of services for a fee) by the purchaser only, but does not convey rights to use any other products, reagents, assays or methods such as the 5' nuclease assay process. The sale of this instrument is expressly conditioned on the purchaser not reselling, repackaging, or distributing this instrument, or any of its components, and no such rights are conveyed expressly, by implication, or by estoppel. For information on obtaining additional rights, please contact [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com) or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

### **Limited Use Label License No. 572: Cy<sup>®</sup> and CyDye<sup>®</sup>**

Notice to Purchaser: Cy<sup>®</sup> and CyDye<sup>®</sup> are registered trademarks of GE Healthcare.

## **Corporate entity**

Life Technologies | Carlsbad, CA 92008 USA | Toll Free in USA 1.800.955.6288

## **TRADEMARKS**

Cy is a trademark of GE Healthcare UK Limited. Adobe, Acrobat, and Reader are trademarks of Adobe Systems, Inc. Microsoft, Internet Explorer, and Windows are trademarks of Microsoft Corporation. Apple, Safari, and Macintosh are trademarks of Apple, Inc. Mozilla and Firefox are trademarks of Mozilla Foundation. TaqMan is a registered trademark of Roche Molecular Systems, Inc., used under permission and license. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

©2015 Thermo Fisher Scientific Inc. All rights reserved.

# Contents

About this guide .....	8
Revision history .....	8
Purpose .....	9
■ <b>CHAPTER 1</b> Product information .....	10
Instrument hardware description .....	10
Instrument overview .....	10
Parts of the instrument .....	13
Software description .....	14
Instrument, desktop, and cloud software features .....	14
Folders, templates, experiments, and projects .....	15
Third-party software installation .....	16
Network connection options .....	16
Experiment types .....	17
■ <b>CHAPTER 2</b> Start, sign on, and configure the instrument .....	19
Installation and instrument verification .....	19
Precautions for use .....	20
Power on the instrument .....	20
Parts of the home screen .....	21
Use the instrument without signing in .....	23
Create a new instrument profile .....	23
Sign In .....	23
Sign out .....	24
Link your instrument profile to your Thermo Fisher Cloud account .....	24
Link your instrument profile to Thermo Fisher Cloud .....	24
Create an instrument profile and link to Thermo Fisher Cloud .....	24
Configure instrument settings .....	25

■ <b>CHAPTER 3</b>	<b>Create and run experiments on the instrument .....</b>	<b>28</b>
Workflow .....		28
Options for running an experiment .....		29
Create and run an experiment from a template .....		29
Run an experiment from a saved file .....		29
Repeat your last experiment run .....		30
Edit an experiment before starting a run .....		31
Define experiment properties .....		31
Edit a method .....		32
Define plate wells .....		35
Load and unload a plate in the instrument .....		36
View, pause, or stop a run .....		37
View well details .....		37
Pause or stop a run .....		37
Adjust the graphical view of an experiment .....		37
Lock the touchscreen during a run .....		38
Transfer, view, or manage files and results .....		38
Transfer experiment results .....		38
View run history .....		38
Manage template (.edt) files .....		39
■ <b>CHAPTER 4</b>	<b>Calibrate and verify instrument performance .....</b>	<b>40</b>
Calibration and verification schedule .....		40
Calibration descriptions .....		41
View the calibration status and set reminders .....		41
View calibration status and set reminders in the instrument .....		41
View calibration status and set reminders in the Cloud .....		42
Perform ROI/Uniformity, Background, and Dye calibrations .....		42
Workflow: Calibration .....		42
Prepare a calibration plate .....		43
Perform calibrations .....		44
View calibration images and transfer results to USB .....		46
Troubleshoot calibration failure .....		47
Identify contamination .....		47
Create a background plate <i>(optional)</i> .....		48
Perform instrument verification using RNase P plates .....		49
Instrument verification description and schedule .....		49
RNase P Instrument Verification Plates .....		49
Analytical performance .....		50
Installation specification .....		50
Prepare an RNase P plate .....		50
Perform RNase P verification .....		51
Troubleshoot verification failure .....		53

Calibrate custom dyes .....	54
Custom dyes overview .....	54
Dilute the custom dye to an optimal concentration .....	54
Calibrate the custom dye .....	57
Calibrate for a custom melt curve experiment .....	60
<b>■ CHAPTER 5 Maintain the instrument .....</b>	<b>62</b>
Backup or restore the instrument .....	62
Decontaminate the sample block .....	63
Materials required .....	63
Clean the sample block .....	63
Solvents for cleaning the sample block .....	65
Replace the instrument fuses .....	65
Materials required .....	65
Replace the fuses .....	66
Power on or off, store, and move .....	66
Enable sleep mode .....	66
Power on the instrument .....	67
Power off the instrument .....	67
Prepare the instrument to ship, move, or store .....	67
Move the instrument .....	68
Return the instrument for service .....	68
<b>■ CHAPTER 6 Profile and instrument configuration tasks .....</b>	<b>70</b>
Initial start up .....	70
Manage profiles .....	71
Create an administrator profile .....	71
View all user profiles (Administrator only) .....	71
Manage all instrument profiles (Administrator only) .....	71
Edit a user profile .....	72
Enable SAE mode (Administrator only) .....	72
Require sign-in (Administrator only) .....	73
Enable remote instrument monitoring (Administrator only) .....	73
Update instrument software (Administrator only) .....	73
Restore factory defaults .....	74
Configure the network .....	74
Set up a wired connection .....	74
Set up a wireless connection .....	75
Select a Cloud region (Administrator only) .....	75
Manage the instrument name (Administrator only) .....	75
Set the date and time .....	76
Manage the Sign Out Timer (Administrator only) .....	76

■	<b>APPENDIX A</b>	<b>Install and connect the instrument to a network .....</b>	<b>77</b>
		Workflow: Install and connect to a network .....	77
		Before you begin .....	78
		Unpack and install the instrument .....	78
		Power on and follow the startup wizard .....	79
		Connect the instrument and the computer directly or to a LAN .....	79
		Instrument and computer connections .....	81
		Download and install the QuantStudio™ Design and Analysis desktop Software .....	82
		Computer requirements for the desktop software .....	82
		Download the desktop software .....	82
		Install the software .....	82
		Networking .....	84
		Supported network configuration options .....	84
		Control and monitor networked instruments .....	85
		About the Ethernet port .....	85
		Firewall ports that must be open .....	85
		Networking guidelines and best practices .....	85
■	<b>APPENDIX B</b>	<b>Troubleshooting .....</b>	<b>87</b>
■	<b>APPENDIX C</b>	<b>Parts and materials .....</b>	<b>88</b>
		Kits, consumables, and accessories .....	88
		96-Well (0.2-mL) consumables .....	88
		384-well consumables .....	89
		96-Well Fast (0.1-mL) consumables .....	90
		QuantStudio™ 3 and 5 Systems accessories .....	91
		General-use materials and consumables .....	91
■	<b>APPENDIX D</b>	<b>Instrument specification and layout .....</b>	<b>92</b>
		Configured system dimensions .....	93
		Instrument clearances .....	94
		Electrical requirements .....	94
		Environmental requirements .....	95
		Network requirements .....	95

■ <b>APPENDIX E Safety</b> .....	<b>97</b>
Symbols on this instrument .....	98
Conformity symbols .....	99
Safety alerts on this instrument .....	100
Location of safety labels on the instrument .....	100
Safety information for instruments not manufactured by Thermo Fisher Scientific ....	101
Instrument safety .....	101
General .....	101
Physical injury .....	102
Electrical .....	103
Cleaning and decontamination .....	103
Safety and electromagnetic compatibility (EMC) standards .....	103
Safety compliance .....	104
EMC .....	104
Environmental design .....	105
Chemical safety .....	105
Biological hazard safety .....	106
 <b>Documentation and support</b> .....	 <b>107</b>
Related documentation .....	107
Obtain information from the Help system .....	108
Customer and technical support .....	108
Limited product warranty .....	109
 <b>Index</b> .....	 <b>110</b>

# About this guide

## Revision history

Revision	Date	Description
C.0	December 2015	Updates – describes new features in Firmware v1.2.x, including: <ul style="list-style-type: none"><li>• Updated workflow to create an instrument profile</li><li>• Streamlined workflow for connectivity to the Thermo Fisher Cloud</li><li>• Improved display of VeriFlex™ Zones</li><li>• Plate insert reminder before starting run</li><li>• Support for 96-well Fast (0.1 mL) plates</li></ul>
B.0	September 2015	Updates – describes new features in Firmware v1.1.x, including: <ul style="list-style-type: none"><li>• <b>Experiment runs:</b> monitor real-time data, edit cycle number and lock screen during a run, view end plot, and support for 384-well plates</li><li>• <b>File management:</b> access Network folders, navigate folder structures and save templates on instrument, and perform batch actions for file management</li><li>• <b>Instrument configuration:</b> support for Security, Audit, and e-Signature (SAE), smart monitoring, and choice of server region</li></ul> Includes information on software feature comparison, definitions of terms, parts of method, network connection options, experiment types, desktop software installation.
A.0	April 2015	New document. Describes installation, operation, and maintenance of the QuantStudio™ 3 and 5 Real-Time PCR Systems.

## Purpose

This guide provides information about installing, using, and maintaining the QuantStudio™ 3 and 5 Real-Time PCR Systems.

**Note:** For information and instructions on performing experiments on these systems, refer to the *QuantStudio™ Design and Analysis desktop Software User Guide* (Pub. no. MAN0010408).

# Product information

- Instrument hardware description ..... 10
- Software description ..... 14
- Network connection options ..... 16
- Experiment types ..... 17

## Instrument hardware description

### Instrument overview

- The QuantStudio™ 3 and 5 Real-Time PCR Systems use fluorescent-based polymerase chain reaction (PCR) reagents to perform:
- Quantitative detection of target nucleic acid sequences (targets).
  - Qualitative detection of targets (endpoint analysis).
  - Qualitative analysis of the PCR product (post-PCR melt curve analysis).

The following fixed-block configurations are available:

QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 5 Real-Time PCR System
<ul style="list-style-type: none"> <li>• 96-Well VeriFlex™ 0.2-mL Block (4 Color)</li> <li>• 96-Well VeriFlex™ 0.1-mL Block (4 Color)</li> </ul>	<ul style="list-style-type: none"> <li>• 96-Well VeriFlex™ 0.2-mL Block (6 Color De-coupled)</li> <li>• 384-Well Block (5 Color)</li> <li>• 96-Well VeriFlex™ 0.1-mL Block (6 Color De-coupled)</li> </ul>

The instrument can be run directly from the touchscreen to create and start experiments. To design experiments or to analyze data, the instrument can be integrated with QuantStudio™ Design and Analysis Software (desktop or cloud).

An optional barcode scanner and optional wireless adapter can be purchased separately.

## Instrument filters and supported dyes

### System dyes


The QuantStudio™ 3 and 5 Systems use a coupled four-color, coupled five-color, or decoupled six-color filter set that supports the dyes shown in the following table and figure. For more information about the spectral dye calibration kits available for the QuantStudio™ 3 and 5 Systems, contact Support.

Peak channel	Color	Filter wavelength (nm) <sup>[1]</sup>		Pre-calibrated dyes	Example custom dyes
		Excitation	Emission		
x1-m1	Blue	470 ± 15	520 ± 15	FAM™ and SYBR™ Green	SYT09
x2-m2	Green	520 ± 10	558 ± 12	VIC™	HEX™, TET™, and JOE™ <sup>[2]</sup>
x3-m3	Yellow	550 ± 10	587 ± 10	ABY™, NED™, and TAMRA™	Cy <sup>fi</sup> 3
x4-m4	Orange	580 ± 10	623 ± 14	JUN™ and ROX™	Texas Red™
x5-m5	Red	640 ± 10	682 ± 14	Cy <sup>fi</sup> 5 and MUSTANG PURPLE™	LIZ™
x6-m6	Deep-Red	662 ± 10	711 ± 12	None <sup>[3]</sup>	Cy <sup>fi</sup> 5.5

<sup>[1]</sup> The central wavelengths are the optimized wavelengths.

<sup>[2]</sup> The HEX™ and TET™ dyes from Thermo Fisher Scientific fall within the emission wavelength range of the system, therefore they can be added and adapted for use in experiments on the system. To add any of these dyes to the Dye Library, perform a custom dye calibration.

<sup>[3]</sup> This filter set currently does not support any dyes supplied by Thermo Fisher Scientific.

① x1-m1 — FAM™, SYBR™ Green

② x2-m2 — VIC™

③ x3-m3 — ABY™, NED™, Cy<sup>fi</sup> 3, and TAMRA™

④ x4-m4 — JUN™, ROX™, and Texas Red™

⑤ x5-m5 — Cy<sup>fi</sup> 5 and MUSTANG PURPLE™

## Custom dyes

The QuantStudio™ 3 and 5 Systems can run assays designed with custom dyes (dyes not supplied by Thermo Fisher Scientific or dyes not pre-calibrated with the instrument). Custom dyes must excite between 455 – 672 nm and emit between 505 – 723 nm. Select a custom dye that does not overlap with other dyes used in an experiment (see the filter-wavelength table in “System dyes” on page 11).

To use a custom dye and add it to the Dye Library, perform a custom dye calibration (see “Calibrate custom dyes” on page 54).

## About data collection



The instrument collects raw fluorescence data at different points during the PCR cycle, depending on the type of run performed.

When you create an experiment template (.edt file) in the desktop or cloud software, you can customize the optical filter channels through which the instrument collects data. You can specify a filter channel set for all PCR thermal protocols and, optionally, a different filter set for the melt curve stage(s).

Run type	Experiment type	Data collection point
Real-time	<ul style="list-style-type: none"> <li>Standard curve</li> <li>Relative standard curve</li> <li>Comparative C<sub>T</sub> (<math>\Delta\Delta C_T</math>)</li> <li>Melt curve</li> </ul>	<p>During the thermal cycling protocol.</p> <p>Typical timing is to collect data at each cycle of a PCR stage or continuously during a melt stage.</p>
Post-PCR (endpoint)	<ul style="list-style-type: none"> <li>Genotyping</li> <li>Presence/ Absence</li> </ul>	<ul style="list-style-type: none"> <li>After thermal cycling is completed. For Presence/Absence and Genotyping experiments, data collection before the PCR cycle is optional, but recommended.</li> <li><i>[Optional]</i> Before thermal cycling starts. Collecting data during the run can confirm genotyping results by viewing traces in allelic discrimination plots or viewing genotyping calls at earlier cycles.</li> </ul>

## Blocks with VeriFlex™ Zones


Applied Biosystems™ VeriFlex™ technology provides independent temperature zones that offer enhanced functionality and precise control over your qPCR runs.

QuantStudio™ 3 System	QuantStudio™ 5 System
Three programmable VeriFlex™ Zones.	Six programmable VeriFlex™ Zones. <sup>[1]</sup>
	

<sup>[1]</sup> Only applicable for 96-well 0.2-mL and 0.1-mL blocks.

These independent zones are ideal for qPCR optimization or the ability to run multiple experiments in the same run. Unlike standard gradients which give a sigmoidal temperature curve across the columns, blocks with VeriFlex™ Zones help deliver accurate temperatures across every zone.

## Parts of the instrument

- ① **Touchscreen** – Controls the instrument.
- ② **USB port** – For connection to an external network drive, jump drive, or other external data storage device.
- ③ **Instrument drawer** – Contains sample plate.

The instrument includes three additional USB ports on the back of the instrument. The instrument recognizes only one external storage device at a time for data transfer.

## Software description


### Instrument, desktop, and cloud software features

The instrument software and the QuantStudio™ Design and Analysis Software (desktop and cloud) include the features described below.

Feature	Instrument	Desktop	Cloud
Use as guest (no sign in)	✓	✓ <sup>[1]</sup>	—
Create templates (unlocked or locked)	—	✓	✓
Edit unlocked templates	✓	✓	✓
Edit locked templates (password assigned by template creator required)	—	✓	✓
Load system or user-created templates (.edt file) to instrument	✓	✓	—
Change experiment settings in template (.edt file) loaded on the instrument  Settings that can be changed in a locked template (no password required): <ul style="list-style-type: none"> <li>• Properties: All settings</li> <li>• Method: No changes allowed</li> <li>• Plate: Sample names</li> </ul> Settings that can be changed in an unlocked template: <ul style="list-style-type: none"> <li>• Properties: All settings</li> <li>• Method: All settings</li> <li>• Plate: Sample names</li> </ul>	✓	—	—
Load plate in instrument	✓	—	—
Start run	✓	✓	—
View real-time data during a run	✓	✓	✓
View instrument status (running, calibration needed, and so on)	✓	—	✓
Analyze results	—	✓	✓
Set calibration reminders	✓	—	✓
Review exported calibration or RNase P verification results	—	✓	✓ (only RNase P)

<sup>[1]</sup> Login required if the Security, Audit, and e-Signature (SAE) module is enabled (QuantStudio™ 5 Systems only)

## Folders, templates, experiments, and projects

Term	Definition	Supported in
Folders in Load Experiment 	Location in which you can store templates (.edt files) on the instrument: <ul style="list-style-type: none"> <li>• <b>My Instrument</b> – Displayed if you are signed in</li> <li>• <b>Public</b> – Location in which all experiments run by guest users are stored</li> <li>• <b>USB</b> – USB for manual transfer to and from a computer</li> <li>• <b>Post Read</b> – If you run an endpoint experiment, the pre-read experiment is automatically saved in this folder for post-read analysis</li> </ul>	Instrument
Experiment Template File (.edt) 	Default settings for an instrument run, can be modified before instrument run.  Two types: <ul style="list-style-type: none"> <li>• Factory-provided, accessed from Open Template (instrument) or Create Experiment (desktop or cloud)</li> <li>• User-created, accessed from Load Experiment</li> </ul>	Instrument Desktop Cloud
Experiment Run File (.eds)	Settings and data for a completed instrument run.	Instrument Desktop Cloud
Project 	Function in the Cloud Data Manager that is used for secondary analysis applications. Does not apply to the QuantStudio™ Design and Analysis cloud Software, which is a primary analysis application.	Cloud

Template and experiment components	Instrument	Desktop	Cloud
Properties	<ul style="list-style-type: none"> <li>• Experiment file name</li> <li>• Reagent information (Reagent barcode, Lot#)</li> <li>• Plate barcode</li> <li>• Data destination (location for auto-transfer of data)</li> <li>• Comments (Tags are not used at this time)</li> </ul>	<ul style="list-style-type: none"> <li>• Experiment file name</li> <li>• Plate barcode</li> <li>• User name</li> <li>• Instrument type</li> <li>• Block type</li> <li>• Experiment type</li> <li>• Chemistry (reagent information)</li> <li>• Run mode</li> <li>• Comments</li> <li>• {Cloud software only} Notifications</li> </ul>	
Method	Thermal cycling conditions	Thermal cycling conditions	

Template and experiment components	Instrument	Desktop	Cloud
Plate	Sample names You cannot edit targets/SNP assays or tasks on the instrument.	Define and assign samples, targets or SNP assays, and tasks in the Quick Setup and Advanced Setup panes of the Plate tab.	
Run	Start and monitor a run in progress View: Time remaining and Temperature, Method, Plots Pause, Resume, Stop a run	Start and monitor a run in progress View: Time remaining and Temperature, Method, Plots	Monitor a run in progress (link to Thermo Fisher Cloud Instrument Details) View: Time remaining and Temperature, Method, Plots
Results	Not applicable	Review plots	
Export	Not applicable	Export results	

### Third-party software installation

Before you install third-party software on the computer running the QuantStudio™ desktop Software, confirm that the third-party software will not:

- Restrict Ethernet communication.
- Interfere with instrument or computer operation.

## Network connection options

You can connect a QuantStudio™ 3 or 5 Instrument to a network or computer in the configurations listed below. For specific information on networking, see “Networking” on page 84.

Direct connection	Local area network (LAN) connection	Cloud connection
wired	wired <i>or</i> wireless	wired <i>or</i> wireless
		

## Experiment types

Purpose	Description
<b>Standard Curve experiment</b>	
Determines absolute target quantity in samples.	<ol style="list-style-type: none"> <li>1. The software measures amplification of the target in a standard dilution series and in test samples.</li> <li>2. The software generates a standard curve using data from the standard dilution series.</li> <li>3. Using the standard curve, the software interpolates the absolute quantity of target in the test samples.</li> </ol>
<b>Relative Standard Curve experiment</b>	
Determines relative target quantity in samples.	<ol style="list-style-type: none"> <li>1. The software measures amplification of the target of interest and of an endogenous control target in a standard dilution series, in a reference (calibrator) sample, and in test samples.  The endogenous control is a target that is expressed equally in all samples; examples of endogenous controls are <math>\beta</math>-actin, GAPDH, and 18S ribosomal RNA. The software can algorithmically incorporate multiple endogenous control targets in relative quantification calculations.  The reference sample is used as the basis for relative quantification results (or 1× sample). For example, in a study of drug effects on gene expression, an untreated control is an appropriate reference sample.</li> <li>2. The software generates standard curves for the target of interest and the endogenous control using data from the corresponding standard dilution series.</li> <li>3. Using the standard curves, the software interpolates the quantities of the target of interest and the endogenous control in each sample. The target quantity in each sample is then normalized to the endogenous control quantity in the sample.</li> <li>4. To determine the relative quantity of the target in test samples, the software divides the normalized target quantity in the sample by the normalized target quantity in the reference sample.</li> </ol>
<b>Comparative <math>C_T</math> (<math>\Delta\Delta C_T</math>) experiment</b>	
Determines relative target quantity in samples.	<ol style="list-style-type: none"> <li>1. The software measures amplification of the target of interest and of an endogenous control target in a reference (calibrator) sample and in test samples.  The endogenous control is a target that is expressed equally in all samples; examples of endogenous controls are <math>\beta</math>-actin, GAPDH, and 18S ribosomal RNA. The software can algorithmically incorporate multiple endogenous control targets in relative quantification calculations.  The reference sample is used as the basis for relative quantification results (or 1× sample). For example, in a study of drug effects on gene expression, an untreated control is an appropriate reference sample.</li> <li>2. The measurements for the target of interest are normalized to the endogenous control.</li> <li>3. To determine the relative quantity of the target in test samples, the software compares the normalized <math>C_T</math> (<math>\Delta C_T</math>) for the sample to the normalized <math>C_T</math> (<math>\Delta C_T</math>) for the reference sample.</li> </ol>

Purpose	Description
<b>Genotyping experiment</b>	
Detects single nucleotide polymorphism (SNP) variants of a target nucleic acid sequence.	<p>Genotyping experiments use preformulated TaqMan™ SNP Genotyping Assays that include:</p> <ul style="list-style-type: none"> <li>Two sequence-specific primers for amplification of sequences containing the SNP of interest.</li> <li>Two allele-specific TaqMan™ probes for Allele 1 and Allele 2.</li> </ul> <ol style="list-style-type: none"> <li>The software normalizes the fluorescence of the reporter dyes to the fluorescence of the passive reference dye in each well.</li> <li>The software plots the normalized intensities (Rn) of the reporter dyes in each sample well on an Allelic Discrimination Plot, which contrasts the reporter dye intensities of the allele-specific probes.</li> <li>The software algorithmically clusters the sample data, and assigns a genotype call to the samples of each cluster according to its position on the plot.</li> </ol>
<b>Presence/Absence experiment</b>	
Determines the presence or absence of a target nucleic acid sequence in a sample.	The software calls the target present or absent based on an algorithmically determined call threshold. (The call threshold is different from the C <sub>T</sub> threshold; the C <sub>T</sub> threshold is not used to make calls.)
<b>Melt Curve experiment</b>	
Determines the melting temperature (T <sub>m</sub> ) of the amplification products of a PCR reaction that used SYBR™ Green dye.	In the software, Melt Curve analysis is included in the default run method for any experiment type that uses SYBR™ Green reagents.

# Start, sign on, and configure the instrument

■ Installation and instrument verification .....	19
■ Precautions for use .....	20
■ Power on the instrument .....	20
■ Parts of the home screen .....	21
■ Use the instrument without signing in .....	23
■ Create a new instrument profile .....	23
■ Sign In .....	23
■ Sign out .....	24
■ Link your instrument profile to your Thermo Fisher Cloud account .....	24
■ Configure instrument settings .....	25

## Installation and instrument verification

Before the first use of the instrument:

- Install the instrument (see “Unpack and install the instrument” on page 78).
- Verify instrument performance (see “Perform instrument verification using RNase P plates” on page 49).

**Note:** Instruments are factory calibrated, so no calibration is necessary at installation. However, we recommend that you verify instrument performance before using the instrument.

**Note:** Regular calibration and verification should be performed according to the “Calibration and verification schedule” on page 40.

## Precautions for use

**CAUTION! PHYSICAL INJURY HAZARD.** Do not remove the instrument cover. There are no components inside the instrument that you can safely service yourself. If you suspect a problem, contact technical support.

**CAUTION! FIRE HAZARD.** For continued protection against the risk of fire, replace fuses only with listed and certified fuses of the same type and rating as those currently in the instrument.

**CAUTION! PHYSICAL INJURY HAZARD.** During instrument operation, the sample block temperature can reach 100°C. Allow it to cool to room temperature before handling.

**CAUTION!** Before using a cleaning or decontamination method other than those recommended by Thermo Fisher Scientific, confirm with Thermo Fisher Scientific that the proposed method will not damage the equipment.


**CAUTION!** Use flat caps for 0.2 mL tubes and 0.1 mL tubes. Rounded caps can damage the heated cover.

## Power on the instrument

To power on the instrument from a powered-off state:

1. Touch anywhere on the touchscreen to determine if the instrument is in sleep mode. If the home screen displays, the instrument is already powered on.
2. If the home screen does not display, power on the instrument by pressing the switch on the rear panel.

The instrument is ready to use when the home screen is displayed.

If left unattended (for about two hours), the instrument automatically enters sleep mode (enabled by default) to conserve power. Refer to the touchscreen Help system for step-by-step instructions for changing the sleep mode setting.

## Parts of the home screen

- ① Avatar and Instrument name
- ② Eject icon
- ③ Help icon
- ④ Status dial
- ⑤ Current user name; instrument block type
- ⑥ Settings button
- ⑦ Buttons for accessing .edf files
- ⑧ Connectivity icons
- ⑨ Sign In (or My Profile) button

**Table 1** Parts of the home screen

Element of the home screen	Function	For more information, see...
Avatar and Instrument name	Set by the administrator to uniquely identify instrument.	"Manage the instrument name (Administrator only)" on page 75.
Eject icon	Touch to open or close the instrument drawer.	—
Help icon	Touch to launch the touchscreen Help system to access step-by-step instructions.	—

Element of the home screen	Function	For more information, see...
Status dial	<ul style="list-style-type: none"> <li>When the instrument is in use – Displays the sample block temperature, the elapsed run time, and the run status. You can swipe the dial to the left or touch  to access real-time views of the run.</li> <li>When the instrument is not in use – Displays Set up run. You can start a run by touching the status dial.</li> <li>When the instrument is locked – Displays a lock icon () within the status dial.</li> </ul>	—
Current user name and block type	<p>Displays the user name of the current signed-in user and the instrument block type.</p> <p><b>Note:</b> If no user is signed-in, the instrument defaults to the Guest profile.</p>	—
Settings button	Touch  <b>Settings</b> to configure, calibrate, or learn about the instrument.	“Configure instrument settings” on page 25.
Buttons for accessing experiment and template files	<b>Load Experiment</b> – Touch to open a user-created .edt file from a Cloud account, USB, instrument folder, or network drive.	“Run an experiment from a saved file” on page 29.
	<b>Open Template</b> – Touch to open a system .edt template file.	“Create and run an experiment from a template” on page 29.
	<p><b>Run Last</b> – Touch to open the last .edt file run on the instrument.</p> <ul style="list-style-type: none"> <li>If you are signed-in, the file opened will be the last file you ran when signed-in.</li> <li>If you are not signed-in, the file opened will be the last file the guest profile ran.</li> </ul>	“Repeat your last experiment run” on page 30.
Connectivity icons	<ul style="list-style-type: none"> <li> – The instrument is connected to a wired network.</li> <li> – The instrument is wirelessly connected.</li> <li> – A USB drive is plugged into the instrument.</li> <li> – The instrument is linked to a Cloud account.</li> </ul>	—
Sign In button (My Profile button when a user is signed in)	<ul style="list-style-type: none"> <li>Touch  <b>Sign In</b> to sign into an instrument profile and link to a Cloud account.</li> <li>Touch  <b>My Profile</b> to change instrument profile settings, link to a Cloud account, or lock the instrument during a run.</li> </ul>	“Sign In” on page 23.

## Use the instrument without signing in

If the instrument is configured by an Administrator to allow guest access ( **Settings** ▶ **Manage Users** ▶ **Sign In Required** set to off), you can use the instrument without signing in.

If you do not sign in to the instrument:

- All actions are logged to Guest user profile.
- You have access only to the Public folder for selecting and storing experiments.
- You cannot transfer data to the Cloud (only to USB or network drive).

## Create a new instrument profile

In the home screen:

1. Touch  **Sign In**.
2. Touch **Get Started**.
3. Touch **Name**, then enter a username and touch **Done**.
4. Touch **PIN Code**, then enter a four-digit numerical password and touch **Enter**.  
**Note:** Touch the **Show PIN** checkbox to switch PIN display on or off.
5. Touch **Confirm PIN** and repeat the previous step.
6. Touch **Create profile**.  
**Note:** “Link your instrument profile to your Thermo Fisher Cloud account” on page 24 at this point, or you can link after signing in.
7. “Sign In” on page 23 to the profile you just created.

## Sign In

“Create a new instrument profile” on page 23 before signing into the instrument.

In the home screen:

1. Touch  **Sign In**.
2. Touch **Sign In**, then select your username.
3. Enter your PIN, then touch **Enter**.

“Link your instrument profile to your Thermo Fisher Cloud account” on page 24 once signed in.

## Sign out

In the home screen:

1. Touch  **My Profile**.
2. Touch **Sign Out**.

## Link your instrument profile to your Thermo Fisher Cloud account

There are two ways to link an instrument profile to a Thermo Fisher Cloud account:

“Link your instrument profile to Thermo Fisher Cloud” on page 24

“Create an instrument profile and link to Thermo Fisher Cloud” on page 24

Linking your instrument profile to your Cloud account allows you to:

- View instrument status from the Cloud.
- Download templates stored in your Cloud account to the instrument.
- Transfer results from the instrument to your Cloud account.

### Link your instrument profile to Thermo Fisher Cloud

1. “Sign In” on page 23 to your instrument profile.
2. Touch  **My Profile**.
3. Touch **Cloud**.
4. Enter your Thermo Fisher Cloud username and password, then touch **Link Account**.
5. Touch **Done** to exit the confirmation screen.

Your instrument profile will link to your Thermo Fisher Cloud account.


### Create an instrument profile and link to Thermo Fisher Cloud

1. Complete step 1 to step 6 in “Create a new instrument profile”, or continue from “Initial start up” on page 70.
2. In the Get Started – Cloud screen, enter your Thermo Fisher Cloud username and password.
3. Touch **Link Account**.

Your instrument profile will link to your Thermo Fisher Cloud account.

## Configure instrument settings

Touch  **Settings** in the home screen to configure settings as needed.

Touch  for step-by-step instructions on configuring settings.

Options	Description
<b>Instrument Settings</b>	
Instrument Name (Administrator only)	Instrument name. If connected to a network, the instrument name must be unique. <i>(Optional)</i> Add an avatar image for the instrument (.jpg, .png, or .gif) from a USB.
Sleep Mode	Enable the instrument to enter a standby mode after a set length of inactivity.
Heated Cover Temperature	Set the temperature of the heated cover for instrument operation and standby mode. Enable heated cover to automatically turn off during standby mode.
Network Drive	Specify a default network location for the signed-in user.
Insert Plate Reminder	Enable a reminder to insert a plate before starting a run from the instrument.
OEM Connection Only (Administrator only)	Required for API access to the instrument. When enabled, the QuantStudio™ Design and Analysis Software cannot connect to the instrument. API access to the instrument is exclusive to authorized OEM partners.
Cloud Region (Administrator only)	Specify the regional server location to access the Thermo Fisher Cloud Dashboard. <b>Note:</b> Once you set a Cloud region, restore factory defaults to change the region.
Date/Time	Set time zone and date and time formats.
Network Connection	Set wired or wireless network connection for the instrument.
Restore Factory Defaults	Restore the instrument to the factory settings.  <b>IMPORTANT!</b> Back up the instrument before restoring factory defaults (see “Backup or restore the instrument” on page 62).  <b>Note:</b> When you restore an instrument to its factory defaults: <ul style="list-style-type: none"> <li>• User profiles and files stored on the instrument are deleted, including all user-created .edt and .eds files and any custom dye and custom melt calibrations.</li> <li>• System templates and factory calibrations remain on the instrument.</li> </ul>

Options	Description
SAE Mode (Administrator only) (QuantStudio™ 5 Systems only)	When enabled: <ul style="list-style-type: none"> <li>• The home screen displays SAE mode.</li> <li>• Runs must be started from the desktop software. The instrument locks out access to load experiments, open template, run last, start a run, change instrument name and date/time.</li> <li>• Only an administrator can perform calibrations.</li> <li>• A user can monitor a run and view the method from the instrument touchscreen.</li> </ul>
<b>About Instrument</b>	
About Instrument	Displays the instrument Model Name, Block Serial Number, and Firmware Version.
License Agreement	Displays the End User Software License Agreement and the Limited Product Warranty. You can export the License Agreement to a USB drive.
<b>Notifications</b>	
—	Enable home screen notifications of instrument errors and software updates. This function is not related to the Notifications function in the cloud software.  The number of new, unviewed notifications displays over Settings in the home screen.
<b>Maintenance and Service</b>	
Software Update (Administrator only)	Update the instrument software.
Monitoring (Administrator only)	Enable: <ul style="list-style-type: none"> <li>• Remote Monitoring Service to automatically send critical device statistics to Thermo Fisher Scientific (does not monitor or send data).</li> <li>• Thermo Fisher Scientific Cloud Monitor to allow real-time monitoring of instrument runs from a Cloud account.</li> </ul>
Instrument Statistics	Displays instrument usage information: Block Cycle Count and LED Life.
Calibrations	<ul style="list-style-type: none"> <li>• Perform calibrations <ul style="list-style-type: none"> <li>– ROI and Uniformity</li> <li>– Dye</li> <li>– Custom (including Background calibration)</li> </ul> </li> <li>• View calibration history and set calibration reminders in <b>History and Reminders</b></li> </ul>
RNase P Verification	Perform instrument verification using an RNase P plate.
Self Verification Test	Check the instrument hardware functions.
Log	View and export Instrument Log.
Backup / Restore	<ul style="list-style-type: none"> <li>• Backup Instrument</li> <li>• Restore a Backup (Administrator only)</li> </ul>
Ship Prep Mode	Place the instrument in a safe state for shipping, moving, or long-term storage.


Options	Description
<b>Run History</b>	
—	<p>Displays the experiments run on the instrument and whether the data was transferred.</p> <p><b>Note:</b> Experiments run while a user is signed-in can only be viewed or transferred by that user or an administrator.</p> <p>Touch an experiment to view experiment run details and to transfer or delete .eds file.</p>
<b>Manage Users</b>	
Sign In Required (Administrator only)	Enable the restriction of instrument use to only signed-in users (disables Guest profile use).
Sign Out Timer (Administrator only)	Set the duration of inactivity before a user is automatically signed out.
Manage Profiles (Administrator only)	<ul style="list-style-type: none"> <li>• Instrument</li> <li>• Cloud</li> <li>• All Profiles</li> </ul> <p>Displays the profile information for the instrument and the associated Cloud account.</p>

# 3

## Create and run experiments on the instrument

■ Workflow .....	28
■ Options for running an experiment .....	29
■ Edit an experiment before starting a run .....	31
■ Load and unload a plate in the instrument .....	36
■ View, pause, or stop a run .....	37
■ Transfer, view, or manage files and results .....	38

### Workflow

## Options for running an experiment

### Create and run an experiment from a template

In the home screen:

1. Touch  **Open Template**.
  2. (Optional) Touch an experiment category in the left column.
  3. Touch the experiment file name.
  4. (Optional) “Define experiment properties”, including Run File Name (.eds file name), Plate Barcode, Reagent Information, and Data Destination.
  5. (Optional) “Edit a method” on page 32.
    - Add, remove, or edit a step, stage, melt curve, or data collection location.
    - Adjust the cover temperature, sample volume, or number of cycles
    - Configure VeriFlex™ Zones, ramp rate, and pause settings.
  6. (Optional) “Define plate wells” with sample names and view Well ID, Targets, or Dyes.
  7. Load a plate in the instrument (see “Load and unload a plate in the instrument” on page 36).
  8. Touch **Start Run**.  
When prompted, confirm that you inserted a plate.
- Note:** To disable this reminder, select **Do not show again** or select  **Settings ▶ Insert Plate Reminder** in the home screen.

### Run an experiment from a saved file

In the home screen:

1. Touch  **Load Experiment**.
2. Touch  **Cloud**,  **USB**,  **My Instrument**, or  **Network Drive** to navigate to your file location.
  - For files saved to the guest profile, touch  **My Instrument ▶ Public**.
  - For pre and post read files, touch  **My Instrument ▶ Post Read**.
3. Touch the experiment file name.
4. (Optional) “Manage template (.edt) files” in either  **USB** or  **My Instrument**.
5. (Optional) “Define experiment properties”, including Run File Name (.eds file name), Plate Barcode, Reagent Information, and Data Destination.

6. (Optional) “Edit a method” on page 32.
  - Add, remove, or edit a step, stage, melt curve, or data collection location.
  - Adjust the cover temperature, sample volume, or number of cycles
  - Configure VeriFlex™ Zones, ramp rate, and pause settings.
7. (Optional) “Define plate wells” with sample names and view Well ID, Targets, or Dyes.
8. Load a plate in the instrument (see “Load and unload a plate in the instrument” on page 36).
9. Touch **Start Run**.  
 When prompted, confirm that you inserted a plate.  
**Note:** To disable this reminder, select **Do not show again** or select  **Settings ▶ Insert Plate Reminder** in the home screen.

## Repeat your last experiment run

This feature is only applicable to runs started from the instrument and not available for runs started from the desktop software. If you are signed-in, this function applies to the last run from your profile.

In the home screen:

1. Touch  **Run Last**.
2. (Optional) “Define experiment properties”, including Run File Name (.eds file name), Plate Barcode, Reagent Information, and Data Destination.
3. (Optional) “Edit a method” on page 32.
  - Add, remove, or edit a step, stage, melt curve, or data collection location.
  - Adjust the cover temperature, sample volume, or number of cycles
  - Configure VeriFlex™ Zones, ramp rate, and pause settings.
4. (Optional) “Define plate wells” with sample names and view Well ID, Targets, or Dyes.
5. Load a plate in the instrument (see “Load and unload a plate in the instrument” on page 36).
6. Touch **Start Run**.  
 When prompted, confirm that you inserted a plate.  
**Note:** To disable this reminder, select **Do not show again** or select  **Settings ▶ Insert Plate Reminder** in the home screen.

## Edit an experiment before starting a run

### Define experiment properties

Open an experiment (see “Run an experiment from a saved file” on page 29) or template (see “Create and run an experiment from a template” on page 29) (.edt files), then touch **Properties** ► **Edit** to define the experiment properties.

- Edit the experiment name.
  - a. Touch the **Run File Name** text field.
  - b. Enter the run file name (.eds file), then touch **Done**.
- Add a plate barcode to your experimental record.
  - a. Touch the **Plate Barcode** text field.
  - b. Enter the **Plate Barcode**, then touch **Done**.
- Record reagents and their expiration dates.
  - a. Touch  **Reagent Information**.
  - b. Touch **Add**, or touch an existing reagent, then touch **Edit** or **Delete**.
  - c. Touch the **Name**, **Type**, **Lot #**, **Reagent Barcode**, **Part #**, or **Expiration Date** field to enter individual reagent information.
  - d. Touch **Done**.
- Automatically transfer data after a run is completed.
  - a. Touch  **Data Destination**.
  - b. Touch , or  to choose a data destination.
  - c. Under your desired data destination, select **Automatically transfer experiment**.
  - d. Touch **Done**.  
 Your selection displays on the Properties screen.
- Add a comment or tag to your experiment record.
  - a. Touch **Comments**.
  - b. Enter text, then touch **Done**.

### Scan a barcode using the optional barcode scanner

The instrument is compatible with an optional Handheld Barcode Scanner (Cat. no. 448842, purchased separately). The barcode scanner reads Code 128 (alphanumeric), which supports 128 ASCII character barcodes.

**WARNING! LASER.** Exposure to direct or reflected laser light can burn the retina and leave permanent blind spots. Never look into the laser beam. Remove jewelry and anything else that can reflect the beam into your eyes. Protect others from exposure to the beam.

To scan a barcode using the hand-held barcode scanner:

1. Select the field in the touchscreen or software where you want to enter the barcode.
2. Hold the hand-held barcode scanner 20–30 cm away from a plate or container label and aim at the center of the barcode, then press the trigger.
3. Slowly move the scanning beam across the barcode until the scanner emits a high-pitched tone.

When the scanner scans a barcode, it automatically:

- Transmits the alphanumeric equivalent of the barcode to the touchscreen or software.
- Transmits a carriage-return character (the equivalent of pressing the Enter key).
- Transmits the reagent information for other fields (Lot #, Part #, Expiration Date, etc.)

For more information on the hand-held barcode scanner, refer to the user documentation shipped with the barcode scanner.

## Edit a method

Access a method for editing.

- “Create and run an experiment from a template” on page 29.
- “Run an experiment from a saved file” on page 29.
- “Repeat your last experiment run” on page 30.

See “Method elements” on page 33 for an overview of the method as it is graphically represented in the software.

In the Method screen:


1. Touch **Edit**.
2. Touch in a field, enter any desired change, then touch **Enter**.  
**Note:** Touch and drag  to quickly increase or decrease a step temperature.
3. Touch **Manage Steps** to add or remove a step, stage, melt curve, or data collection location.

**Note:** Manage Steps also provides access to the Advanced Options: **VeriFlex™ Zones**, **Ramp Rates**, and **Add Pause**.

## Method elements

A method consists of one or more stages, each with one or more steps. Within each stage, the number of times the steps repeat is indicated in the Number of Cycles field. The method also specifies the sample volume and heated cover temperature.

See the figure below for a representation of method elements.

## Manage Steps

In the Method screen, touch **Edit ▶ Manage Steps**.

Add or remove a step, stage, melt curve, or data collection point using the **Manage Steps** option. You can also access “Advanced Options” on page 34 from the Manage Steps screen.

- Add a melt curve.
  - a. Touch **Melt curves ▶ Add melt curve**.
  - b. Touch the **+** on the left or right border of a step to add a melt curve before or after, respectively.
  - c. Choose **Continuous** or **Step and hold**.
  - d. Touch melt curve parameters to edit if necessary, then touch **Done**.
  - e. Touch **Done**.

Depending on the experiment type, there may be restrictions on the addition of melt curves.

- Remove a melt curve.
  - a. Touch **Remove melt curve**.
  - b. Touch **—** on the melt curve to be removed.
  - c. Touch **Done**.
- Add a step.
  - a. Touch **Add/Remove steps ▶ Add steps**.
  - b. Touch the **+** on the left or right border of a step to add a step before or after, respectively.
  - c. Enter parameters for the new step, then touch **Enter**.
  - d. Touch **Done**.
- Remove a step.
  - a. Touch **Add/Remove steps ▶ Remove steps**.
  - b. Touch the **—** on any step to remove it, then touch **Done**.
- Add a stage.
  - a. Touch **Add/Remove stages ▶ Add stages**.
  - b. Touch the **+** on the left or right border of a stage to add a stage before or after, respectively.
  - c. Touch **Done**.
  - d. Edit parameters of the new stage in the Method screen (see “Edit a method” on page 32).
- Remove a stage.
  - a. Touch **Add/Remove stages ▶ Remove stages**.
  - b. Touch the **—** on any stage to remove it, then touch **Done**.
- Add or remove data collection locations.
  - a. Touch **Data collection location**.
  - b. Touch  to switch data collection on or off.

### Advanced Options

In the Method screen, touch **Edit ▶ Manage Steps ▶ Advanced Options**.

- To use VeriFlex™ Zones (96-well blocks only):
  - a. Touch **VeriFlex™ Zones**.
  - b. Touch  on the step which you want to apply VeriFlex™ Zones.

- c. Touch each zone to edit the temperature, then touch **Enter**.  
 The background colors of VeriFlex™ Zones change with the temperature edit. Zones with higher temperatures display a background brighter than those with lower temperatures.
- d. (Optional) Apply additional VeriFlex™ Zones.
- e. Touch **Done**.
- To edit Ramp Rates:
  - a. Touch **Ramp Rates**.
  - b. Touch the ramp rates fields.
  - c. Touch **Enter**.
  - d. (Optional) Edit additional steps, then touch **Done**.
- To add a pause into the method:
  - a. Touch **Add Pause**.
  - b. Touch **||** in a stage.
  - c. Enter the pause temperature, and the cycle after which you want the pause to occur.
  - d. Touch **Enter**.  
 The pause is represented by a **P** in the corner of the stage.
  - e. Touch **Done**.

## Define plate wells

Touch the **Plate** tab of an open experiment or template file to access these functions.

Touch  to view the plate layout or  to view the well table.

In the plate layout:


1. Touch **Manage** or touch an individual well.
  - a. Touch the **Samples** tab, select one or more wells, then touch **Edit**.
  - b. Enter a sample name for the well(s) you have selected, then touch **Done**.
  - c. Touch the **Targets** tab to view the targets assigned to the wells.
  - d. In the Targets tab, touch **Details** to view target information, then touch **Done**.
2. Touch **Done** to return to the Plate tab.
3. Touch  to access the well table.
  - a. Touch **Edit** to edit sample names for individual wells.
  - b. Touch a sample name text field, enter a new name, then touch **Done**.

## Load and unload a plate in the instrument

**CAUTION!** Use flat caps for 0.2-mL tubes and 0.1-mL tubes. Rounded caps can damage the heated cover.

1. Load the plate.
  - a. Touch  to eject the instrument drawer.
  - b. Load the plate onto the plate adaptor so that:
    - Well A1 of the plate is in the top-left corner of the plate adaptor.
    - The barcode faces the front of the instrument.

**IMPORTANT!** Plates should be loaded and unloaded by trained operators who have been warned of the moving parts hazard.

**Note:** (For 96-well 0.2-mL blocks only) Do not remove the black plate adaptor before loading a plate or strip tubes. Strip tubes may fit loosely in the adaptor, but when the drawer closes, the heated cover will apply the appropriate pressure to seat the tube strips securely in the adaptor.

**Note:** The 384-well and 96-well Fast (0.1-mL) block configurations do not require a plate adaptor.

- c. Touch  to close the instrument drawer.
2. Unload the plate.
  - a. Touch  to eject the instrument drawer.
  - b. Remove the plate.
  - c. Touch  to close the instrument drawer.

**CAUTION! PHYSICAL INJURY HAZARD.** During instrument operation, the plate temperature can reach 100°C. Allow it to cool to room temperature before handling.

**Note:** If the instrument does not eject the plate, contact Support.

## View, pause, or stop a run

During an instrument run, you can:

- “View well details” on page 37
- “Pause or stop a run” on page 37
- “Adjust the graphical view of an experiment” on page 37
  - Touch  or swipe left once to access real-time views of the run method or to edit the number of cycles.
  - Touch  or swipe left twice to view the real-time data plot.

### View well details

In the home screen, during an instrument run:

1. Touch the right arrow  or swipe left twice.
2. Touch **Well details**.
3. Touch **Samples**, **Targets**, or **Tasks** to select a graphical representation of each.
4. Touch **Close** to return to the home screen.

### Pause or stop a run

In the home screen, during an instrument run:

1. Touch the right arrow  or swipe left to access real-time views of the run.
2. Stop or pause the run.
  - Touch **Stop Run**.
  - Touch **Pause**, then enter a pause temperature.

**CAUTION! PHYSICAL INJURY HAZARD.** During instrument operation, the plate temperature can reach 100°C. If you want to access the plate during a run pause, enter room temperature as the pause temperature and allow the plate to cool to room temperature before handling.

3. (Optional) After pausing a run, touch **Edit** to change the number of cycles.

### Adjust the graphical view of an experiment

In the home screen, during an instrument run:

1. Touch the right arrow  or swipe left twice to access real-time views of the run.
2. Touch **Zoom**.
3. Touch  or  to zoom in or out.
4. Touch the arrows to pan left, right, up, or down on the graph.
5. Touch **Close** to return to the default view.

## Lock the touchscreen during a run

After you have started a run, you can lock the touchscreen so that other users cannot interfere with instrument operation.

**Note:** You must be signed-in to use this feature.

In the home screen:

1. To lock the touchscreen:
  - a. Touch  **My Profile**.
  - b. Touch **Lock Screen**, then touch **Lock**.
2. To unlock the touchscreen:
  - a. Touch anywhere on the touchscreen.
  - b. Touch the **PIN Code** field, and enter your PIN.

**Note:** The touchscreen automatically unlocks when the run is complete.

## Transfer, view, or manage files and results

### Transfer experiment results

In the home screen, when a run is complete:

1. Touch **Transfer File**.
2. Touch  **Cloud**,  **USB**, or  **Network** to select the data destination.
3. Navigate to and select your folder destination.
4. Touch **OK**.
5. Touch **Transfer**.

“View run history” on page 38 to transfer completed run files (.eds files) at any time.

### View run history

In the instrument home screen:

Touch  **Settings** ▶ **Run History**.

- Touch an individual run to view details, then touch **Delete** to delete the record, or **Transfer** to export the experiment .eds file.
- Touch **Manage** to select multiple items for simultaneous viewing, deletion, or transfer.

Guests (users not signed-in) can only view guest records. Users with instrument profiles can also view their own records. Administrators can view all experiment records.

## Manage template (.edt) files

This feature applies to experiment templates (.edt files) on your  **USB** or  **My Instrument**. To manage completed run files (.eds files), see “Transfer experiment results” on page 38.

In the home screen:

1. Touch  **Load Experiment**.
2. Touch  **USB** or  **My Instrument** to select your file location and navigate to your file.
3. Touch **Manage Files**.
4. (Optional) Navigate through the **My Instrument**, **Public**, **USB**, and **Post Read** folders to access an .edt file. Folder availability depends on your sign-in status.
5. Touch to select one or more files.
6. Delete or copy files.
  - Touch **Delete Files** and confirm deletion.
  - Touch **Copy Files**, choose a file destination, then touch **Paste Files**.
7. Touch **Done**.

# Calibrate and verify instrument performance

■ Calibration and verification schedule .....	40
■ Calibration descriptions .....	41
■ View the calibration status and set reminders .....	41
■ Perform ROI/Uniformity, Background, and Dye calibrations .....	42
■ Perform instrument verification using RNase P plates .....	49
■ Calibrate custom dyes .....	54
■ Calibrate for a custom melt curve experiment .....	60

## Calibration and verification schedule

The instrument is factory calibrated and does not require calibration at installation. To ensure optimal performance, perform calibrations at the recommended frequency.

**Note:** After instrument installation, we recommend performing instrument verification using the provided RNase P plate.

---

**IMPORTANT!** Perform calibrations and run experiments under the environmental conditions specified in “Environmental requirements” on page 95. Exposure to extreme temperatures can have adverse effects on the run results, as well as shortening the life span of the instrument components.

---

Calibration type	Recommended frequency
ROI/Uniformity	<ul style="list-style-type: none"><li>Every two years (recommended)</li></ul> <p><b>Note:</b> Set the calibration interval for your instrument: <b>Settings ▶ Maintenance and Service ▶ Calibrations ▶ History and Reminders ▶ Edit ▶ Exp interval field.</b></p> <ul style="list-style-type: none"><li>Background calibration can be performed as needed to check for contamination (depends on usage and laboratory conditions).</li></ul>
Dye	
Background	
RNase P instrument verification	<ul style="list-style-type: none"><li>After installing or moving the instrument</li><li>As needed thereafter to confirm instrument performance</li></ul>

**Note:** For preparing custom dye plates and performing custom calibrations, see “Calibrate custom dyes” on page 54.

## Calibration descriptions

Calibration type	Purpose	Description	Pass criteria
ROI/Uniformity	The software uses calibration data to map the increase in fluorescence to the plate wells during subsequent runs and to evaluate well-to-well consistency of the signals.	The software captures a plate image for each optical filter.	<ul style="list-style-type: none"> <li>The image for each filter distinguishes all wells of the plate.</li> <li>Each well in the image must be distinct and visible at the same luminosity relative to the other wells in the image.</li> </ul>
Background	The software uses calibration data to remove background fluorescence during a run. <b>Note:</b> You can also run this calibration to determine if contamination is related to the sample block or the plate.	The software captures a background image for each optical filter in the absence of sample and reagent, and compares the fluorescence from each well to the average for the plate.	The plate images for all filters are free of abnormal fluorescence.
Dye	The software uses calibration data to characterize and distinguish the individual contribution of each dye in the total fluorescence signals collected by the instrument.	The software extracts a spectral profile for each dye standard, then produces a set of spectral profiles plotted as fluorescence vs filter.	<ul style="list-style-type: none"> <li>Dye spectra peak within the same filter as their group.</li> <li>Dye spectra can diverge slightly at other wavelengths.</li> </ul>

## View the calibration status and set reminders

### View calibration status and set reminders in the instrument

**Note:** The calibration reminders feature requires a connection between the instrument and a computer network.

In the home screen:

1. Touch  **Settings** ▶ **Maintenance and Service** ▶ **Calibrations** ▶ **History and Reminders**.
2. In the Calibration Reminders screen, view the status of each calibration type.
3. (Optional) Touch a calibration row to view the history of that specific calibration type, then touch **Done**.
4. Touch **Edit** to configure the calibration reminder settings. For each calibration type:
  - a. Slide the control **On** or **Off** to enable or disable the calibration reminder.
  - b. Touch the **Exp interval** and **Remind me** text fields to set the calibration time tables.

- c. Touch **Save** to save the settings or **Cancel** to exit the screen without saving.
5. (Optional) Touch **Export** to transfer the calibration report to a  Cloud account,  USB, or  Network drive.
6. Touch **Done** to return to the Calibrations screen.

### View calibration status and set reminders in the Cloud

**Note:** The calibration reminders feature requires a connection between the instrument and a computer network.

1. In the Thermo Fisher Cloud, click  for the Instrument Connect page.
2. Select any of your registered instruments.
3. In the Summary tab:
  - Click **+ Calibrations** to view the status of each calibration type.
  - (Optional) Click **+ Calibration Reminders** to set the calibration reminder time table and enter the notification email address(es).  
**Note:** More than one email address can receive the calibration reminders.  
**Note:** The settings are automatically saved.
  - (Optional) In the Downloads section, click **Maintenance Summary.pdf** to download the calibration status report.
4. In the Calibrations History tab:
  - View the history of each calibration type.
  - (Optional) Click  to download the calibration history report.

## Perform ROI/Uniformity, Background, and Dye calibrations

### Workflow: Calibration

If you start an ROI/Uniformity calibration, the instrument automatically prompts for calibrations to be performed in this order.

Perform an ROI/Uniformity calibration  
 (Always followed by the other calibrations below.)

Perform a Background calibration  
 (Perform any time that the ROI/Uniformity calibration is current.)

Perform Dye calibrations  
 (Perform any time that the ROI/Uniformity and Background calibrations are current.)

## Prepare a calibration plate

### Materials required for calibration plate preparation

- Plate(s) for the calibration you are performing:
  - ROI/Uniformity plate (one ROI plate needed)
  - Background calibration plate
  - Dye calibration plates

**Note:** We recommend calibrating with all Spectral Dye Calibrations Plates available for your block configuration even if you are not using all the dyes in the plates.

**Note:** Do not discard the packaging for the calibration plates. Each calibration plate can be used up to 3 times if the plate is:

- Stored in its packing sleeve at  $-15$  to  $-25^{\circ}\text{C}$
- Used within 6 months after opening
- Used before the plate's expiry date
- Centrifuge with plate adapter; buckets cleaned before use
- Powder-free gloves
- Safety glasses

### Thaw, vortex, and centrifuge a calibration plate


1. Remove the calibration plate from the freezer, then thaw the plate in its packaging. Keep plates protected from light until you perform the calibration.
  - Thaw each plate for 30 minutes.
  - Use each plate within 2 hours of thawing.

---

**IMPORTANT!** Do not remove the plate from its packaging until you are ready to use it. The fluorescent dyes in the wells of calibration plates are photosensitive. Prolonged exposure to light can diminish the fluorescence of the dyes.

---

2. While wearing powder-free gloves, remove the calibration plate from its packaging and retain the packaging. Do not remove the optical film.
3. Vortex the plate for 5 seconds, then centrifuge for 2 minutes at  $750$  to  $1000 \times g$ .
4. Confirm that the liquid in each well is at the bottom of the well and free of bubbles. If it is not, centrifuge the plate again.

---

**IMPORTANT!** Keep the bottom of the plate clean. Fluids and other contaminants on the bottom of the plate can contaminate the sample block and cause an abnormally high background signal.

---

## Perform calibrations

1. In the home screen:

Calibration type	Touch
ROI/Uniformity <sup>[1]</sup>	⚙️ <b>Settings ▶ Maintenance and Service ▶ Calibrations ▶ ROI and Uniformity</b>
Background	⚙️ <b>Settings ▶ Maintenance and Service ▶ Calibrations ▶ Custom ▶ Background</b>
Dye	⚙️ <b>Settings ▶ Maintenance and Service ▶ Calibrations ▶ Dye</b>


<sup>[1]</sup> Automatically followed by Background and Dye calibrations.

2. Follow the instructions on the screen to start the calibration.

**Note:** *Dye calibration only:* Select the Dye Plate to run, then touch **Next**.

3. Load the plate.

- a. Touch  to eject the instrument drawer.
- b. Load the plate onto the plate adaptor so that:
  - Well A1 of the plate is in the top-left corner of the plate adaptor.
  - The barcode faces the front of the instrument.

---

**IMPORTANT!** Plates should be loaded and unloaded by trained operators who have been warned of the moving parts hazard.

---

**Note:** (For 96-well 0.2-mL blocks only) Do not remove the black plate adapter before loading a plate or strip tubes. Strip tubes may fit loosely in the adapter, but when the drawer closes, the heated cover will apply the appropriate pressure to seat the tube strips securely in the adapter.

**Note:** The 384-well and 96-well Fast (0.1-mL) block configurations do not require a plate adapter.

- c. Touch  to close the instrument drawer.

4. Touch **Start**.

5. When the run is complete and the screen displays Calibration Complete, touch **View Results** to check the calibration status.

Calibration status	Action
Passed	Touch <b>Next</b> to proceed to the next required calibration.
Failed	See "Troubleshoot calibration failure" on page 47.

**Note:** You can view the calibration images only after the ROI/Uniformity and Background calibrations pass.

6. Unload the plate.
  - a. Touch  to eject the instrument drawer.
  - b. Remove the plate.
  - c. Touch  to close the instrument drawer.

**CAUTION! PHYSICAL INJURY HAZARD.** During instrument operation, the plate temperature can reach 100°C. Allow it to cool to room temperature before handling.

**Note:** If the instrument does not eject the plate, contact Support.

7. Return the plate to its original packaging.





**Note:** Each calibration plate can be used up to 3 times if the plate is:

- Stored in its packing sleeve at -15 to -25°C
- Used within 6 months after opening
- Used before the plate's expiry date

## View calibration images and transfer results to USB

The instrument performs the ROI, Uniformity, and Background calibrations in sequence. You can view the calibration images after the Background calibration is complete.

1. In the Calibration Status screen, touch **Details**.
2. In the Details screen, touch a calibration type to view its images and plots.

Calibration	Example results indicating successful calibration
ROI  <b>Note:</b> Select the desired filter combination from the Filter Set drop-down list.	Green circles around all wells and bright well centers.  
Uniformity	Signals from each well following a uniform trend.  
Background	Few, if any, signals with abnormally high fluorescence.  
Dye	Signals from each well following a uniform trend.  

3. In the Calibration Status screen:

a. Touch **Accept Results** or **Reject Results**.

Accepting the results saves the calibration data to the instrument and overwrites existing data.

b. (Optional) Touch **Transfer EDS** to transfer the calibration data to a USB.

## Troubleshoot calibration failure

Observation	Possible cause	Recommended action
Calibration failed	The plate was improperly prepared.	Ensure the following: <ul style="list-style-type: none"> <li>The correct plate was used for the calibration performed.</li> <li>The plate was properly thawed.</li> <li>The plate was properly spun down.</li> <li>The heat seal of the plate was properly sealed.</li> </ul>
	The plate is damaged or contaminated.	Check for damage, improper heat seal, or contamination.
		Order a replacement plate. If the replacement plate fails, contact Support.
High fluorescent signal	Signals that exceed the limit of normal fluorescence may indicate fluorescent contaminants on the plate or the sample block.	See "Identify contamination" on page 47.
Calibration failed but plate is undamaged	The incorrect plate was used for calibration performed.	Use the plate that matches the calibration performed.
	The plate was improperly prepared.	Repeat the calibration with the plate properly prepared.
		If the calibration fails again, order a replacement plate. If the replacement plate fails, contact Support.

## Identify contamination

Signals that exceed the limit of normal fluorescence may indicate fluorescent contaminants on the calibration plate or the sample block. Common contaminants include ink residue from permanent pens, powder from disposable gloves, and dust.

To identify and resolve a possible contamination problem:

1. View the calibration data and note the wells that failed the quality check.
2. Remove the plate from the instrument, rotate the plate 180°, then perform the calibration again.

3. Determine the location of the failed wells again as in step 1.

If the position(s) of the failed well(s) are...	Action
Identical	The sample block is contaminated. Decontaminate the sample block (see “Decontaminate the sample block” on page 63).
Reversed	The plate is contaminated. Discard the plate, then perform the calibration using a new calibration plate.

4. If the calibration fails after you decontaminate the sample block and replace the plate, contact Support.

### Create a background plate (optional)

Whenever possible, use a background plate listed in Appendix C, “Parts and materials”. These plates contain a buffer that accurately simulates the reagents used for PCR, and, therefore, produces high-quality calibration data.

If a background plate is not available, you can create one as described below.

#### Required materials

- MicroAmp™ optical 96-well reaction plate
- Optical adhesive cover or optical flat caps
- Pipettor, 200-μL (with pipette tips)
- Powder-free gloves
- Safety glasses
- Deionized water

---

**IMPORTANT!** Wear powder-free gloves while creating the background plate.

---

1. Remove a reaction plate from its box and place it on a clean, dry surface.
2. Aliquot the appropriate volume of deionized water to each well of the reaction plate.  
Recommended volume for a 96-well plate is 10–20 μL per well.
3. Seal the plate using an optical adhesive cover or optical flat caps.
4. Use the plate for background calibration as you would a background plate from the spectral calibration kit.

## Perform instrument verification using RNase P plates

Instruments are factory calibrated, so calibration is not necessary for initial installation. However, we recommend that you verify the instrument performance before using the instrument.

After future calibrations, we highly recommend verifying the instrument performance. The instrument must have valid ROI/Uniformity, Background, and Dye calibrations to perform instrument verification.

### Instrument verification description and schedule

Purpose	Description	Pass criteria	Frequency
Confirms the performance of the instrument.	Quantifies the number of copies of the human RNase P gene in samples with known concentrations of the corresponding DNA template.	$[(C_{TA}) - 3(\sigma_{CTA})] > [(C_{TB}) + 3(\sigma_{CTB})]$ where: <ul style="list-style-type: none"> <li><math>C_{TA}</math> = Average <math>C_T</math> of unknown population A</li> <li><math>\sigma_{CTA}</math> = Standard deviation of unknown population A</li> <li><math>C_{TB}</math> = Average <math>C_T</math> of unknown population B</li> <li><math>\sigma_{CTB}</math> = Standard deviation of unknown population B</li> </ul> The instrument successfully distinguishes between unknown populations A and B with a statistical confidence level of 99.7%	<ul style="list-style-type: none"> <li>After installing or moving the instrument</li> <li>As needed to confirm instrument performance</li> </ul>

### RNase P Instrument Verification Plates

The RNase P plate contains the reagents necessary for the detection and quantitation of genomic copies of the human RNase P gene (a single-copy gene encoding the RNase moiety of the RNase P enzyme). Each well contains: PCR master mix, RNase P primers, FAM™ dye-labeled probe, and a known concentration of human genomic DNA template.



Figure 1 96-Well RNase P plate

- |                             |                     |
|-----------------------------|---------------------|
| ① Unknown A (5000)          | ⑤ STD 5000 copies   |
| ② NTC (no template control) | ⑥ STD 10000 copies  |
| ③ STD 1250 copies           | ⑦ STD 20000 copies  |
| ④ STD 2500 copies           | ⑧ Unknown B (10000) |

**Figure 2** 384-Well RNase P plate

- |                             |                     |
|-----------------------------|---------------------|
| ① NTC (no template control) | ⑤ STD 10000 copies  |
| ② STD 1250 copies           | ⑥ STD 10000 copies  |
| ③ STD 2500 copies           | ⑦ Unknown A (5000)  |
| ④ STD 5000 copies           | ⑧ Unknown B (10000) |

### Analytical performance

After the run, the software calculates average copy number values and standard deviation values. The instrument passes verification if the following inequality is true:

$$[(C_{TA}) - 3(\sigma_{CTA})] > [(C_{TB}) + 3(\sigma_{CTB})]$$

where:

- $C_{TA}$  = Average  $C_T$  of unknown population A
- $\sigma_{CTA}$  = Standard deviation of unknown population A
- $C_{TB}$  = Average  $C_T$  of unknown population B
- $\sigma_{CTB}$  = Standard deviation of unknown population B

### Installation specification

The instrument passes the installation specification if the inequality holds and the instrument successfully distinguishes between unknown populations A and B with a statistical confidence level of 99.7%.

The software automatically adjusts the threshold and omits a defined number of wells from the unknown populations to meet the installation specifications. To view any of the omitted well(s), open the verification file in the desktop software or cloud software.

### Prepare an RNase P plate


#### Materials required for RNase P plate preparation

- RNase P instrument verification plate
- Centrifuge with plate adapter; buckets cleaned before use
- Powder-free gloves
- Safety glasses

#### Thaw, vortex, and centrifuge an RNase P plate

1. Remove the RNase P plate from the freezer, then thaw the plate in its packaging.
  - Thaw the plate for approximately 5 minutes.
  - Use the plate within 30 minutes of thawing.
2. While wearing powder-free gloves, remove the plate from its packaging.


3. Vortex the plate for 5 seconds, then centrifuge for 2 minutes at 750 to 1000 × g.
4. Confirm that the liquid in each well is at the bottom of the well and free of bubbles. If it is not, centrifuge the plate again.

**IMPORTANT!** Keep the bottom of the plate clean. Fluids and other contaminants on the bottom of the plate can contaminate the sample block and cause an abnormally high background signal.

## Perform RNase P verification

1. In the home screen, touch **Settings** ▶ **Maintenance and Service** ▶ **RNase P Verification**.
2. Load the plate.
  - a. Touch to eject the instrument drawer.
  - b. Load the plate onto the plate adaptor so that:
    - Well A1 of the plate is in the top-left corner of the plate adaptor.
    - The barcode faces the front of the instrument.

**IMPORTANT!** Plates should be loaded and unloaded by trained operators who have been warned of the moving parts hazard.

**Note:** (For 96-well 0.2-mL blocks only) Do not remove the black plate adapter before loading a plate or strip tubes. Strip tubes may fit loosely in the adapter, but when the drawer closes, the heated cover will apply the appropriate pressure to seat the tube strips securely in the adapter.

**Note:** The 384-well and 96-well Fast (0.1-mL) block configurations do not require a plate adapter.

- c. Touch to close the instrument drawer.

3. Touch **Start**.
4. When the run is complete and the screen displays Verification Complete, touch **View Results** to confirm the status of the run.

Calibration status	Action
Passed	Instrument is ready for use.
Failed	See "Troubleshoot verification failure" on page 53.

5. In the RNase P Verification Status screen, touch:
  - **Accept Results** to save the results to the instrument
  - **Reject Results** to delete the RNase P verification results
  - **Export Results** to export the calibration results to a USB
6. Unload the plate.
  - a. Touch  to eject the instrument drawer.
  - b. Remove the plate.
  - c. Touch  to close the instrument drawer.

**CAUTION! PHYSICAL INJURY HAZARD.** During instrument operation, the plate temperature can reach 100°C. Allow it to cool to room temperature before handling.

**Note:** If the instrument does not eject the plate, contact Support.

## Troubleshoot verification failure

Observation	Possible cause	Recommended action
Verification failed	The plate was improperly prepared.	Ensure the following: <ul style="list-style-type: none"> <li>The correct plate was used for the verification performed.</li> <li>The plate was properly thawed.</li> <li>The plate was properly spun down.</li> <li>The heat seal of the plate was properly sealed.</li> </ul>
		Open the verification file in the desktop software or cloud software to view the flags and troubleshooting details for failed wells.
	The plate is damaged or contaminated.	Check for damage, improper heat seal, or contamination.
		Order a replacement plate. If the replacement plate fails, contact Support.
High fluorescent signal	Signals that exceed the limit of normal fluorescence may indicate fluorescent contaminants on the plate or the sample block.	See “Identify contamination” on page 47.
Verification failed but plate is undamaged	The incorrect plate was used for verification.	Use the correct RNase P plate for verification.
	The plate was improperly prepared.	Repeat the verification with a new properly prepared plate. <b>Note:</b> The verification procedure is an experiment run, so each RNase P plate can only be used once.
		Open the verification file in the desktop software or cloud software to view the flags and troubleshooting details for failed wells.
		If the verification fails again, order a replacement plate. If the replacement plate fails, contact Support.

## Calibrate custom dyes

### Custom dyes overview

The QuantStudio™ 3 and 5 Real-Time PCR Systems can be used to run assays designed with custom dyes (dyes not manufactured by Thermo Fisher Scientific or dyes not pre-calibrated with the instrument). Custom dyes must excite between 455–672 nm and emit between 505–723 nm.

### Custom dye calibration

For each custom dye, determine the optimal dye concentration. Use this concentration for preparing all subsequent dye calibration plates.

#### Dilute the custom dye to an optimal concentration

Prepare a custom dye dilution plate

Run the dilution plate as an experiment

Determine the optimal dye concentration

#### Calibrate the custom dye for each plate

Create a custom dye calibration plate

Add a custom dye to the software

Perform a custom dye calibration

### Dilute the custom dye to an optimal concentration

#### Custom dye dilution guidelines

Prepare a dilution series for each custom dye.

- Target several dye concentrations within a range of 100–2000 nM.
- Choose a 2- or 3-fold difference in dilution points.
- Dispense 10–20 µL per well for 96-well plates or 5 µL per well for 384-well plates.
- Dilute the dye in buffer compatible with your master mix.
- *(Intercalating dyes only)* Add the appropriate amount of amplified PCR product to generate fluorescence.

#### Prepare a custom dye dilution plate



---

**IMPORTANT!** Wear powder-free gloves throughout the procedure.


---

1. Prepare a 2 or 3-fold dilution series of the custom dye.
2. Dispense aliquots of each dilution into the center of a reaction plate, then seal the plate.

A full plate is not needed. See the figures below for suggested replicates.

3. Vortex the plate for 5 seconds, then centrifuge for 2 minutes at 750 to 1000  $\times$  g.
4. Confirm that the liquid in each well is at the bottom of the well and free of bubbles. If it is not, centrifuge the plate again.

**IMPORTANT!** Keep the bottom of the plate clean. Fluids and other contaminants on the bottom of the plate can contaminate the sample block and cause an abnormally high background signal.

## Run the dilution plate as an experiment

1. Load the plate into the instrument.
2. Set up a genotyping experiment on the instrument.
  - Load an experiment created in either the desktop or cloud software.
    - a. In the home screen, touch **Load Experiment**.
    - b. Select the experiment location, then the experiment file.
  - Create a new experiment on the instrument touchscreen.
    - a. In the home screen, touch **Open Template ▸ Genotyping ▸ Genotyping Post**.
    - b. (Optional) In the Properties tab, edit the experiment properties.
    - c. In the Method tab, set the hold temperature to 60°C with a 2 minute hold and enter the appropriate reaction volume.
    - d. In the Plate tab, enter the dilution series information for the appropriate wells.
3. Touch **Start Run**.
4. When the run is complete, download the results for analysis.
5. Unload the plate from the instrument.

**CAUTION! PHYSICAL INJURY HAZARD.** During instrument operation, the plate temperature can reach 100°C. Allow it to cool to room temperature before handling.

---

## Determine the optimal dye concentration

Review the dye signal data and select the dilution to calibrate.

1. In the **Results** tab of either the desktop or cloud software, select **Raw Data Plot**.

This plot displays the raw fluorescent signal of each dye, for individual wells and at individual cycles over the duration of the PCR run.

**Note:** The Raw Data Plot is not available for viewing on the instrument touchscreen.
2. For each replicate population of dilutions, select the wells in the Plate Layout tab to view in the plot.

- Examine the raw data and identify the well(s) yielding signals according to the ranges shown in the following table.

Plate type	Acceptable signal range
96-well	800,000 to 3,200,000
384-well	400,000 to 2,000,000

**Note:** You can also export the raw data and average for the various concentrations.

- Select the lowest (optimal) dye concentration that falls within the acceptable signal range.


## Calibrate the custom dye

### Create a custom dye calibration plate

**IMPORTANT!** Wear powder-free gloves while creating the dye plate.

Create a full plate of the custom dye diluted to the optimal concentration:

- Dilute the custom dye to the optimal concentration in buffer. Prepare an adequate volume, using a recommended volume range of 10–20  $\mu\text{L}$ /well for a 96-well plate.
- Pipet the appropriate volume of the diluted custom dye to all wells of an optical reaction plate.
- Seal the plate.
- Vortex the plate for 5 seconds, then centrifuge for 2 minutes at 750 to 1000  $\times$  g.
- Confirm that the liquid in each well is at the bottom of the well and free of bubbles. If it is not, centrifuge the plate again.

**IMPORTANT!** Keep the bottom of the plate clean. Fluids and other contaminants on the bottom of the plate can contaminate the sample block and cause an abnormally high background signal.

### Add a custom dye to the software

- In the home screen, touch **Settings** **► Maintenance and Service** **► Calibrations** **► Custom** **► Custom Dye**.
- Touch **Add Custom Dye** to enter information about the new custom dye.

3. Enter the dye information:

Field/option	Action
Custom Dye Name	Enter a name for the custom dye. <b>IMPORTANT!</b> <ul style="list-style-type: none"> <li>• Do not use a system dye name for a custom dye name.</li> <li>• Dye names are case and spacing sensitive.</li> </ul>
Type	Select: <ul style="list-style-type: none"> <li>• <b>Reporter</b> – The dye works in conjunction with a quencher dye to report an increase of PCR product.</li> <li>• <b>Quencher</b> – The dye suppresses the fluorescence of a reporter dye until amplification of PCR product.</li> <li>• <b>Both</b> – The dye reports an increase of PCR product without the aid of a quencher dye.</li> </ul>

4. Touch **Save** to add the custom dye or **Cancel** to exit the screen without saving the new dye details.

**Note:** You must add the custom dye to the desktop or cloud software dye libraries before creating, running, or analyzing experiments that use custom dyes. Refer to the *QuantStudio™ Design and Analysis desktop Software User Guide* (Pub. no. MAN0010408) or the *QuantStudio™ Design and Analysis cloud Software Help* (Pub. no. MAN0010414).

### Perform a custom dye calibration


**IMPORTANT!** If you added a custom dye to the dye library in the desktop software or cloud software, you will need to re-enter the custom dye information in the instrument touchscreen before you perform the custom dye calibration (see “Add a custom dye to the software” on page 57).

1. Load the plate (see “Load and unload a plate in the instrument” on page 36).
2. In the home screen, touch  **Settings** ▶ **Maintenance and Service** ▶ **Calibrations** ▶ **Custom** ▶ **Custom Dye**.
3. Touch the custom dye you wish to calibrate.
4. Review the custom dye information, make changes as necessary, then touch **Update**.
5. Enter the calibration temperature.
6. (Optional) Touch **Reagents**, then enter reagent information.
7. Touch **Start**.
8. When the run is complete and the screen displays Calibration Complete, touch **View Results** ▶ **Details**.

9. Review the plot. Passing calibration results show uniform signals with peaks aligned with the dye's wavelength.

Peak channel	Filter wavelength (nm) <sup>[1]</sup>	
	Excitation	Emission
x1-m1	470 ± 15	520 ± 15
x2-m2	520 ± 10	558 ± 12
x3-m3	550 ± 10	587 ± 10
x4-m4	580 ± 10	623 ± 14
x5-m5	640 ± 10	682 ± 14
x6-m6	662 ± 10	711 ± 12

<sup>[1]</sup> The central wavelengths are the optimized wavelengths.

**Note:** Example dye calibration plot.  
The peaks for your dye may align with a different filter set.

10. Select an action depending on whether the custom dye calibration passed or failed.

Calibration status	Action
Passed	<ul style="list-style-type: none"> <li>Touch <b>Accept Results</b> or <b>Reject Results</b>.</li> </ul> <p><b>Note:</b> Accepting the results saves the calibration data to the instrument and overwrites existing data.</p> <ul style="list-style-type: none"> <li>[Optional] Touch <b>Transfer EDS</b> to transfer the calibration data to a USB.</li> </ul>
Failed	<ul style="list-style-type: none"> <li>Create another custom dye plate using the next dye concentration greater than the concentration determined in “Determine the optimal dye concentration”, then perform the calibration again.</li> <li>See “Troubleshoot calibration failure” on page 47.</li> </ul>

11. Unload the plate (see “Load and unload a plate in the instrument” on page 36).

## Calibrate for a custom melt curve experiment

**Note:** A custom melt calibration calibrates a custom dye and a melt calibration at the same time.

Before running the custom melt calibration:

- Ensure that all calibrations are current (touch  **Settings** ▶ **Maintenance and Service** ▶ **Calibrations** ▶ **History and Reminders**).
- “Add a custom dye to the software” on page 57.

In the home screen:

1. Touch  **Settings** ▶ **Maintenance and Service** ▶ **Calibrations** ▶ **Custom** ▶ **Custom Melt**.
2. Touch **PCR + Melt** or **Melt only** as appropriate for the kit you are using.
3. Load the plate (see “Load and unload a plate in the instrument” on page 36).
4. Select or add a dye, then select a filter set appropriate for your dye's wavelength (see filter-wavelength table below).

**Note:** Refer to your reagent kit documentation for dye name and wavelength information.

Peak channel	Filter wavelength (nm) <sup>[1]</sup>	
	Excitation	Emission
x1-m1	470 ± 15	520 ± 15
x2-m2	520 ± 10	558 ± 12
x3-m3	550 ± 10	587 ± 10
x4-m4	580 ± 10	623 ± 14
x5-m5	640 ± 10	682 ± 14
x6-m6	662 ± 10	711 ± 12

<sup>[1]</sup> The central wavelengths are the optimized wavelengths.


---

**IMPORTANT!** If the selected filter set does not match your reagent kit documentation, then the incorrect wavelength may be collected during a run.

---

5. (Optional) Touch **Reagents**, then enter reagent information.
6. Touch **Start**.
7. When the run is complete and the screen displays Calibration Complete, touch **View Results** ▶ **Details**.

8. Review the plot. Passing calibration results show uniform signals with peaks aligned with the dye's wavelength.

**Note:** Example dye calibration plot.  
The peaks for your dye may align with a different filter set.

9. Select an action depending on whether the custom dye calibration passed or failed.

Calibration status	Action
Passed	<ul style="list-style-type: none"> <li>Touch <b>Accept Results</b> or <b>Reject Results</b>.</li> </ul> <p><b>Note:</b> Accepting the results saves the calibration data to the instrument and overwrites existing data.</p> <ul style="list-style-type: none"> <li>(Optional) Touch <b>Transfer EDS</b> to transfer the calibration data to a USB.</li> </ul>
Failed	<ul style="list-style-type: none"> <li>Create another custom dye plate using the next dye concentration greater than the concentration determined in "Determine the optimal dye concentration", then perform the calibration again.</li> <li>See "Troubleshoot calibration failure" on page 47.</li> </ul>

10. Unload the plate (see "Load and unload a plate in the instrument" on page 36).

**Note:** You must add the custom dye to the desktop or cloud software dye libraries before creating, running, or analyzing experiments that use custom dyes. Refer to the *QuantStudio™ Design and Analysis desktop Software User Guide* (Pub. no. MAN0010408) or the *QuantStudio™ Design and Analysis cloud Software Help* (Pub. no. MAN0010414).

**Note:** To perform a custom melt experiment, you can either create a Standard Curve or a Custom experiment with melt, then specify the data points per degree in the method.

# Maintain the instrument

■ Backup or restore the instrument .....	62
■ Decontaminate the sample block .....	63
■ Replace the instrument fuses .....	65
■ Power on or off, store, and move .....	66

---

**IMPORTANT!** This chapter contains user maintenance procedures for the instrument. Procedures other than those described in this document must be performed by a qualified Service Engineer.

---

## Backup or restore the instrument

In the home screen:

Touch  **Settings** ▶ **Maintenance and Service** ▶ **Backup/Restore**.

To	Action
Backup	<ol style="list-style-type: none"> <li>1. Touch <b>Backup Instrument</b>.</li> <li>2. Touch  <b>USB</b> or  <b>Cloud</b> to choose where to backup your instrument.</li> <li>3. Enter a backup file name, then touch <b>Done</b>.</li> <li>4. Select which elements you will backup, or leave them all selected.</li> <li>5. Touch <b>Backup</b>.</li> </ol>
Restore (Administrator only)	<ol style="list-style-type: none"> <li>1. Touch <b>Restore a Backup</b>.</li> <li>2. Touch  <b>USB</b> or  <b>Cloud</b> to select where you will restore from.</li> <li>3. Select your backup file, then touch <b>Restore</b>.</li> </ol>

## Decontaminate the sample block

Perform this procedure to eliminate fluorescent contaminants from the instrument sample block. Contamination is generally evident in failed background calibrations where one or more wells consistently exhibit abnormally high signals.

**CAUTION! PHYSICAL INJURY HAZARD.** Do not remove the instrument cover. There are no components inside the instrument that you can safely service yourself. If you suspect a problem, contact Support.

**CAUTION! PHYSICAL INJURY HAZARD.** During instrument operation, the sample block temperature can reach 100°C. Allow it to cool to room temperature before handling.

**CAUTION!** Before using a cleaning or decontamination method other than those recommended by Thermo Fisher Scientific, confirm with Thermo Fisher Scientific that the proposed method will not damage the equipment.

### Materials required

- Safety glasses
- Powder-free gloves
- Tissue, lint-free
- Cotton or nylon swabs and lint-free cloths
- Pipette (100-μL) with pipette tips
- Deionized water
- Ethanol, 95% solution
- Bleach, 10% solution

### Clean the sample block


**CAUTION! PHYSICAL INJURY HAZARD.** During instrument operation, the sample block temperature can reach 100°C. Allow it to cool to room temperature before handling.

**IMPORTANT!** Wear powder-free gloves when you perform this procedure.

1. Identify the contaminated wells of the sample block (see “Identify contamination” on page 47).
2. Prepare the instrument and access the sample block:
  - a. Power off and unplug the instrument, then allow it to cool for 15 minutes.
  - b. Pull the instrument drawer forward to expose the sample block.

3. Rinse the contaminated wells of the sample block with deionized water (see “Solvents for cleaning the sample block” on page 65).

4. Close the drawer and test the block for contamination:
  - a. Push the instrument drawer back in to the instrument.
  - b. Plug in, then power on the instrument.
  - c. Perform a background calibration to confirm that you have eliminated the contamination.
5. If the contamination remains:
  - a. Repeat step 2 – step 3.
  - b. Clean the contaminated wells of the sample block using a 95% ethanol solution (see “Solvents for cleaning the sample block” on page 65).
  - c. Repeat step 3 – step 4 to rinse the sample block and to confirm that you have eliminated the contamination.

---

**IMPORTANT!** Always use deionized water to rinse wells after cleaning with bleach or ethanol solution.

---

6. If the contamination still remains:
  - a. Repeat step 2 – step 3.
  - b. Clean the contaminated wells of the sample block using a 10% bleach solution (see “Solvents for cleaning the sample block” on page 65).
  - c. Repeat step 3 – step 4 to rinse the sample block and to confirm that you have eliminated the contamination.

---

**IMPORTANT!** Always use deionized water to rinse wells after cleaning with bleach or ethanol solution.

---

7. If the contamination continues to remain, contact Support.

## Solvents for cleaning the sample block

**IMPORTANT!** Use these solvent cleaning procedures *only* in conjunction with the "Clean the sample block" procedures.

### Rinse the sample block with deionized water

1. Pipet a small volume of deionized water into each contaminated well.
2. In each well, pipet the water up and down several times to rinse the well.
3. Pipet the water to a waste beaker.
4. Use a cotton swab to scrub inside of each contaminated well.
5. Use a lint-free cloth to absorb the excess deionized water.

### Clean the sample block with 95% ethanol

1. Pipet a small volume of 95% ethanol solution into each contaminated well.
2. In each well, pipet the solution up and down several times to rinse the well.
3. Pipet the ethanol solution to a waste beaker.

**IMPORTANT!** Always use deionized water to rinse wells after cleaning with bleach or ethanol solution.

### Clean the sample block with 10% bleach

1. Pipet a small volume of 10% bleach solution into each contaminated well.
2. In each well, pipet the solution up and down several times to rinse the well.
3. Pipet the bleach solution to a waste beaker.

**IMPORTANT!** Always use deionized water to rinse wells after cleaning with bleach or ethanol solution.

## Replace the instrument fuses

**CAUTION! FIRE HAZARD.** For continued protection against the risk of fire, replace fuses only with listed and certified fuses of the same type and rating as those currently in the instrument.

### Materials required

- Fuses (2) – 10A, Time-Lag T, 250VAC, 5 × 20mm
- Safety glasses
- Powder-free gloves
- Screwdriver, flathead

## Replace the fuses

1. Power off and unplug the instrument, then allow it to cool for 15 minutes.
2. Using a flat-head screwdriver, unscrew and remove the fuse holder.

3. Remove each fuse from its fuse holder and inspect it for damage. Carbon typically coats the inside of failed fuses.

Good	Failed
	

4. Replace each failed fuse.  
**Note:** The voltage and amperage ratings are on the fuse holder.
5. Install the fuse holder back into the instrument.
6. Plug in, then power on the instrument.  
The installation is successful if the instrument powers on.

**Note:** Fuse failure can result from fluctuations in the supplied power to the system. To prevent further failures, consider installing an electrical protective device, such as a UPS or a surge protector. If issues with the fuse persist, contact Support.

## Power on or off, store, and move

### Enable sleep mode


In the home screen:

1. Touch  **Settings** ▶ **Instrument Settings** ▶ **Sleep Mode**.
2. Slide the control **On** or **Off**.
3. Touch **Edit Time**, then enter the period of inactivity until the software goes into sleep mode.
4. Touch **Enter**.
5. Touch **OK**.

## Power on the instrument

To power on the instrument from a powered-off state:

1. Touch anywhere on the touchscreen to determine if the instrument is in sleep mode. If the home screen displays, the instrument is already powered on.
2. If the home screen does not display, power on the instrument by pressing the switch on the rear panel.

The instrument is ready to use when the home screen is displayed.

If left unattended (for about two hours), the instrument automatically enters sleep mode (enabled by default) to conserve power. Refer to the touchscreen Help system for step-by-step instructions for changing the sleep mode setting.

## Power off the instrument

The instrument operates in low-power mode when not in use. However, the instrument can be powered off completely so that the components draw no power.

**Note:** If you will be shutting down the instrument for >1 week, see “Prepare the instrument to ship, move, or store” on page 67.

1. Power off the instrument from the power switch on the rear of the instrument.
2. Power off the computer.

## Prepare the instrument to ship, move, or store

In the home screen:

1. Touch  **Settings** ▶ **Maintenance and Service** ▶ **Ship Prep Mode** ▶ **Next**.
2. Touch  to eject the instrument drawer.
3. Load the packing plate or an empty plate, then touch  to close the drawer.
4. Touch **Lock Block**.
5. Power off the instrument using the power switch on the back of the instrument.

The instrument is now ready to ship, move, or store.

## Move the instrument

**CAUTION! PHYSICAL INJURY HAZARD.** Do not attempt to lift the instrument or any other heavy objects unless you have received related training. Incorrect lifting can cause painful and sometimes permanent back injury. Use proper lifting techniques when lifting or moving the instrument. At least two people are required to lift it.

**IMPORTANT!** Moving your instrument can create subtle changes in the alignment of the instrument optics. Recalibrate the instrument if necessary.

- Ensure that the surface on which you place the instrument can support at least 35 kg (77 lbs).
- Ensure that the path to transport the instrument is clear of obstructions.
- At least two people are needed to lift and carry the instrument.
- Keep your spine in a good neutral position.
- Bend at the knees and lift with your legs.
- Do not lift an object and twist your torso at the same time.
- Coordinate your intentions with your assistant before lifting and carrying.

**IMPORTANT!** After moving the instrument, perform an RNase P instrument verification run. If the run fails, perform ROI/ uniformity, background, and dye calibrations.

## Return the instrument for service

The service process requires 2 to 3 weeks.

Before returning the instrument for service:

1. Back up the instrument (see “Backup or restore the instrument” on page 62).
2. In the home screen, touch  **Settings** ▶ **Instrument Settings** ▶ **Reset Factory defaults** to unlink from all Cloud accounts.
3. Set the instrument to Ship Prep Mode (see “Prepare the instrument to ship, move, or store” on page 67).

To return the instrument for service:

1. Contact your local customer care center or technical support group to obtain a copy of the Certificate of Instrument Decontamination, a service notification, a service call number, and packaging materials (if required).
2. Follow the instructions in the form to decontaminate the instrument.

**IMPORTANT!** The instrument must be decontaminated before packing it for shipping.

3. Complete and sign a copy of the Certificate of Instrument Decontamination.
4. Fax the Certificate of Instrument Decontamination to the customer care center.

5. Pack the instrument in the provided packaging and follow the instructions in the table below.

Prepare	Include	Exclude
<ol style="list-style-type: none"> <li>1. Transfer any data files from the instrument.</li> <li>2. Load an empty plate in the sample block.</li> <li>3. Use the touchscreen to place the instrument in ship mode.</li> </ol> <p><b>Note:</b> The empty plate and ship mode protects the internal components of the instrument during transport.</p>	<ul style="list-style-type: none"> <li>• Instrument</li> <li>• Completed and signed Certificate of Instrument Decontamination</li> </ul> <p><b>Note:</b> The instrument will not be accepted for service without a hard copy of the Certificate of Instrument Decontamination.</p>	<p>Any accessories, including:</p> <ul style="list-style-type: none"> <li>• Power cord</li> <li>• Ethernet cable</li> <li>• USB drive</li> <li>• Wireless adapter</li> </ul> <p><b>Note:</b> If included with the instrument, these items will be disposed of during service and not returned.</p>

6. Attach the postage provided with the Certificate of Instrument Decontamination to the box, then ship the instrument to the designated facility.

# Profile and instrument configuration tasks

■ Initial start up .....	70
■ Manage profiles .....	71
■ Enable SAE mode (Administrator only) .....	72
■ Require sign-in (Administrator only) .....	73
■ Enable remote instrument monitoring (Administrator only) .....	73
■ Update instrument software (Administrator only) .....	73
■ Restore factory defaults .....	74
■ Configure the network .....	74
■ Select a Cloud region (Administrator only) .....	75
■ Manage the instrument name (Administrator only) .....	75
■ Set the date and time .....	76
■ Manage the Sign Out Timer (Administrator only) .....	76

## Initial start up

Perform these tasks during the initial start up or after you “Restore factory defaults” on page 74.

1. (Optional) “Configure the network” on page 74, or touch **Neither/Decide Later**.
2. (Optional) “Manage the instrument name (Administrator only)” on page 75.
3. “Set the date and time” on page 76.
4. “Select a Cloud region (Administrator only)” on page 75, or touch **Next** to continue.
5. “Create an administrator profile” on page 71.
6. “Link your instrument profile to Thermo Fisher Cloud” on page 24 or touch **Skip**.

## Manage profiles

To	See
Create a profile	<ul style="list-style-type: none"> <li>• “Create an administrator profile” on page 71</li> <li>• “Create a new instrument profile” on page 23</li> </ul>
Configure a profile	<ul style="list-style-type: none"> <li>• “Link your instrument profile to your Thermo Fisher Cloud account” on page 24</li> <li>• “Edit a user profile” on page 72</li> </ul>
View or manage all profiles	<ul style="list-style-type: none"> <li>• “View all user profiles (Administrator only)” on page 71</li> <li>• “Manage all instrument profiles (Administrator only)” on page 71</li> </ul>

### Create an administrator profile

“Initial start up” on page 70 of the instrument automatically prompts for the creation of an administrator profile. The first profile created during instrument startup is given administrator privileges. Administrators can assign additional users administrative privileges (see “Manage all instrument profiles (Administrator only)” on page 71).

1. Touch **Name**, then enter a username and touch **Done**.
2. Touch **PIN**, then enter a four-digit numerical password and touch **Enter**.  
**Note:** Touch the **Show PIN** checkbox to switch PIN display on or off.
3. Touch **Confirm PIN** and repeat the previous step.
4. Touch **Create profile**.

### View all user profiles (Administrator only)

In the home screen:

1. Navigate to a list of all profiles.
  - Touch  **My Profile**, then touch the **All Profiles** tab.
  - Touch  **Settings** ▶ **Manage Users** ▶ **Manage Profiles**.

A list of Users, the date their profile was created, and User type is displayed.

2. Touch **Done**.

### Manage all instrument profiles (Administrator only)

In the home screen:

1. Touch  **My Profile**, then touch the **All Profiles** tab.  
Alternatively, touch  **Settings** ▶ **Manage Users** ▶ **Manage Profiles**, then touch the **All Profiles** tab.
2. Select the instrument profile to edit.

3. Edit the profile.
  - To delete the profile, touch **Delete profile** ▶ **Delete**.
  - To reset the PIN, touch **Reset PIN** ▶ **Reset**.
    - The PIN will be deleted, and the user will be directed to enter a new PIN upon the next sign in.
  - To enable or disable administrative privileges, slide the control to Administrator or Standard.
4. Touch **Done**.

## Edit a user profile

“Sign In” on page 23 to access these functions. “Manage all instrument profiles (Administrator only)” on page 71 from this screen if you are an administrator.

In the home screen:

1. Touch  **My Profile**.  
 Administrators can also navigate to this screen by touching  **Settings** ▶ **Manage Users** ▶ **Manage Profiles**.
2. Touch **Edit**.
3. Select the fields to edit, then make changes.
4. Touch **Done**.

## Enable SAE mode (Administrator only)

**Note:** The Security, Audit, and e-Signature (SAE) module is available for QuantStudio™ 5 Systems only.

The SAE module works with the instrument and desktop software to record and restrict various user activities. Audits and e-Signatures are tracked within the desktop software.

In the home screen:

1. Touch  **Settings** ▶ **Instrument Settings** ▶ **SAE Mode**.
2. Slide the control **On** or **Off**.
3. Touch **OK**.

The following functions are disabled in SAE mode.

- Start a run from the instrument.
- Edit an instrument name.
- Edit an instrument date or time.
- Restore to factory defaults.
- Perform a software update.

## Require sign-in (Administrator only)

In the home screen:

1. Touch  **Settings** ▶ **Manage Users** ▶ **Sign In Required**.
2. Slide the control **On** or **Off**.
3. Touch **Done**.

Requiring a user to sign in disables guest profiles.

## Enable remote instrument monitoring (Administrator only)

In the home screen:

1. Touch  **Settings** ▶ **Maintenance and Service** ▶ **Monitoring**.
2. Configure your instrument monitoring settings.

Option	Description
Remote Monitoring Service	Slide the control <b>On</b> or <b>Off</b> to enable or disable the instrument to automatically notify Thermo Fisher Scientific support teams in real time of potential instrument issues. Data and experiments are not monitored.
Thermo Fisher Cloud Monitor	Slide the control <b>On</b> or <b>Off</b> to enable or disable real-time monitoring of amplification plots from your Thermo Fisher Cloud dashboard.

3. Touch **OK**.

## Update instrument software (Administrator only)

In the home screen:

1. Touch  **Settings** ▶ **Maintenance and Service** ▶ **Software Update**.
2. Touch  **USB** or  **Cloud** to select the location of the update files.
3. When prompted, confirm your request to update the software.

## Restore factory defaults

---

**IMPORTANT!** Back up the instrument before restoring factory defaults (see “Backup or restore the instrument” on page 62).

---

In the home screen:

1. Touch  **Settings** ▶ **Instrument Settings** ▶ **Restore Factory Defaults**.
2. Touch **Restore Factory Defaults**.
3. Power **Off**, then power **On** the instrument to effect the change in settings.  
All profile and personal information will be permanently removed from the instrument.

Following a factory default restore, follow the procedures described in “Initial start up” on page 70.

## Configure the network

What do you want to do?

 “Set up a wired connection” on page 74

 “Set up a wireless connection” on page 75

### Set up a wired connection

In the home screen:

1. Touch  **Settings** ▶ **Instrument Settings** ▶ **Network Connection** ▶  **Wired**.
2. Touch a radio button to choose to connect to a network by DHCP or a Static IP address on the Network Configuration screen.

Option	Action
<b>Automatic</b>	Touch <b>DHCP</b>
<b>Manual</b>	<ol style="list-style-type: none"> <li>1. Touch <b>Static IP</b></li> <li>2. Enter the appropriate IP addresses for the instrument, the Subnet Mask, and, optionally, the Default Gateway, the Primary DNS Server, and the Secondary DNS Server using the numeric editor. Addresses are in the form of X.X.X.X, where each X is a 3-digit number, from 001 to 255.</li> </ol>

**Note:** Ask your system administrator if the IP address is assigned statically or dynamically. For static addresses, you need to know the IP address for the instrument, the subnet mask, and the default gateway.

3. Touch **OK** to save the changes and return to the Network Connection screen.

“Create an administrator profile” on page 71 to continue the initial start up of the instrument.

## Set up a wireless connection

In the home screen:

1. Touch  **Settings** ▶ **Instrument Settings** ▶ **Network Connection** ▶  **Wireless** to display a list of the available networks.
2. Select a network or touch **Join other network**, then enter the network password and touch **Enter**.  
If you choose Join other network, enter the network name and security type.
3. Touch **Join** to continue or **Cancel** to exit.
4. In the Network Connection Complete screen, touch **Next** to continue.
5. Touch **OK**.
6. Edit the Network details in the Network Configuration screen or touch **Done**.

Proceed to “Create an administrator profile” on page 71 to continue the initial start up of the instrument.

## Select a Cloud region (Administrator only)

During the initial startup, you will be automatically prompted to select a Cloud region. Return to the selection through  **Settings**, if necessary.

In the home screen:

1. Touch  **Settings** ▶ **Instrument Settings** ▶ **Cloud Region**.
2. Touch a region.  
The server used to store your data is located in this region.
3. Touch **OK**.

## Manage the instrument name (Administrator only)

In the home screen:

1. Touch  **Settings** ▶ **Instrument Settings** ▶ **Instrument Name**.
2. Touch the Instrument Name field, enter an instrument name, then touch **Done**.
3. (Optional) Touch **Add Avatar** to associate an avatar with your instrument. You must insert a USB drive with images to use this option.
4. Touch **OK**.

## Set the date and time

In the home screen:

1. Touch  **Settings** ▶ **Instrument Settings** ▶ **Date/Time**.
2. Select a time zone from the drop-down list.
3. Select a date format.
  - a. Touch **Date Format**.
  - b. Select the radio button of your preferred date format.
  - c. Touch **Next**, touch the date field, and enter the date.
  - d. Touch **Enter**, then touch **Done**.
4. Select a time format.
  - a. Touch **Time Format**.
  - b. Slide the control left or right to choose a 12-hour or 24-hour clock.
  - c. Touch **Next**, touch the time field, and enter the time.
  - d. Touch **Enter**, then touch **Done**.
5. Touch **Done**.

## Manage the Sign Out Timer (Administrator only)

In the home screen:

1. Touch  **Settings** ▶ **Manage Users** ▶ **Sign Out Timer**.
2. Touch the **Edit Time** field, then enter the desired duration of inactivity before automatic user sign out.
3. Touch **Enter**, then touch **Done**.

# Install and connect the instrument to a network

■ Workflow: Install and connect to a network .....	77
■ Before you begin .....	78
■ Unpack and install the instrument .....	78
■ Power on and follow the startup wizard .....	79
■ Connect the instrument and the computer directly or to a LAN .....	79
■ Instrument and computer connections .....	81
■ Download and install the QuantStudio™ Design and Analysis desktop Software .....	82
■ Networking .....	84

## Workflow: Install and connect to a network

Perform all steps in “Before you begin” on page 78

“Unpack and install the instrument” on page 78

“Power on and follow the startup wizard” on page 79

“Perform instrument verification using RNase P plates” on page 49

“Connect the instrument and the computer directly or to a LAN” on page 79

“Download and install the QuantStudio™ Design and Analysis desktop Software” on page 82

## Before you begin

This section provides instructions for customer installation of the computer and the instrument. For installation by a field service engineer (FSE), contact Support to order a service call.

Before starting the installation::

- Review “Networking” on page 84 to determine the configuration for your instrument and obtain the network information you need
- Review site requirements in *QuantStudio™ 3 and 5 Real-Time PCR Systems Site Preparation Guide* (Pub. no. MAN0010405)
- Review the connections in “Instrument and computer connections” on page 81

## Unpack and install the instrument

1. Prepare the installation site as described in *QuantStudio™ 3 and 5 Real-Time PCR Systems Site Preparation Guide* (Pub. no. MAN0010405).
2. Follow the pre-printed instructions on the instrument box to unpack the instrument, accessories, and reference documentation. Save the packing material for future use or recycle it.

The instrument box contains:

- One instrument
  - Accessories: power cable, Ethernet cable, USB drive, reaction tube retainer
  - Shipping plate
- Note:** Save the shipping plate but do not use it to operate the instrument.
- Reference documentation: Welcome note, unpacking and set up instructions card, system documentation insert

If you ordered the wireless adapter, it is provided separately.

3. Place the instrument on the bench.
4. Plug the power cable into the instrument power port at the rear of the instrument and the other end into a receptacle.
5. Connect your instrument as required by your network configuration (see “Supported network configuration options” on page 84 and “Instrument and computer connections” on page 81).
  - Connect an Ethernet cable to the connector on the rear of the instrument and to a computer or a LAN.
  - or
  - Connect wirelessly via the wireless adapter.

**Note:** Do not connect the High Power USB WiFi Module (Cat. no. A26774) to the instrument if it is connected to a network by an Ethernet cable. Configuring the instrument for both wired and wireless connection can interfere with instrument operation.

## Power on and follow the startup wizard

1. Power on the instrument.
2. Follow the startup wizard through the following tasks:
  - Accept license agreement
  - Choose networking option  
Select **Wired** or **Wireless** as needed for your configuration (see “Supported network configuration options” on page 84)
  - Configure instrument  
Specify time zone, date format, and time format.
  - Create administrator profile

You can perform any of the steps above at a later time if you do not have the information needed to complete the startup screens. See “Configure instrument settings” on page 25.

---

**IMPORTANT!** Before using the instrument for the first time, we recommend that you “Perform instrument verification using RNase P plates” on page 49.

---

## Connect the instrument and the computer directly or to a LAN

This section describes direct wired connection of the computer provided by Thermo Fisher Scientific to the instrument or to a LAN.

Do not connect a customer-provided computer to the instrument.

1. Connect an Ethernet cable from the instrument or a LAN to the computer.
2. Power on the computer, then log in using a Windows™ Administrator account.
3. In the Windows™ desktop, right-click **My Network Places ▶ Properties**.
4. Right-click **Local Area Connection**, then select **Properties**.
5. Select **Internet Protocol (TCP/IP) ▶ Properties**.


6. Set the Internet Protocol (TCP/IP) Properties for either DHCP or Static IP communication:

Network configuration	Action
DHCP	<ol style="list-style-type: none"> <li>1. Select <b>Obtain an IP address automatically</b>.</li> <li>2. Set the DNS address. If the computer obtains DNS addresses: <ul style="list-style-type: none"> <li>• Automatically – Select <b>Obtain DNS server address automatically</b>.</li> <li>• Statically – Select <b>Use the following DNS address</b>, then enter the address of the preferred and alternate DNS servers (if available).</li> </ul> </li> </ol>
Static IP	<ol style="list-style-type: none"> <li>1. Select <b>Use the following IP address</b>.</li> <li>2. In the IP Address field, enter the static IP address.</li> <li>3. If necessary, enter a subnet mask.</li> <li>4. If necessary, enter a static gateway address in the Default Gateway field.</li> </ol>


7. If your network requires advanced TCP/IP setup (such as WINS), define the settings:
- a. Click **Advanced** in the Internet Protocol (TCP/IP) Properties dialog box.
  - b. Define the IP Settings, DNS, and WINS tabs as instructed by your systems administrator, then click **OK**.
8. Close all dialog boxes by clicking **OK**, then re-start the computer.  
The computer is now visible to other computers on the network.

## Instrument and computer connections

**Figure 3** Instrument back panel

- |   |                                     |
|---|-------------------------------------|
| ① USB ports   | ④ RS232 Port – For service use only |
| ② WiFi USB port – Connect USB wireless adapter for wireless network access (ordered separately) | ⑤ Fuse Cover                        |
| ③ Ethernet Port – RJ45 port for 10/100 Mbps Ethernet communication with the instrument          | ⑥ Power Switch                      |
|   | ⑦ Power Port – 100 to 240 VAC       |

**Figure 4** Instrument-to-computer connections

- |   |   |
|---|---|
| ① Detachable power supply cord compatible with local power supply receptacle. | ③ Connection between the computer and the monitor, keyboard, and mouse. |
| ② Connection between the computer and the instrument (direct connection).     | ④ Connection between the computer and the handheld barcode scanner.     |

## Download and install the QuantStudio™ Design and Analysis desktop Software

### Computer requirements for the desktop software

If you purchased a computer provided by Thermo Fisher Scientific, you can install the QuantStudio™ Design and Analysis desktop Software and use it to control the instrument.


We do not support the use of customer-provided computers to control the instrument.

However, you can install the desktop software on a customer-provided computer and use the software to create templates and analyze data. Minimum requirements for a customer-provided computer are:

- Windows™ 7 operating system (32-bit or 64-bit)
- Pentium® 4 processor or compatible
- 4 GB RAM
- 500 GB hard drive
- Monitor resolution 1280x1024


### Download the desktop software

1. Sign in to your account at **thermofisher.com**.

**Note:** If you do not have an account, create one.


2. Go to <http://www.thermofisher.com/quantstudio3-5softwaredownloads>.
3. Download the software and the example files.

### Install the software

1. Log in to the computer on which you are installing the software with a Windows™ Administrator account.
2. Unzip the downloaded software and example files.

3. Double-click **setup.exe**.

4. Follow the InstallShield Wizard prompts to install the software.
5. Accept the License Agreement.
6. Select **Typical** as the setup preference, then click **Next**.
7. If you are installing the software on a computer provided by Thermo Fisher Scientific, install the software on the D:\ drive. If you are installing on a customer-provided computer, install the software in your preferred location.

8. Click **Finish**.

## Networking

**IMPORTANT!** This section provides general networking information. It *does not* provide adequate detail to integrate the instrument into all possible network architectures. Because a network may contain advanced features (such as a firewall or network domains), we recommend that you consult a network administrator before connecting the instrument to your laboratory network.

### Supported network configuration options

We support the following direct, LAN (local area network) and Cloud network configurations. Configurations other than those listed are not recommended.

DO NOT connect the High Power USB WiFi Module (Cat. no. A26774) to the instrument if it is connected to a network by an Ethernet cable. Configuring the instrument for both wired and wireless connection can interfere with instrument operation.

Connection	Requirements
<p>Direct</p> 	<ul style="list-style-type: none"> <li>A computer provided by Thermo Fisher Scientific installed with the QuantStudio™ Design and Analysis desktop Software</li> <li>Direct Ethernet connection instrument-to-computer</li> </ul>
<p>Local area network (LAN) connection</p> 	<ul style="list-style-type: none"> <li>A computer provided by Thermo Fisher Scientific installed with the QuantStudio™ Design and Analysis desktop Software</li> <li>Direct Ethernet connection <i>or</i> wireless connection instrument-to-LAN (Wi-Fi module connected to instrument is required for wireless connection)</li> <li>Direct Ethernet connection <i>or</i> wireless connection computer-to-LAN</li> <li>Open network ports: 5353 and 7000</li> </ul> <p>Multiple computers can access the same instrument through the network, but only one computer can start a run at any given time.</p>
<p>Cloud</p> 	<ul style="list-style-type: none"> <li>A Microsoft™ Windows™ or Apple™ Macintosh™ computer with internet connection and web browsing access (for Cloud access only, no instrument control)</li> <li>Direct Ethernet connection <i>or</i> wireless connection instrument-to-Cloud (Wi-Fi module connected to instrument is required for wireless connection)</li> <li>Direct Ethernet connection <i>or</i> wireless connection computer-to-Cloud</li> <li>Open network ports: 5353 and 7000</li> </ul>

## Control and monitor networked instruments

When the instrument is connected to a network:

- Computers on the network that are running the desktop software can control the instrument. Networked instruments can be controlled by only one computer at a time.
- Instruments linked to the Cloud cannot be controlled remotely. However, you can:
  - Remotely access the cloud software to create (and analyze) experiments.
  - From the instrument, download the experiments and start a run.
  - Monitor a run from the Cloud in real time.

## About the Ethernet port

The Ethernet port of the instrument supports:

- Static IP network service with subnet mask, primary and secondary data network service (DNS), and default gateway settings, or dynamic host configuration protocol (DHCP) network service.
- mDNS/DNS for local domains.

**Note:** Because mDNS is limited to direct network connections, an instrument configured for mDNS may not be visible to other nodes that are separated by a router, hub, or another network device.

- IPv4 linklocal (IPv4LL) in the RFC (also known as Automatic Private IP Addressing [APIPA] or Internet Protocol Automatic Configuration [IPAC])

**Note:** When an instrument is set for DHCP, APIPA is automatically enabled, and the instrument provides an IP address when no address is supplied by the DHCP server.

## Firewall ports that must be open

Ports	Condition
80/443	Standard ports for instrument-to-Cloud and computer-to-Cloud connection
mDNS, 7000	Instrument-to-computer connection
mDNS, 5353	Instruments discovery

## Networking guidelines and best practices

- Consult a network administrator before connecting the instrument to a network.
- To enable the full functionality of the software, the computer requires a network connection.
- Open the firewall port for the instruments to be discovered.  
See “Firewall ports that must be open” on page 85 for information on relevant ports to open on the computer and router.
- Observe the restrictions to mDNS and Autodiscovery.  
The instrument supports mDNS but only when the instrument and computer share a direct network connection and are within the same subnet. Network computers that are separated from the instrument by a router, hub, or another network device may not be able to access the instrument by its host name.

- Confirm the uniqueness of the instrument name.
  - The instrument name must be unique within the subnet. The desktop software can automatically discover instruments on the link-local network.
  - The instrument does not test the uniqueness of the instrument name within the subnet when it is set.

# Troubleshooting

Observation	Possible cause	Recommended action
Inconsistent communication between instrument-computer or instrument-Cloud	Instrument is configured for both wired and wireless network connection.	Ensure only one connectivity option is plugged into the instrument (either an Ethernet cable <i>or</i> wireless adapter, but not both).
		Configure for wired <i>or</i> wireless network connection.
	Weak or unstable internet connection, especially if configured for wireless.	Change configuration to wired connection.
		Use a wireless network with a stronger or more consistent signal.
The connection between instrument and computer not recognized	The connection is not fully established.	Power the instrument off, then power on again.
Insufficient disk space message	Insufficient disk space to save a run.	<ol style="list-style-type: none"> <li>1. In the home screen, touch  <b>Settings ▶ Run History ▶ Manage</b>.</li> <li>2. Delete or transfer experiments from the instrument.</li> </ol>
Touchscreen is black	Instrument is in sleep mode.	Touch anywhere on the instrument touchscreen.
	Instrument is not powered on.	If you touch the instrument touchscreen and it remains black, check if the instrument is powered on. The power switch is located on the rear panel of the instrument.
		If the instrument does not power on, check that the power supply is properly connected.
		If the instrument does not power on and the power supply is properly connected, contact Support.
Forgot PIN for instrument profile	Non-administrator forgot instrument profile PIN.	See “Manage all instrument profiles (Administrator only)” on page 71.
	Administrator forgot instrument profile PIN.	Have another administrator reset the PIN for the forgotten-PIN profile (see “Manage all instrument profiles (Administrator only)” on page 71).
		If there is not another administrator profile on the instrument, you must restore factory defaults (see “Restore factory defaults” on page 74).

# Parts and materials

■ Kits, consumables, and accessories .....	88
■ 96-Well (0.2-mL) consumables .....	88
■ 384-well consumables .....	89
■ 96-Well Fast (0.1-mL) consumables .....	90
■ QuantStudio™ 3 and 5 Systems accessories .....	91
■ General-use materials and consumables .....	91

## Kits, consumables, and accessories

The following kits, consumables, and reagents are used with the QuantStudio™ 3 and 5 Systems.

Unless otherwise indicated, all materials are available through **thermofisher.com**.

**Note:** Store ROI/Uniformity, Background, Dye, and RNase P plates at –20°C and use them by the expiration date on the packaging. All other consumables can be stored at room temperature.

## 96-Well (0.2-mL) consumables

Consumable	Contents	Cat. no.
MicroAmp™ Optical 8-Cap Strip	300 strips	4323032
MicroAmp™ Optical 8-Tube Strip (0.2 mL)	125 strips	4316567
MicroAmp™ Optical Tube without Cap (0.2 mL)	2000 tubes	N8010933
MicroAmp™ 96-Well Tray/Retainer Set (blue) (for 0.2 mL)	10 sets	4381850
MicroAmp™ Optical 96-Well Reaction Plate (0.2 mL)	10 plates	N8010560
	500 plates	4316813
MicroAmp™ EnduraPlate™ Optical 96-Well Reaction Plate with Barcode (blue) (0.2 mL)	20 plates	4483343
MicroAmp™ Optical Adhesive Film Kit	1 kit	4313663

Calibration or instrument verification plate	Cat. no.
96-Well Region of Interest (ROI) and Background Plates (2 plates)	4432364
QuantStudio™ 3 or 5 10-Dye Spectral Calibration Kit, 96-Well 0.2-mL	A26343
96-Well 0.2-mL Spectral Calibration Plate 1 (containing FAM™, VIC™, ROX™, and SYBR™ dyes)	A26331
96-Well 0.2-mL Spectral Calibration Plate 2 (containing ABY™, JUN™, and MUSTANG PURPLE™ dyes)	A26332
96-Well 0.2-mL Spectral Calibration Plate 3 (containing TAMRA™, NED™, and Cy <sup>fi</sup> 5 dyes)	A26333
96-Well 0.2-mL TaqMan™ RNase P Instrument Verification Plate	4432382

## 384-well consumables

Consumable	Contents	Cat. no.
MicroAmp™ Optical 384-Well Reaction Plate with Barcode	50 plates	4309849
	500 plates	4326270
	1000 plates	4343814
MicroAmp™ EnduraPlate™ Optical 384-Well Reaction Plate with Barcode (clear)	20 plates	4483285
	500 plates	4483273
MicroAmp™ Optical Adhesive Film Kit	1 kit	4313663

Calibration or instrument verification plate	Cat. no.
384-Well Region of Interest (ROI) and Background Plates (2 plates)	4432320
QuantStudio™ 3 or 5 10-Dye Spectral Calibration Kit, 384-Well	A26341
384-Well Spectral Calibration Plate 1 (containing FAM™, VIC™, ROX™, TAMRA™, and SYBR™ dyes)	A26334
384-Well Spectral Calibration Plate 2 (containing ABY™, JUN™, MUSTANG PURPLE™, NED™, and Cy <sup>fi</sup> 5 dyes)	A26335
384-Well TaqMan™ RNase P Instrument Verification Plate	4455280

## 96-Well Fast (0.1-mL) consumables

Consumable	Contents	Cat. no.
MicroAmp™ Optical 8-Cap Strip	300 strips	4323032
MicroAmp™ Optical Fast 8-Tube Strip (0.1 mL)	125 strips	4358293
MicroAmp™ Optical Fast Tube with Cap (0.1 mL)	1000 tubes	4358297
MicroAmp™ 96-Well Tray (blue) (for 0.1 mL)	10 trays	4379983
MicroAmp™ Optical 96-Well Fast Reaction Plate (0.1 mL)	10 plates	4346907
MicroAmp™ EnduraPlate™ Optical 96-Well Fast Reaction Plate with Barcode (clear) (0.1 mL)	20 plates	4481194
	500 plates	4483494
MicroAmp™ Optical Adhesive Film Kit	1 kit	4313663

Calibration or instrument verification plate	Cat. no.
96-Well Fast Region of Interest (ROI) and Background Plates (2 plates)	4432426
QuantStudio™ 3 or 5 10-Dye Spectral Calibration Kit, 96-Well 0.1-mL	A26342
96-Well 0.1-mL Spectral Calibration Plate 1 (containing FAM™, VIC™, ROX™, and SYBR™ dyes)	A26336
96-Well 0.1-mL Spectral Calibration Plate 2 (containing ABY™, JUN™, and MUSTANG PURPLE™ dyes)	A26337
96-Well 0.1-mL Spectral Calibration Plate 3 (containing TAMRA™, NED™, and Cy <sup>fi</sup> 5 dyes)	A26340
96Well Fast TaqMan™ RNase P Instrument Verification Plate	4351979

## QuantStudio™ 3 and 5 Systems accessories

Accessory	Contents	Cat. no.
MicroAmp™ Multi-Removal Tool	1 tool	4313950
MicroAmp™ Cap Installing Tool (handle style)	1 tool	4330015
MicroAmp™ Optical Adhesive Film Kit	1 kit	4313663
MicroAmp™ Optical Adhesive Film	25 films	4360954
	100 films	4311971
MicroAmp™ Adhesive Film Applicator	5 applicators	4333183
Real Time PCR Grade Water	10 × 1.5 mL tubes	AM9935
Handheld Barcode Scanner	1 scanner	448842
High Power USB WiFi Module	1 module	A26774

## General-use materials and consumables

The following general-use materials and consumables are required to calibrate, maintain, and operate the instrument. Unless indicated otherwise, all materials shown below are available from major laboratory suppliers (MLS).

Material/consumable	Source
Bleach, 10% solution	MLS
Centrifuge with 96-well plate buckets	MLS
Cotton or nylon swabs and lint-free cloths	MLS
Ethanol, 95% solution	MLS
Optical clear adhesive film for PCR	MLS
Pipettors, 100-µL and 200-µL (with pipette tips)	MLS
Powder-free gloves	MLS
Safety glasses	MLS
Screwdriver, flathead	MLS
Tissue, lint-free	MLS
Deionized water	MLS



# Instrument specification and layout

■ Configured system dimensions .....	93
■ Electrical requirements .....	94
■ Environmental requirements .....	95
■ Network requirements .....	95

**Note:** For information on instrument components and connections, see “Instrument and computer connections” on page 81.

## Configured system dimensions

Allow space for the configured instrument. A typical setup is shown below (dimensions are rounded to nearest unit).

## Instrument clearances

During instrument setup and maintenance, it is necessary to access the back of the instrument. If the back of the instrument faces a wall, it will be necessary to have enough space to rotate the instrument on the bench for access.

**IMPORTANT!** For safety, the power outlet used for powering the instrument must be accessible at all times.

Component	Top	Front	Sides	Back
QuantStudio™ 3 and 5 Real-Time PCR Systems	30.5 cm (12 in.)	30.5 cm (12 in.)	15.25 cm (6 in.)	15.25 cm (6 in.)
Computer <i>(Optional)</i>	—	15.25 cm (6 in.)	—	15.25 cm (6 in.)

## Electrical requirements

**WARNING!** For safety, the power outlet used for powering the instrument must be accessible at all times. See “Instrument clearances” for information about the space needed between the wall and the instrument. In case of emergency, you must be able to immediately disconnect the main power supply to all the equipment. Allow adequate space between the wall and the equipment so that the power cords can be disconnected in case of emergency.

- Electric receptacle with grounding capability
- Maximum power dissipation: 960 W (approximate, not including computer and monitor)
- Mains AC line voltage tolerances must be up to  $\pm 10$  percent of nominal voltage

Device	Rated voltage	Circuit required	Rated frequency	Rated power
QuantStudio™ 3 and 5 Systems	100–240 $\pm 10\%$ VAC <sup>[1]</sup>	10 A	50/60 Hz	960 W
Computer (desktop)	100–240 $\pm 10\%$ VAC	10 A	50/60 Hz	125 VA
Monitor				65 VA
Computer (laptop)	100–240 $\pm 10\%$ VAC	10 A	50/60 Hz	90 VA

<sup>[1]</sup> If the supplied power fluctuates beyond the rated voltage, a power line regulator may be required. High or low voltages can adversely affect the electronic components of the instrument.

## Environmental requirements

Table 2 Environmental requirements

Condition	Acceptable range
Installation site	Indoor use only
Electromagnetic interference	Do not use this device in close proximity to sources of strong electromagnetic radiation (for example, unshielded intentional RF sources). Strong electromagnetic radiation may interfere with the proper operation of the device.
Altitude	Between sea level and 2000 m (6500 ft.) above sea level
Operating conditions	<ul style="list-style-type: none"> <li>Temperature: 15 to 30°C (60 to 85°F)</li> </ul> <p><b>Note:</b> The room temperature must not fluctuate more than <math>\pm 2^\circ\text{C}</math> over a 2-hour period.</p> <ul style="list-style-type: none"> <li>Humidity: 15–80% relative humidity (noncondensing)</li> </ul>
Thermal output	During operation, the net thermal output, based on the actual current draw of the instrument, is expected to be approximately 960 W (3275 Btu/h).
Vibration	Ensure that the instrument is not adjacent to strong vibration sources, such as a centrifuge, pump, or compressor. Excessive vibration will affect instrument performance.
Pollution degree	<p>The instrument has a Pollution Degree rating of II. The instrument may only be installed in an environment that has nonconductive pollutants such as dust particles or wood chips. Typical environments with a Pollution Degree II rating are laboratories and sales and commercial areas.</p> <p>The noise output of the instrument is <math>\leq 60</math> dB when running.</p>
Other conditions	<p>Ensure the instrument is located away from any vents that could expel particulate material onto the instrument components.</p> <p>Avoid placing the instrument and computer adjacent to heaters, cooling ducts, or in direct sunlight.</p>

## Network requirements

The instrument:

- Is factory-configured for IPv4 TCP/IP communication and includes a fast Ethernet adapter (10/100 Mbps) with an RJ45-type connector for integrating the device into a local area network (LAN).
- Can alternatively be configured for wireless networking (wireless dongle required, sold separately as an optional accessory).

The instrument should be configured for *either* wired or wireless networking, not both.

If a field service representative is to install your instrument:

- If the instrument will be connected to a LAN, an active, tested network jack must be in place before the scheduled installation date.
- A representative from your information technologies department must be available during the installation to help connect the instrument to your network.

**Required materials:**

- Wired: CAT6 Ethernet cable of sufficient length with RJ45 connectors (for a 1000 Mbps network connection or a CAT5 for a 100 Mbps connection)  
*or*
- Wireless: 802.11b/g/n Single-Band Wireless Dongle (Cat. no. A26774, ordered separately)

# Safety

■ Symbols on this instrument .....	98
■ Safety alerts on this instrument .....	100
■ Safety information for instruments not manufactured by Thermo Fisher Scientific .....	101
■ Instrument safety .....	101
■ Safety and electromagnetic compatibility (EMC) standards .....	103
■ Chemical safety .....	105
■ Biological hazard safety .....	106

---

 **WARNING! GENERAL SAFETY.** Using this product in a manner not specified in the user documentation may result in personal injury or damage to the instrument or device. Ensure that anyone using this product has received instructions in general safety practices for laboratories and the safety information provided in this document.

- Before using an instrument or device, read and understand the safety information provided in the user documentation provided by the manufacturer of the instrument or device.
  - Before handling chemicals, read and understand all applicable Safety Data Sheets (SDSs) and use appropriate personal protective equipment (gloves, gowns, eye protection, etc). To obtain SDSs, see the “Documentation and Support” section in this document.
-


## Symbols on this instrument

Symbols may be found on the instrument to warn against potential hazards or convey important safety information. In this document, the hazard symbol is used along with one of the following user attention words:

- **CAUTION!** – Indicates a potentially hazardous situation that, if not avoided, may result in minor or moderate injury. It may also be used to alert against unsafe practices.
- **WARNING!** – Indicates a potentially hazardous situation that, if not avoided, could result in death or serious injury.
- **DANGER!** – Indicates an imminently hazardous situation that, if not avoided, will result in death or serious injury.

Symbol	English	Français
	Caution, risk of danger Consult the manual for further safety information.	Attention, risque de danger Consulter le manuel pour d'autres renseignements de sécurité.
	Caution, risk of electrical shock	Attention, risque de choc électrique
	Moving parts	Parties mobiles
	Caution, hot surface	Attention, surface chaude
	Potential biohazard	Danger biologique potentiel
	Ultraviolet light	Rayonnement ultraviolet
	Potential slipping hazard	Danger de glisser potentiel
I	On	On (marche)
O	Off	Off (arrêt)
ⓘ	On/Off	On/Off (marche/arrêt)
⏻	Standby	En attente

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## Conformity symbols


Conformity mark	Description
	Indicates conformity with safety requirements for Canada and U.S.A.
	Indicates conformity with European Union requirements for safety and electromagnetic compatibility.
	Indicates conformity with Australian standards for electromagnetic compatibility.

## Safety alerts on this instrument

Additional text may be used with one of the symbols described above when more specific information is needed to avoid exposure to a hazard. See the following table for safety alerts found on the instrument.

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### Location of safety labels on the instrument

## Safety information for instruments not manufactured by Thermo Fisher Scientific

Some of the accessories provided as part of the instrument system are not designed or built by Thermo Fisher Scientific. Consult the manufacturer's documentation for the information needed for the safe use of these products.

### Instrument safety

#### General

**CAUTION! Do not remove instrument protective covers.** If you remove the protective instrument panels or disable interlock devices, you may be exposed to serious hazards including, but not limited to, severe electrical shock, laser exposure, crushing, or chemical exposure.

**CAUTION! Solvents and Pressurized fluids.** Wear eye protection when working with any pressurized fluids. Use caution when working with any polymeric tubing that is under pressure:

- Extinguish any nearby flames if you use flammable solvents.
- Do not use polymeric tubing that has been severely stressed or kinked.
- Do not use polymeric tubing with tetrahydrofuran or nitric and sulfuric acids.
- Be aware that methylene chloride and dimethyl sulfoxide cause polymeric tubing to swell and greatly reduce the rupture pressure of the tubing.
- Be aware that high solvent flow rates (~40mL/min) may cause a static charge to build up on the surface of the tubing and electrical sparks may result.

## Physical injury

---

**CAUTION! Moving and Lifting Injury.** The instrument is to be moved and positioned only by the personnel or vendor specified in the applicable site preparation guide. Improper lifting can cause painful and permanent back injury.

Things to consider before lifting or moving the instrument or accessories:

- Depending on the weight, moving or lifting may require two or more persons.
- If you decide to lift or move the instrument after it has been installed, do not attempt to do so without the assistance of others, the use of appropriate moving equipment, and proper lifting techniques.
- Ensure you have a secure, comfortable grip on the instrument or accessory.
- Make sure that the path from where the object is to where it is being moved is clear of obstructions.
- Do not lift an object and twist your torso at the same time. Keep your spine in a good neutral position while lifting with your legs.
- Participants should coordinate lift and move intentions with each other before lifting and carrying.
- For smaller packages, rather than lifting the object from the packing box, carefully tilt the box on its side and hold it stationary while someone else slides the contents out of the box.

---

**CAUTION! Moving Parts.** Moving parts can crush, pinch and cut. Keep hands clear of moving parts while operating the instrument. Disconnect power before servicing.

---

---

**WARNING!** Do not attempt to lift or move the instrument without the assistance of others. Use appropriate moving equipment and proper lifting technique, improper lifting may result in serious injury.

---

## Electrical

**WARNING! Fuse Installation.** Before installing the instrument, verify that the fuses are properly installed and the fuse voltage matches the supply voltage. Replace fuses only with the type and rating specified for the unit. Improper fuses can damage the instrument wiring system and cause a fire.

**WARNING! Ensure appropriate electrical supply.** For safe operation of the instrument:

- Plug the system into a properly grounded receptacle with adequate current capacity.
- Ensure the electrical supply is of suitable voltage.
- Never operate the instrument with the ground disconnected. Grounding continuity is required for safe operation of the instrument.

**WARNING! Power Supply Line Cords.** Use properly configured and approved line cords for the power supply in your facility.

**WARNING! Disconnecting Power.** To fully disconnect power either detach or unplug the power cord, positioning the instrument such that the power cord is accessible.

## Cleaning and decontamination

**CAUTION! Cleaning and Decontamination.** Use only the cleaning and decontamination methods specified in the manufacturer's user documentation. It is the responsibility of the operator (or other responsible person) to ensure the following requirements are met:

- No decontamination or cleaning agents are used that could cause a HAZARD as a result of a reaction with parts of the equipment or with material contained in the equipment.
- The instrument is properly decontaminated a) if hazardous material is spilled onto or into the equipment, and/or b) prior to having the instrument serviced at your facility or sending the instrument for repair, maintenance, trade-in, disposal, or termination of a loan (decontamination forms may be requested from customer service).
- Before using any cleaning or decontamination methods (except those recommended by the manufacturer), users should confirm with the manufacturer that the proposed method will not damage the equipment.

## Safety and electromagnetic compatibility (EMC) standards

The instrument design and manufacture complies with the standards and requirements for safety and electromagnetic compatibility as noted in the following table:

**Safety compliance**

Reference	Description
EU Directive 2014/35/EU	European Union “Low Voltage Directive”
IEC 61010-1 EN 61010-1 UL 61010-1 CSA C22.2 No. 61010-1	<i>Safety requirements for electrical equipment for measurement, control, and laboratory use – Part 1: General requirements</i>
IEC 61010-2-010 EN 61010-2-010	<i>Safety requirements for electrical equipment for measurement, control and laboratory use – Part 2-010: Particular requirements for laboratory equipment for the heating of materials</i>
IEC 61010-2-081 EN 61010-2-081	<i>Safety requirements for electrical equipment for measurement, control and laboratory use – Part 2-081: Particular requirements for automatic and semi-automatic laboratory equipment for analysis and other purposes</i>

**EMC**

Reference	Description
Directive 2014/30/EU	European Union “EMC Directive”
EN 61326-1/ IEC 61326-1	<i>Electrical Equipment for Measurement, Control and Laboratory Use – EMC Requirements – Part 1: General Requirements</i>
AS/NZS CISPR 11	<i>Limits and Methods of Measurement of Electromagnetic Disturbance Characteristics of Industrial, Scientific, and Medical (ISM) Radiofrequency Equipment</i>
ICES-001, Issue 4	<i>Industrial, Scientific and Medical (ISM) Radio Frequency Generators</i>
FCC Part 15 Subpart B (47 CFR)	<i>U.S. Standard Radio Frequency Devices</i>

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Environmental  
design

Reference	Description
Directive 2012/19/EU	European Union “WEEE Directive” – Waste electrical and electronic equipment
Directive 2011/65/EU	European Union “RoHS Directive” – Restriction of hazardous substances in electrical and electronic equipment
Directive 2006/66/EC	European Union “Battery Directive”

## Chemical safety

**WARNING! GENERAL CHEMICAL HANDLING.** To minimize hazards, ensure laboratory personnel read and practice the general safety guidelines for chemical usage, storage, and waste provided below, and consult the relevant SDS for specific precautions and instructions:

- Read and understand the Safety Data Sheets (SDSs) provided by the chemical manufacturer before you store, handle, or work with any chemicals or hazardous materials. To obtain SDSs, see the “Documentation and Support” section in this document.
- Minimize contact with chemicals. Wear appropriate personal protective equipment when handling chemicals (for example, safety glasses, gloves, or protective clothing).
- Minimize the inhalation of chemicals. Do not leave chemical containers open. Use only with adequate ventilation (for example, fume hood).
- Check regularly for chemical leaks or spills. If a leak or spill occurs, follow the manufacturer's cleanup procedures as recommended in the SDS.
- Handle chemical wastes in a fume hood.
- Ensure use of primary and secondary waste containers. (A primary waste container holds the immediate waste. A secondary container contains spills or leaks from the primary container. Both containers must be compatible with the waste material and meet federal, state, and local requirements for container storage.)
- After emptying a waste container, seal it with the cap provided.
- Characterize (by analysis if necessary) the waste generated by the particular applications, reagents, and substrates used in your laboratory.
- Ensure that the waste is stored, transferred, transported, and disposed of according to all local, state/provincial, and/or national regulations.
- **IMPORTANT!** Radioactive or biohazardous materials may require special handling, and disposal limitations may apply.

**WARNING! HAZARDOUS WASTE (from instruments).** Waste produced by the instrument is potentially hazardous. Follow the guidelines noted in the preceding General Chemical Handling warning.

**WARNING! 4L Reagent and Waste Bottle Safety.** Four-liter reagent and waste bottles can crack and leak. Each 4-liter bottle should be secured in a low-density polyethylene safety container with the cover fastened and the handles locked in the upright position.

## Biological hazard safety

---

**WARNING! Potential Biohazard.** Depending on the samples used on this instrument, the surface may be considered a biohazard. Use appropriate decontamination methods when working with biohazards.

---

---

**WARNING! BIOHAZARD.** Biological samples such as tissues, body fluids, infectious agents, and blood of humans and other animals have the potential to transmit infectious diseases. All work should be conducted in properly equipped facilities using the appropriate safety equipment (for example, physical containment devices). Safety equipment also may include items for personal protection, such as gloves, coats, gowns, shoe covers, boots, respirators, face shields, safety glasses, or goggles. Individuals should be trained according to applicable regulatory and company/ institution requirements before working with potentially biohazardous materials. Follow all applicable local, state/provincial, and/or national regulations. The following references provide general guidelines when handling biological samples in laboratory environment.

- U.S. Department of Health and Human Services, *Biosafety in Microbiological and Biomedical Laboratories (BMBL)*, 5th Edition, HHS Publication No. (CDC) 21-1112, Revised December 2009; found at:  
[www.cdc.gov/biosafety/publications/bmbl5/BMBL.pdf](http://www.cdc.gov/biosafety/publications/bmbl5/BMBL.pdf)
  - World Health Organization, *Laboratory Biosafety Manual*, 3rd Edition, WHO/CDS/CSR/LYO/2004.11; found at:  
[www.who.int/csr/resources/publications/biosafety/Biosafety7.pdf](http://www.who.int/csr/resources/publications/biosafety/Biosafety7.pdf)
-

# Documentation and support

## Related documentation

Document	Publication number	Description
<i>QuantStudio™ 3 and 5 Real-Time PCR Systems Installation, Use, and Maintenance Guide</i>	MAN0010407	Describes the QuantStudio™ 3 and 5 Real-Time PCR Systems hardware and software and provides information on preparing, using, maintaining, and troubleshooting the system.
<i>QuantStudio™ Real-Time PCR System Help</i>	MAN0010422	Describes the QuantStudio™ 3 and 5 Real-Time PCR Systems touchscreen and provides procedures for configuration, calibration, and performing a run.
<i>QuantStudio™ Design and Analysis desktop Software Command-Line Application Guide</i>	MAN0010409	Describes how to use the command-line interface of the QuantStudio™ Design and Analysis desktop Software and provides the procedure to automate the creation of new experiment files and export data from existing files.
<i>QuantStudio™ Design and Analysis desktop Software User Guide</i>	MAN0010408	Describes how to perform the six different experiments on the QuantStudio™ Design and Analysis desktop Software
<i>QuantStudio™ Design and Analysis desktop Software Help</i>	MAN0010415	Describes the QuantStudio™ Design and Analysis desktop Software and provides procedures for common tasks.
<i>SAE Admin Console Help</i>	MAN0010417	Describes the Security, Audit, and e-Signature (SAE) Administrator Console and provides procedures for common tasks.
<i>SAE Admin Console User Guide</i>	MAN0010410	Describes how to use the Security, Audit, and e-Signature (SAE) Administrator Console.

Document	Publication number	Description
<i>QuantStudio™ Design and Analysis cloud Software Help</i>	MAN0010414	Describes the QuantStudio™ Design and Analysis cloud Software and provides procedures for common tasks.
<i>QuantStudio™ 3 and 5 Real-Time PCR Systems Site Preparation Guide</i>	MAN0010405	Explains how to prepare your site to receive and install the QuantStudio™ 3 and 5 Real-Time PCR Systems. Intended for personnel who schedule, manage, and perform the tasks required to prepare the site for installation of the QuantStudio™ 3 and 5 Real-Time PCR Systems.

**Note:** For additional documentation, see “Customer and technical support” on page 108.

## Obtain information from the Help system

The instrument has a Help system that describes how to use each feature of the touchscreen. Touch  on the instrument touchscreen to access the Help system.

## Customer and technical support

Visit [thermofisher.com/support](https://thermofisher.com/support) for the latest in services and support, including:

- Worldwide contact telephone numbers
- Product support, including:
  - Product FAQs
  - Software, patches, and updates
- Order and web support
- Product documentation, including:
  - User guides, manuals, and protocols
  - Certificates of Analysis
  - Safety Data Sheets (SDSs; also known as MSDSs)

**Note:** For SDSs for reagents and chemicals from other manufacturers, contact the manufacturer.

## Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at **[www.thermofisher.com/us/en/home/global/terms-and-conditions.html](http://www.thermofisher.com/us/en/home/global/terms-and-conditions.html)**. If you have any questions, please contact Life Technologies at **[www.thermofisher.com/support](http://www.thermofisher.com/support)**.

# Index

## A

- accessories 91
- add to method 34
- administrator profile 71

## B

- Background calibration
  - create plate 48
  - description 41
  - plate, prepare 43
  - run 44
  - when to perform 40
- backup instrument 62
- barcode scan 31
- biohazard safety 106

## C

- calibration
  - reminders 41, 42
  - Background 44
  - custom dye. *See* custom dye
  - custom melt 60
  - Dye 44
  - plate, prepare 43
  - reminders 41, 42
  - ROI/Uniformity 44
  - schedule 40
  - status 41, 42
  - transfer results 46
  - view images 46
  - workflow 42
- calibration failure, troubleshoot 47
- Cloud server, select location 75
- Comparative CT experiment 17
- computer
  - connection to instrument or LAN 84
  - connections 81
  - requirements for desktop software 82
  - set up 79
- connect to network, workflow 77
- consumables
  - 384-well 89

- 96-well 0.1-mL 90
- 96-well 0.2-mL 88
- contamination, identify 47
- custom dye
  - add to the software 57
  - calibration plate, create 57
  - calibration, perform 58
  - dilution guidelines 54
  - optimal concentration 56
  - wavelength requirements 54
- custom dye calibration, workflow 54
- custom dye dilution plate, prepare 54
- custom dyes 12
- custom melt calibration 60

## D

- data collection 12
- date and time 76
- desktop software, download and install 82
- disk space message
- documentation, related 107
- Dye calibration
  - description 41
  - plate, prepare 43
  - run 44
  - when to perform 40
- dyes
  - custom 12, 54
  - system 11

## E

- experiment
  - copy or delete 39
  - create from template 29
  - edit 31
  - run from saved file 29
  - run last 30
  - transfer 38
- experiment properties. *See* properties
- experiment types 17

**F**

factory defaults 74  
 files, copy or delete 39  
 firewall ports 85  
 fuses, replace 66

**G**

Genotyping experiment 17  
 guest profile  
   disable 73  
   limitations 23

**H**

Help system, access 108  
 home screen, parts of 21

**I**

installation  
   before you begin 78  
   desktop software 82  
   workflow 77  
 installation specification 50  
 instrument  
   profile, create 23  
   connections 81  
   enable for monitoring in cloud software 73  
   move 68  
   name 75  
   overview 10  
   parts of 13  
   power off 67  
   power on 20, 67  
   prepare to move or ship or store 67  
   return for service 68  
 instrument filters 11  
 instrument layout 92  
 instrument settings 25  
 instrument specification 92  
 instrument verification. *See* RNase P verification

**L**

limited product warranty 109  
 link to a Cloud account 24  
 load plate 36  
 lock the touchscreen 38

**M**

maintenance 62

Melt Curve experiment 17  
 method

  add steps 33  
   edit 32  
   parts of 33  
   pause 34  
   ramp rate 34

**N**

network  
   configurations supported 84  
   guidelines and best practices 85  
   wired connection set up 74  
   wireless connection set up 75  
 networking 84

**O**

online Help. *See* Help system

**P**

pause run 37  
 performance verification. *See* RNase P verification  
 PIN reset 71  
 plate wells, define 35  
 Presence/Absence experiment 17  
 profiles  
   administrator 71  
   delete 71  
   edit 72  
   guest, disable 73  
   instrument 23  
   link to Cloud 24  
   require sign-in 73  
   reset PIN 71  
   sign in 23  
   view all 71  
 properties, define 31  
 protocol. *See* method

**Q**

QuantStudio Design and Analysis desktop software,  
   download and install 82

**R**

related documentation 107  
 Relative Standard Curve experiment 17  
 Replace with index term  
 require sign in 73  
 restore instrument 62

RNase P verification  
     analytical performance 50  
     plate description 49  
     plate, prepare 50  
     run 51  
 ROI/Uniformity calibration  
     description 41  
     plate, prepare 43  
     run 44  
     when to perform 40  
 run  
     pause or stop 37  
     post-PCR 12  
     real-time 12  
 run history 38  
 run last experiment 30  
 run protocol. *See* method

## S

SAE mode (Security, Audit, and e-Signature) 72  
 safety, biohazard 106  
 sample block  
     clean 63  
     decontaminate 63  
 settings 25  
 sign in, require 73  
 sign out  
     profile 24  
     timer, set 76  
 sleep mode 66  
 software features 14  
 software updates 73  
 Standard Curve experiment 17  
 standby mode 66

startup wizard 79  
 steps, add to method 33  
 stop run 37  
 supported dyes 11  
 symbols, safety 98  
 system dyes 11

## T

terms and conditions 109  
 thermal cycling protocol. *See* method  
 third-party software 16  
 transfer files 38  
 troubleshoot 87

## U

unload plate 36  
 unpack instrument 78  
 use the instrument, workflow 28

## V

verification failure, troubleshoot 53

## W

warranty 109  
 well details 37  
 workflow, use the instrument 28

## Z

zoom 37

For support visit [thermofisher.com/support](http://thermofisher.com/support) or email [techsupport@lifetech.com](mailto:techsupport@lifetech.com)  
[thermofisher.com](http://thermofisher.com)

16 December 2015

**ThermoFisher**  
S C I E N T I F I C

# MagMAX™ Viral/Pathogen II Nucleic Acid Isolation Kit

## INSTRUCTIONS FOR USE

**Catalog Number** A48383

**Publication Number** MAN0019746

**Revision** A.0

For In Vitro Diagnostic Use.

**ThermoFisher**  
S C I E N T I F I C

Thermo Fisher Scientific Baltics UAB | V.A. Graiciuno 8, LT-02241 | Vilnius, Lithuania

For descriptions of symbols on product labels or product documents, go to [thermofisher.com/symbols-definition](https://thermofisher.com/symbols-definition).

The customer is responsible for compliance with regulatory requirements that pertain to their procedures and uses of the instrument.

The information in this guide is subject to change without notice.

**DISCLAIMER:** TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

**Revision history:** Pub. No. MAN0019746

Revision	Date	Description
A.0	20 November 2020	New document.

**TRADEMARKS:** All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

©2020 Thermo Fisher Scientific Inc. All rights reserved.

# Contents

■ CHAPTER 1	Product information .....	5
	Intended use .....	5
	Product information .....	5
	Contents and storage .....	6
	Required materials not supplied .....	6
■ CHAPTER 2	Extract nucleic acid (automated method) .....	8
	Before you begin .....	8
	Extract nucleic acid—Automated method (200-μL sample input volume) .....	9
	Set up the instrument (200-μL sample input volume) .....	9
	Prepare the processing plates (200-μL sample input volume) .....	10
	Prepare Binding Bead Mix (200-μL sample input volume) .....	10
	Prepare sample plate .....	10
	Process the samples (200-μL sample input volume) .....	11
	Extract nucleic acid—Automated method (400-μL sample input volume) .....	12
	Set up the instrument .....	12
	Prepare the processing plates .....	12
	Prepare Binding Bead Mix (400-μL sample input volume) .....	13
	Prepare sample plate .....	13
	Process the samples (400-μL sample input volume) .....	14
■ CHAPTER 3	Extract nucleic acid (manual method) .....	15
	Before you begin .....	15
	Extract nucleic acid—Manual method (200-μL sample input volume) .....	16
	Prepare Binding Bead Mix (200-μL sample input volume) .....	16
	Digest with Proteinase K (200-μL sample input volume) .....	16
	Wash the beads (200-μL sample input volume) .....	17
	Elute the nucleic acid (200-μL sample input volume) .....	17
	Extract nucleic acid—Manual method (400-μL sample input volume) .....	18
	Prepare Binding Bead Mix (400-μL sample input volume) .....	18
	Digest with Proteinase K (400-μL sample input volume) .....	18
	Wash the beads (400-μL sample input volume) .....	19
	Elute the nucleic acid (400-μL sample input volume) .....	20

- **APPENDIX A Safety ..... 21**
  - Chemical safety ..... 21
  - Biological hazard safety ..... 22
- **APPENDIX B Documentation and support ..... 23**
  - Related documentation ..... 23
  - Customer and technical support ..... 23
  - Limited product warranty ..... 23

# Product information

■ Intended use .....	5
■ Product information .....	5
■ Contents and storage .....	6
■ Required materials not supplied .....	6

## Intended use

MagMAX™ Viral/Pathogen II Nucleic Acid Isolation Kit is a magnetic bead technology-based nucleic acid purification kit intended for the isolation and purification of viral RNA and DNA from human nasopharyngeal swabs. The kit is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of magnetic bead purification, either manual or automated, and *in vitro* diagnostic procedures.

## Product information

MagMAX™ Viral/Pathogen II Nucleic Acid Isolation Kit (Cat. No. A48383) is specifically designed to recover RNA and DNA from viral particles contained in viral transport media (VTM). The kit utilizes MagMAX™ magnetic-bead technology, ensuring reproducible recovery of high-quality nucleic acid.

This product is intended for *in-vitro* diagnostic use and includes the following features:

- Automated workflow using the KingFisher™ Flex Magnetic Particle Processor with 96 Deep-Well Head allows for 96 nasopharyngeal swab specimens to be processed in <30 minutes
- Protocol options support an automated workflow or a manual workflow
- Flexible protocol accommodates sample volume inputs from 200 to 400 µL of viral transport medium
- No need for carrier RNA
- Elution volume of 50 µL

## Contents and storage

The MagMAX™ Viral/Pathogen II Nucleic Acid Isolation Kit contains sufficient reagents for 1,000 reactions with 400 µL volume input or 2,000 reactions with 200 µL volume input.

Component	Amount	Storage
Binding Solution	550 mL	15°C to 25°C
Wash Solution	1,000 mL	
Binding Beads	20 mL	
Proteinase K	10 mL	
Elution Buffer	100 mL	

## Required materials not supplied

Unless otherwise indicated, all materials are available through [thermofisher.com](https://www.thermofisher.com). "MLS" indicates that the material is available from [fisherscientific.com](https://www.fisherscientific.com) or another major laboratory supplier.

Catalog numbers that appear as links open the web pages for those products.

Item	Source
<b>Automated nucleic acid extraction system and materials</b>	
KingFisher™ Flex Magnetic Particle Processor with 96 Deep-Well Head	<a href="#">5400630</a>
KingFisher™ Flex 96 Deep-Well Heating Block	24075430
KingFisher™ Deep-Well 96 Plate	<a href="#">95040450</a> , A48305, A48424, 95040455
96-well plate for the tip comb, one of the following: <ul style="list-style-type: none"> <li>KingFisher™ 96 KF microplate</li> <li>Tip Comb Presenting Plate for KF 96</li> <li>Nunc™ MicroWell™ 96-Well Microplate, Flat Bottom</li> <li>Nunc™ MicroWell™ 96-Well Microplate, barcoded</li> <li>ABgene™ 96-Well Polypropylene Storage Microplate</li> <li>ABgene™ 96-Well 1.2-mL Polypropylene Deepwell Storage Plate</li> <li>Nunc™ F96 MicroWell™ Black Polystyrene Plate</li> <li>Nunc™ F96 MicroWell™ White Polystyrene Plate</li> <li>KingFisher™ Deep-Well 96 Plate</li> </ul>	<ul style="list-style-type: none"> <li><a href="#">97002540</a></li> <li>267600</li> <li>167008</li> <li>269787</li> <li>AB0796</li> <li><a href="#">AB1127</a></li> <li>137101</li> <li>136101</li> <li><a href="#">95040450</a>, A48305, A48424, 95040455</li> </ul>
KingFisher™ 96 tip comb for DW magnets	97002534, A48438, A48414

(continued)

Item	Source
<b>Manual nucleic acid extraction system and materials</b>	
Magnetic Stand-96	AM10027 AM10050
Compact Digital Microplate Shaker	88882005
Incubator capable of reaching 65°C with slatted shelves	MLS
KingFisher™ Deep-Well 96 Plate	95040450, A48305, A48424, 95040455
Standard 96-well plate for the eluate, one of the following: <ul style="list-style-type: none"> <li>KingFisher™ 96 KF microplate</li> <li>MicroAmp™ Fast Optical 96-Well Reaction Plate with Barcode, 0.1 mL</li> <li>MicroAmp™ Fast Optical 96-Well Reaction Plate, 0.1 mL</li> <li>MicroAmp™ Optical 96-Well Reaction Plate with Barcode, 0.2 mL</li> <li>MicroAmp™ Optical 96-Well Reaction Plate, 0.2 mL</li> </ul>	<ul style="list-style-type: none"> <li>97002540</li> <li>4346906, 4366932</li> <li>4346907</li> <li>4306737, 4326659</li> <li>N8010560, 4316813</li> </ul>
MicroAmp™ Clear Adhesive Film	4306311
<b>Equipment</b>	
Laboratory mixer, vortex or equivalent	MLS
Single and multichannel adjustable pipettors (1.00 µL to 1,000.0 µL)	MLS
Cold block or ice	MLS
<b>Reagents</b>	
Fisher BioReagents™ Ethanol, Absolute, Molecular Biology Grade <sup>[1]</sup> , or equivalent	BP2818100, BP2818500, BP28184
Nuclease-free Water (not DEPC-Treated)	MLS
(Optional) Extraction control, if required for your assay	See the assay guide for more information
<b>Tubes, plates, and other consumables</b>	
MicroAmp™ Clear Adhesive Film	4306311
MicroAmp™ Adhesive Film Applicator	4333183
Sterile conical tubes for reagent preparation	MLS
Sterile aerosol barrier (filtered) pipette tips	thermofisher.com/pipettetips

<sup>[1]</sup> Available at [fisherscientific.com](https://www.fisherscientific.com).

# 2

## Extract nucleic acid (automated method)

■ Before you begin .....	8
■ Extract nucleic acid—Automated method (200-μL sample input volume) .....	9
■ Extract nucleic acid—Automated method (400-μL sample input volume) .....	12

Automated nucleic acid extraction is performed with the KingFisher™ Flex Magnetic Particle Processor with 96 Deep-Well Head using sample volume input of either 200 μL or 400 μL. See your assay documentation for specific sample volume input recommendations.

### Before you begin

- Ensure that you read and understand the information provided in this guide before you begin the extraction procedure.
- Review your assay documentation to determine if an extraction control is recommended to verify the efficacy of the nucleic acid preparation. Follow the extraction control guidelines provided in the assay documentation.
- Determine the number of required reactions based on the number of patient samples to be processed, plus one Negative Control per plate.
- Prepare fresh 80% Ethanol using Ethanol, Absolute, Molecular Biology Grade and Nuclease-free Water (not DEPC-Treated) for the required number of reactions, sufficient for 1 mL per reaction, plus 10% overage.
- Label the short side of each KingFisher™ Deep-Well 96 Plate (4):

Label	Number of plates
Sample plate	1
Wash 1	1
Wash 2	1
Elution plate	1

- Label the short side of the KingFisher™ 96 KF microplate (1):

Label	Number of plates
Tip comb	1

---

**Note:** The following items can be used to hold the tip comb instead of the KingFisher™ 96 KF microplate:

- Tip Comb Presenting Plate for KF 96
  - Nunc™ MicroWell™ 96-Well Microplate, Flat Bottom
  - Nunc™ MicroWell™ 96-Well Microplate, barcoded
  - ABgene™ 96-Well Polypropylene Storage Microplate
  - ABgene™ 96-Well 1.2-mL Polypropylene Deepwell Storage Plate
  - Nunc™ F96 MicroWell™ Black Polystyrene Plate
  - Nunc™ F96 MicroWell™ White Polystyrene Plate
  - KingFisher™ Deep-Well 96 Plate
- 
- Mark the Negative Control well on the plate.

## Extract nucleic acid—Automated method (200- $\mu$ L sample input volume)

### Set up the instrument (200- $\mu$ L sample input volume)

1. Ensure that the KingFisher™ Flex Magnetic Particle Processor with 96 Deep-Well Head is set up with the KingFisher™ Flex 96 Deep-Well Heating Block.

---

**IMPORTANT!** Failure to use the proper magnetic head and heat block results in lower yields and potential harm to the instrument.

---

2. Ensure that the **MVP\_2Wash\_200\_Flex** program has been downloaded from the MagMAX™ Viral/Pathogen II Nucleic Acid Isolation Kit product page at [www.thermofisher.com](http://www.thermofisher.com) and loaded onto the instrument.

## Prepare the processing plates (200-µL sample input volume)

**Note:** During the wash steps, the Wash Solution may develop inert white or brown particulates that float in solution. This is not a cause for concern and does not negatively affect performance.

Prepare the processing plates according to the following table. Cover the plates with a temporary seal (such as MicroAmp™ Clear Adhesive Film), then store at room temperature for up to 1 hour while you set up the sample plate.

Plate ID	Plate position	Plate type	Reagent	Volume per well
Wash 1 Plate	2	KingFisher™ Deep-Well 96 Plate	Wash Solution	500 µL
Wash 2 Plate	3		80% Ethanol	500 µL
Elution Plate	4		Elution Buffer	50 µL
Tip Comb Plate	5	Place a KingFisher™ 96 tip comb for DW magnets in a KingFisher™ 96 KF microplate or equivalent plate <sup>[1]</sup>		

<sup>[1]</sup> See “Before you begin” on page 8 for equivalent plates.

## Prepare Binding Bead Mix (200-µL sample input volume)

Prepare the required amount of Binding Bead Mix on each day of use.

1. Vortex the Binding Beads to ensure that the bead mixture is homogeneous.
2. For the number of required reactions, prepare the Binding Bead Mix according to the following table:

Component	Volume per well <sup>[1]</sup>
Binding Solution	265 µL
Binding Beads	10 µL
<b>Total volume per well</b>	<b>275 µL</b>

<sup>[1]</sup> Include 10% overage when making the Binding Bead Mix for use with multiple reactions.

3. Mix well by inversion, then store at room temperature.

## Prepare sample plate

1. Invert the Binding Bead Mix 5 times gently to mix, then add 275 µL to each sample well and the Negative Control well in the Sample Plate (KingFisher™ Deep-Well 96 Plate).

**Note:** Remix the Binding Bead Mix by inversion frequently during pipetting to ensure even distribution of beads to all samples or wells. The Binding Bead Mix is viscous, so pipet slowly to ensure that the correct amount is added. DO NOT reuse pipette tips to add Binding Bead Mix to the samples, as the high viscosity will cause variations in the volumes added.

2. Add 200 µL of sample to each sample well.

3. Add 200  $\mu$ L of Nuclease-free Water (not DEPC-Treated) to the Negative Control well.
4. Add 5  $\mu$ L of Proteinase K to each well, including the Negative Control well.
5. (Optional) If using an extraction control, add the required volume to each well, including the Negative Control well. For more information about using an extraction control, see the assay documentation.

---

**Note:** The Proteinase K (see step 4) and extraction control can be pre-mixed on each day of use, then kept on ice. Add the combined required volume of Proteinase K and extraction control to each well of the Sample Plate.

For example, if your assay recommends 5  $\mu$ L of the extraction control per reaction, add 10  $\mu$ L of pre-mixed Proteinase K and extraction control to each well during step 4.

---

## Process the samples (200- $\mu$ L sample input volume)

1. Select the **MVP\_2Wash\_200\_Flex** on the KingFisher™ Flex Magnetic Particle Processor with 96 Deep-Well Head.
2. Start the run, then load the prepared plates into position when prompted by the instrument.
3. After the run is complete (~22 minutes after start), immediately remove the Elution Plate from the instrument, then cover the plate with MicroAmp™ Clear Adhesive Film.

---

**IMPORTANT!** To prevent evaporation, seal the plate containing the eluate immediately.

---

The samples are eluted in 50  $\mu$ L of Elution Buffer (see “Prepare the processing plates (200- $\mu$ L sample input volume)” on page 10).

---

**Note:**

- Significant bead carry over may adversely impact RT-PCR performance. If bead carry over is observed, place elution plate on a magnetic stand to pellet the beads, then pipette the eluate to a new 96-well plate for use in real-time PCR. Review real-time PCR results to determine if re-extraction is required.
  - To ensure reliable performance of the KingFisher™ Flex Magnetic Particle Processor, perform preventive maintenance as instructed by the manufacturer.
- 

Place the Elution Plate on ice for immediate use in real-time RT-PCR.

## Extract nucleic acid—Automated method (400-µL sample input volume)

### Set up the instrument

1. Ensure that the KingFisher™ Flex Magnetic Particle Processor with 96 Deep-Well Head is set up with the KingFisher™ Flex 96 Deep-Well Heating Block.

---

**IMPORTANT!** Failure to use the proper magnetic head and heat block results in lower yields and potential harm to the instrument.

---

2. Ensure that the **MVP\_2Wash\_400\_Flex** program has been downloaded from the MagMAX™ Viral/Pathogen II Nucleic Acid Isolation Kit product page at [www.thermofisher.com](http://www.thermofisher.com) and loaded onto the instrument.

### Prepare the processing plates

---

**Note:** During the wash steps, the Wash Solution may develop inert white or brown particulates that float in solution. This is not a cause for concern and does not negatively affect performance.

---

Prepare the processing plates according to the following table. Cover the plates with a temporary seal (such as MicroAmp™ Clear Adhesive Film), then store at room temperature for up to 1 hour while you set up the sample plate.

Plate ID	Plate position	Plate type	Reagent	Volume per well
Wash 1 Plate	2	KingFisher™ Deep-Well 96 Plate	Wash Solution	1,000 µL
Wash 2 Plate	3		80% Ethanol	1,000 µL
Elution Plate	4		Elution Buffer	50 µL
Tip Comb Plate	5	Place a KingFisher™ 96 tip comb for DW magnets in a KingFisher™ 96 KF microplate or equivalent plate <sup>[1]</sup>		

<sup>[1]</sup> See “Before you begin” on page 8 for equivalent plates.

## Prepare Binding Bead Mix (400-µL sample input volume)

Prepare the required amount of Binding Bead Mix on each day of use.

1. Vortex the Binding Beads to ensure that the bead mixture is homogeneous.
2. For the number of required reactions, prepare the Binding Bead Mix according to the following table:

Component	Volume per well <sup>[1]</sup>
Binding Solution	530 µL
Binding Beads	20 µL
<b>Total volume per well</b>	<b>550 µL</b>

<sup>[1]</sup> Include 10% overage when making the Binding Bead Mix for use with multiple reactions.

3. Mix well by inversion, then store at room temperature.

## Prepare sample plate

1. Invert the Binding Bead Mix 5 times gently to mix, then add 550 µL to each sample well and the Negative Control well in the Sample Plate (KingFisher™ Deep-Well 96 Plate).

---

**Note:** Remix the Binding Bead Mix by inversion frequently during pipetting to ensure even distribution of beads to all samples or wells. The Binding Bead Mix is viscous, so pipet slowly to ensure that the correct amount is added. DO NOT reuse pipette tips to add Binding Bead Mix to the samples, as the high viscosity will cause variations in the volumes added.

---

2. Add 400 µL of sample to each sample well.
3. Add 400 µL of Nuclease-free Water (not DEPC-Treated) to the Negative Control well.
4. Add 10 µL of Proteinase K to each well, including the Negative Control well.
5. (Optional) If using an extraction control, add the required volume to each well, including the Negative Control well. For more information about using an extraction control, see the assay documentation.

---

**Note:** The Proteinase K (see step 4) and extraction control can be pre-mixed on each day of use, then kept on ice. Add the combined required volume of Proteinase K and extraction control to each well of the Sample Plate.

For example, if your assay recommends 10 µL of the extraction control per reaction, add 20 µL of pre-mixed Proteinase K and extraction control to each well during step 4.

---

## Process the samples (400-µL sample input volume)

1. Select the **MVP\_2Wash\_400\_Flex** on the KingFisher™ Flex Magnetic Particle Processor with 96 Deep-Well Head.
2. Start the run, then load the prepared plates into position when prompted by the instrument.
3. After the run is complete (~24 minutes after start), immediately remove the Elution Plate from the instrument, then cover the plate with MicroAmp™ Clear Adhesive Film.

---

**IMPORTANT!** To prevent evaporation, seal the plate containing the eluate immediately.

---

The samples are eluted in 50 µL of Elution Buffer (see “Prepare the processing plates” on page 12).

---

**Note:**

- Significant bead carry over may adversely impact RT-PCR performance. If bead carry over is observed, place elution plate on a magnetic stand to pellet the beads, then pipette the eluate to a new 96-well plate for use in real-time PCR. Review real-time PCR results to determine if re-extraction is required.
  - To ensure reliable performance of the KingFisher™ Flex Magnetic Particle Processor, perform preventive maintenance as instructed by the manufacturer.
- 

Place the Elution Plate on ice for immediate use in real-time RT-PCR.

## Extract nucleic acid (manual method)

■ Before you begin .....	15
■ Extract nucleic acid—Manual method (200- $\mu$ L sample input volume) .....	16
■ Extract nucleic acid—Manual method (400- $\mu$ L sample input volume) .....	18

Manual nucleic acid extraction can be performed using sample volume input of either 200  $\mu$ L or 400  $\mu$ L. See your assay documentation for specific sample volume input recommendations.

### Before you begin

- Ensure that you read and understand the information provided in this guide before you begin the extraction procedure.
- Review your assay documentation to determine if an extraction control is recommended to verify the efficacy of the nucleic acid preparation. Follow the extraction control guidelines provided in the assay documentation.
- Determine the number of required reactions based on the number of patient samples to be processed, plus one Negative Control per plate.
- Prepare fresh 80% Ethanol using Ethanol, Absolute, Molecular Biology Grade and Nuclease-free Water (not DEPC-Treated) for the required number of reactions, plus 10% overage.

Sample input volume	Volume of 80% Ethanol per reaction
200 $\mu$ L	0.75 mL
400 $\mu$ L	1.5 mL

- Mark the Negative Control well on the plate.

## Extract nucleic acid—Manual method (200- $\mu$ L sample input volume)

### Prepare Binding Bead Mix (200- $\mu$ L sample input volume)

Prepare the required amount of Binding Bead Mix on each day of use.

1. Vortex the Binding Beads to ensure that the bead mixture is homogeneous.
2. For the number of required reactions, prepare the Binding Bead Mix according to the following table:

Component	Volume per well <sup>[1]</sup>
Binding Solution	265 $\mu$ L
Binding Beads	10 $\mu$ L
<b>Total volume per well</b>	<b>275 <math>\mu</math>L</b>

<sup>[1]</sup> Include 10% overage when making the Binding Bead Mix for use with multiple reactions.

3. Mix well by inversion, then store at room temperature.

### Digest with Proteinase K (200- $\mu$ L sample input volume)

This section provides volumes for the sample plate. Your plate layout will depend on the number of samples you run.

1. Invert the Binding Bead Mix 5 times gently to mix, then add 275  $\mu$ L to each sample well and Negative Control well.

---

**Note:** Remix the Binding Bead Mix by inversion frequently during pipetting to ensure even distribution of beads to all samples or wells. The Binding Bead Mix is viscous, so pipet slowly to ensure that the correct amount is added. DO NOT reuse pipette tips to add Binding Bead Mix to the samples, as the high viscosity will cause variations in the volumes added.

---

2. Add 200  $\mu$ L of sample to each sample well of a KingFisher™ Deep-Well 96 Plate.
3. Add 200  $\mu$ L of Nuclease-free Water (not DEPC-Treated) to the Negative Control well.
4. Add 5  $\mu$ L of Proteinase K to each well, including the Negative Control well.
5. (Optional) If using an extraction control, add the required volume to each well, including the Negative Control well. For more information about using an extraction control, see the assay documentation.

---

**Note:** The Proteinase K (see step 4) and extraction control can be pre-mixed on each day of use, then kept on ice. Add the combined required volume of Proteinase K and extraction control to each well of the Sample Plate.

For example, if your assay recommends 5  $\mu$ L of the extraction control per reaction, add 10  $\mu$ L of pre-mixed Proteinase K and extraction control to each well during step 4.

---

6. Seal the plate with MicroAmp™ Clear Adhesive Film, then shake the sealed plate at 1,050 rpm for 2 minutes.
7. Incubate the sealed plate at 65°C for 5 minutes (ensure the bottom of the plate is uncovered), then shake the plate at 1,050 rpm for 5 minutes.
8. Place the sealed plate on the magnetic stand for 10 minutes or until all of the beads have collected.

## Wash the beads (200-µL sample input volume)

---

**Note:** During the wash steps, the Wash Solution may develop inert white or brown particulates that float in solution. This is not a cause for concern and does not negatively affect performance.

---

1. Keeping the plate on the magnet, carefully remove the cover, then discard the supernatant from each well.

---

**IMPORTANT!** Avoid disturbing the beads.

---

2. Remove the plate from the magnetic stand, then add 500 µL of Wash Solution to each sample.
3. Reseal the plate, then shake at 1,050 rpm for 1 minute.
4. Place the plate back on the magnetic stand for 2 minutes, or until all the beads have collected.
5. Keeping the plate on the magnet, carefully remove the cover, then discard the supernatant from each well.

---

**IMPORTANT!** Avoid disturbing the beads.

---

6. Repeat step 2 to step 5 using 500 µL of 80% Ethanol.
7. Repeat step 2 to step 5 using 250 µL of 80% Ethanol.
8. Dry the beads by shaking the plate (uncovered) at 1,050 rpm for 2 minutes.

## Elute the nucleic acid (200-µL sample input volume)

1. Add 50 µL of Elution Buffer to each sample, then seal the plate with MicroAmp™ Clear Adhesive Film.
2. Shake the sealed plate at 1,050 rpm for 5 minutes.
3. Place the plate in an incubator at 65°C for 10 minutes.
4. Remove the plate from the incubator, then shake the plate at 1,050 rpm for 5 minutes.
5. Place the sealed plate on the magnetic stand for 3 minutes or until clear to collect the beads against the magnets.

6. Keeping the plate on the magnet, carefully remove the seal, transfer the eluates to a fresh standard (not deep-well) 96-well plate, then seal the plate with MicroAmp™ Clear Adhesive Film.

---

**IMPORTANT!** To prevent evaporation, seal the plate containing the eluate immediately after the transfers are complete.

---

**Note:** Significant bead carry over may adversely impact RT-PCR performance. If bead carry over is observed, extend the time on the magnetic stand to further pellet the beads, then pipette the eluate to a new 96-well plate for use in real-time PCR. Review real-time PCR results to determine if re-extraction is required.

---

Place the plate on ice for immediate use in real-time RT-PCR.

## Extract nucleic acid—Manual method (400-µL sample input volume)

### Prepare Binding Bead Mix (400-µL sample input volume)

Prepare the required amount of Binding Bead Mix on each day of use.

1. Vortex the Binding Beads to ensure that the bead mixture is homogeneous.
2. For the number of required reactions, prepare the Binding Bead Mix according to the following table:

Component	Volume per well <sup>[1]</sup>
Binding Solution	530 µL
Binding Beads	20 µL
<b>Total volume per well</b>	<b>550 µL</b>

<sup>[1]</sup> Include 10% overage when making the Binding Bead Mix for use with multiple reactions.

3. Mix well by inversion, then store at room temperature.

### Digest with Proteinase K (400-µL sample input volume)

This section provides volumes for the sample plate. Your plate layout will depend on the number of samples you run.

1. Invert the Binding Bead Mix 5 times gently to mix, then add 550 µL to each sample well and Negative Control well.

---

**Note:** Remix the Binding Bead Mix by inversion frequently during pipetting to ensure even distribution of beads to all samples or wells. The Binding Bead Mix is viscous, so pipet slowly

to ensure that the correct amount is added. DO NOT reuse pipette tips to add Binding Bead Mix to the samples, as the high viscosity will cause variations in the volumes added.

---

2. Add 400 µL of sample to each sample well of a KingFisher™ Deep-Well 96 Plate.
  3. Add 400 µL of Nuclease-free Water (not DEPC-Treated) to the Negative Control well.
  4. Add 10 µL of Proteinase K to each well, including the Negative Control well.
  5. (Optional) If using an extraction control, add the required volume to each well, including the Negative Control well. For more information about using an extraction control, see the assay documentation.
- 

**Note:** The Proteinase K (see step 4) and extraction control can be pre-mixed on each day of use, then kept on ice. Add the combined required volume of Proteinase K and extraction control to each well of the Sample Plate.

For example, if your assay recommends 10 µL of the extraction control per reaction, add 20 µL of pre-mixed Proteinase K and extraction control to each well during step 4.

---

6. Seal the plate with MicroAmp™ Clear Adhesive Film, then shake the sealed plate at 1,050 rpm for 2 minutes.
7. Incubate the sealed plate at 65°C for 5 minutes (ensure the bottom of the plate is uncovered), then shake the plate at 1,050 rpm for 5 minutes.
8. Place the sealed plate on the magnetic stand for 10 minutes or until all of the beads have collected.

## Wash the beads (400-µL sample input volume)

---

**Note:** During the wash steps, the Wash Solution may develop inert white or brown particulates that float in solution. This is not a cause for concern and does not negatively affect performance.

---

1. Keeping the plate on the magnet, carefully remove the cover, then discard the supernatant from each well.
- 

**IMPORTANT!** Avoid disturbing the beads.

---

2. Remove the plate from the magnetic stand, then add 1 mL of Wash Solution to each sample.
  3. Reseal the plate, then shake at 1,050 rpm for 1 minute.
  4. Place the plate back on the magnetic stand for 2 minutes, or until all the beads have collected.
  5. Keeping the plate on the magnet, carefully remove the cover, then discard the supernatant from each well.
- 

**IMPORTANT!** Avoid disturbing the beads.

---

6. Repeat step 2 to step 5 using 1 mL of 80% Ethanol.

7. Repeat step 2 to step 5 using 500  $\mu$ L of 80% Ethanol.
8. Dry the beads by shaking the plate (uncovered) at 1,050 rpm for 2 minutes.

### Elute the nucleic acid (400- $\mu$ L sample input volume)

1. Add 50  $\mu$ L of Elution Buffer to each sample, then seal the plate with MicroAmp™ Clear Adhesive Film.
2. Shake the sealed plate at 1,050 rpm for 5 minutes.
3. Place the plate in an incubator at 65°C for 10 minutes.
4. Remove the plate from the incubator, then shake the plate at 1,050 rpm for 5 minutes.
5. Place the sealed plate on the magnetic stand for 3 minutes or until clear to collect the beads against the magnets.
6. Keeping the plate on the magnet, carefully remove the seal, transfer the eluates to a fresh standard (not deep-well) 96-well plate, then seal the plate with MicroAmp™ Clear Adhesive Film.

---

**IMPORTANT!** To prevent evaporation, seal the plate containing the eluate immediately after the transfers are complete.

---

**Note:** Significant bead carry over may adversely impact RT-PCR performance. If bead carry over is observed, extend the time on the magnetic stand to further pellet the beads, then pipette the eluate to a new 96-well plate for use in real-time PCR. Review real-time PCR results to determine if re-extraction is required.

---

Place the plate on ice for immediate use in real-time RT-PCR.

# Safety

**WARNING! GENERAL SAFETY.** Using this product in a manner not specified in the user documentation may result in personal injury or damage to the instrument or device. Ensure that anyone using this product has received instructions in general safety practices for laboratories and the safety information provided in this document.

- Before using an instrument or device, read and understand the safety information provided in the user documentation provided by the manufacturer of the instrument or device.
- Before handling chemicals, read and understand all applicable Safety Data Sheets (SDSs) and use appropriate personal protective equipment (gloves, gowns, eye protection, and so on). To obtain SDSs, see the “Documentation and Support” section in this document.

## Chemical safety

**WARNING! GENERAL CHEMICAL HANDLING.** To minimize hazards, ensure laboratory personnel read and practice the general safety guidelines for chemical usage, storage, and waste provided below. Consult the relevant SDS for specific precautions and instructions:

- Read and understand the Safety Data Sheets (SDSs) provided by the chemical manufacturer before you store, handle, or work with any chemicals or hazardous materials. To obtain SDSs, see the "Documentation and Support" section in this document.
- Minimize contact with chemicals. Wear appropriate personal protective equipment when handling chemicals (for example, safety glasses, gloves, or protective clothing).
- Minimize the inhalation of chemicals. Do not leave chemical containers open. Use only with sufficient ventilation (for example, fume hood).
- Check regularly for chemical leaks or spills. If a leak or spill occurs, follow the manufacturer cleanup procedures as recommended in the SDS.
- Handle chemical wastes in a fume hood.
- Ensure use of primary and secondary waste containers. (A primary waste container holds the immediate waste. A secondary container contains spills or leaks from the primary container. Both containers must be compatible with the waste material and meet federal, state, and local requirements for container storage.)
- After emptying a waste container, seal it with the cap provided.
- Characterize (by analysis if needed) the waste generated by the particular applications, reagents, and substrates used in your laboratory.
- Ensure that the waste is stored, transferred, transported, and disposed of according to all local, state/provincial, and/or national regulations.
- **IMPORTANT!** Radioactive or biohazardous materials may require special handling, and disposal limitations may apply.

## Biological hazard safety

**WARNING! BIOHAZARD.** Biological samples such as tissues, body fluids, infectious agents, and blood of humans and other animals have the potential to transmit infectious diseases. Conduct all work in properly equipped facilities with the appropriate safety equipment (for example, physical containment devices). Safety equipment can also include items for personal protection, such as gloves, coats, gowns, shoe covers, boots, respirators, face shields, safety glasses, or goggles. Individuals should be trained according to applicable regulatory and company/ institution requirements before working with potentially biohazardous materials. Follow all applicable local, state/provincial, and/or national regulations. The following references provide general guidelines when handling biological samples in laboratory environment.

- U.S. Department of Health and Human Services, *Biosafety in Microbiological and Biomedical Laboratories (BMBL)*, 5th Edition, HHS Publication No. (CDC) 21-1112, Revised December 2009; found at:  
<https://www.cdc.gov/labs/pdf/CDC-BiosafetyMicrobiologicalBiomedicalLaboratories-2009-P.pdf>
- World Health Organization, *Laboratory Biosafety Manual*, 3rd Edition, WHO/CDS/CSR/LYO/2004.11; found at:  
[www.who.int/csr/resources/publications/biosafety/Biosafety7.pdf](http://www.who.int/csr/resources/publications/biosafety/Biosafety7.pdf)

**WARNING! Potential Biohazard.** If you use the kit with the automated nucleic extraction workflow, the surface of the KingFisher™ purification system may be considered a biohazard. Use appropriate decontamination methods when working with biohazards.

# Documentation and support

## Related documentation

Document	Publication Number
<i>Thermo Scientific™ KingFisher™ Flex User Manual</i>	N07669

## Customer and technical support

Visit [thermofisher.com/support](http://thermofisher.com/support) for the latest service and support information.

- Worldwide contact telephone numbers
- Product support information
  - Product FAQs
  - Software, patches, and updates
  - Training for many applications and instruments
- Order and web support
- Product documentation
  - User guides, manuals, and protocols
  - Certificates of Analysis
  - Safety Data Sheets (SDSs; also known as MSDSs)

---

**Note:** For SDSs for reagents and chemicals from other manufacturers, contact the manufacturer.


---

## Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at [www.thermofisher.com/us/en/home/global/terms-and-conditions.html](http://www.thermofisher.com/us/en/home/global/terms-and-conditions.html). If you have any questions, please contact Life Technologies at [www.thermofisher.com/support](http://www.thermofisher.com/support).

**J.T. Baker**

**Ethanol**

**OPENED 8462**

For use in High Performance Liquid Chromatography.

**'BAKER HPLC ANALYZED'**

2.5 l  
Lot: 2115305855

Use Before: 2023-06  
For Laboratory, Research or  
Manufacturing use.

$C_2H_5OH$   
M: 46.07 g/mol  
Absolute

For a Certificate of Analysis, visit: [www.avantormaterials.com](http://www.avantormaterials.com)


**UN 1170**

Cas: 64-17-5  
Country of Origin: POLAND

Avantor Performance Materials Poland S.A.  
Sowińskiego 11, 44-101 Gliwice, Poland. Phone: +48 32 23 92 000

JTB: 2080100 ©198

applied biosystems

by Thermo Fisher Scientific

**MagMAX™ Viral/Pathogen Binding Beads**

Part 5 of 5 MagMAX™ Viral/Pathogen II Nucleic Acid Isolation Kit (REF A48383)

20 mL

**LOT** 01151807

**EXP** 2022-08-30

15 °C — 25 °C


CE IVD

 [www.thermofisher.com/order/catalog/product/A48383](http://www.thermofisher.com/order/catalog/product/A48383)

 Thermo Fisher Scientific Baltics UAB  
V.A. Graiciuno 8, LT-02241 Vilnius, Lithuania

Made in Lithuania  
LBL0020129

appliedbiosystems

REF A48099

by Thermo Fisher Scientific

TaqPath™ COVID-19 Control Dilution Buffer

Positive Control Dilution Buffer part of TaqPath™

COVID-19 RT-PCR Kit

Part 3 of 4 of TaqPath COVID-19 CE-IVD RT-PCR Kit (01) 1 0190302 01548 6

(REF A48067) (17) 220301

10 Tubes LOT 2107055 -30°C -10°C (10) 2107055 (240) A48099

CE IVD i 2022-03-01 EC REF Life Technologies Europe B.V.

Life Technologies Corp. Kwartsweg 2, 2665 NN

6055 Sunol Blvd. Bleiswijk, The Netherlands

Pleasanton, CA 94566 USA Read SDS Product of USA

24-hour Medical or Spill Emergency Response (Chemtrec)  
North America +1 866 536 0631 / Other +1 301 431 8585

1-10-21

invitrogen

UltraPure™  
Distilled Water  
DNase/RNase Free  
500mL


Distilled Water

REF	10977-035
LOT	2375266
	2023-08


Life Technologies Limited  
3 Fountain Drive  
Inchinnan Business Park  
Paisley, PA4 9RF, UK

Made in UK

applied biosystems

by Thermo Fisher Scientific

REF A42364LT

MagMAX™ Viral/Pathogen  
Elution Buffer

Part 4 of 5 MagMAX™ Viral/Pathogen II  
Nucleic Acid Isolation Kit (REF A48383)

100 mL

LOT 01153026

2022-09-03


[www.thermofisher.com/order/catalog/product/A48383](http://www.thermofisher.com/order/catalog/product/A48383)

Thermo Fisher Scientific Baltics UAB  
V.A. Graiciuno 8, LT-02241 Vilnius, Lithuania  
Made in Lithuania

22/11

2021

15 °C 25 °C

(01)10190302015806

(17)220903

(10)01153026

(240)A42364LT

CE IVD

Read SDS

LBL0020130 Rev.02.00

**applied biosystems**  
by Thermo Fisher Scientific

**TaqPath™ COVID-19 RT-PCR Kit**

Real-time RT-PCR test for in vitro qualitative  
detection of nucleic acid from the SARS-CoV-2

Part 1 of 4 of TaqPath™ COVID-19 CE-IVD RT-PCR Kit (REF A48067)

**REF** A48102

Σ 1000

**LOT**

2107176



2022-07-02

CE

**IVD**

-30°C  -10°C

 Life Technologies Corporation  
6055 Sunol Boulevard  
Pleasanton, CA 94566, USA

**EC REP**

Life Technologies Europe B.V.  
Kwartsweg 2, 2665 NN  
Bleiswijk, The Netherlands

**Read SDS**

Made in UK

100093686 Rev C

**applied biosystems**  
by Thermo Fisher Scientific

**TagPath™ 1-Step Multiplex Master Mix**

No ROX

Part 4 of 4 of TagPath™ COVID-19 CE-IVD RT-PCR Kit (REF A48067)

10 mL

CE IVD LOT 01124599

2022-07-08

Life Technologies Corp.  
6055 Sunol Boulevard  
Pleasanton, CA 94566 USA

Life Technologies Europe B.V.  
Kwartweg 2, 2005, NN  
Bleiswijk, The Netherlands

Read SDS

Made in Lithuania

REF A48111

QR Code

10190302015523 (01)  
17220708 (17)  
01124599 (10)  
A48111 (240)

LBLO019788 Rev.04.00

# MagMAX™ Viral/Pathogen Wash Solution

Part 2 of 5 MagMAX™ Viral/Pathogen Nucleic Acid Isolation Kit (RNA)

1000 mL

Lot 01152606

2022-09-02

Thermo Fisher Scientific  
V.A. Disinfectant & L.T.A.  
Made in Lithuania

Contains Guanidine  
Causes skin irritation. Causes  
eye irritation. Harmful if  
SWALLOWED: Immediate  
POISON CENTER or doctor  
IF ON SKIN: Wash with  
and water. IF IN EYES:  
with water for several  
contact lenses, if possible.  
Continue rinsing.  
Dispose of contents in  
approved waste disposal

Doc 54

**opened 26/11**  
**invitrogen**

by Thermo Fisher Scientific

**REF** AM9930

**LOT** 2107156

**Ambion™**

**Nuclease-Free  
Water**

30 °C

**For RNA use only**

15 °C

Not DEPC Treated

Autoclaved, 0.2 µm filtered

Qty: 500 mL

For Research Use Only. Not for  
use in diagnostic procedures.

**Read SDS**

Product of USA

**Only use in a biosafety  
cabinet cleaned w/ RNase  
zap**

Life Technologies Corporation  
2130 Woodward St. Austin, TX 78744

**MagMAX™ Viral/Pathogen Wash Solution**

Part 2 of 5 MagMAX™ Viral/Pathogen Nucleic Acid Isolation Kit (RNA)

1000 mL

Lot 01152606

2022-09-02

Thermo Fisher Scientific  
V.A. Disinfectant & L.T.A.  
Made in Lithuania

Contains Guanidine  
Causes skin irritation. Causes  
eye irritation. Harmful if  
SWALLOWED: Immediate  
POISON CENTER or doctor  
IF ON SKIN: Wash with  
and water. IF IN EYES:  
with water for several  
contact lenses, if possible  
Continue rinsing  
Dispose of contents  
Approved waste disposal

**opened 26/11**  
**invitrogen**

by Thermo Fisher Scientific

**REF** AM9930

**LOT** 2107156

**Ambion™**

**Nuclease-Free  
Water**

30 °C

**For RNA use only**

15 °C

Not DEPC Treated

Autoclaved, 0.2 µm filtered

Qty: 500 mL

For Research Use Only. Not for  
use in diagnostic procedures.

**Read SDS**

Product of USA

**Only use in a biosafety  
cabinet cleaned w/ RNase  
zap**

Life Technologies Corporation  
2130 Woodward St. Austin, TX 78744

**applied biosystems**

by Thermo Fisher Scientific

**MagMAX™ Viral/Pathogen Wash Solution**

Part 2 of 5 MagMAX™ Viral/Pathogen II Nucleic Acid Isolation Kit (REF A48383)

1000 mL

**LOT** 01152606

**EXP** 2022 - 09 - 02

[www.thermofisher.com/order/catalog/product/A48383](http://www.thermofisher.com/order/catalog/product/A48383)

Thermo Fisher Scientific Baltics UAB  
V.A. Graiciuno 8, LT-02241 Vilnius, Lithuania

Made in Lithuania

**REF** A42360LT

(01)10190302015783  
(17)220902  
(10)01152606  
(240)A42360LT

**CE** **IVD**

**WARNING.**  
Contains Guanidine Thiocyanate.  
Causes skin irritation. Causes eye irritation. Harmful if swallowed.  
**SWALLOWED:** Immediately call a POISON CENTER or doctor.  
**IF ON SKIN:** Wash with plenty of water. **IF IN EYES:** Rinse with water for several minutes. Remove contact lenses, if present, and continue rinsing. Rinse thoroughly. Do. Continue rinsing. Rinse thoroughly. Dispose of contents/containers in an approved waste disposal plan.

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Pure Medical B.V.